

X-PAK Expeditionary Personnel Aid Kit

REV012318

IDEAL FOR:

Individual Sports *Adventure Travel*
Backpacking *Surfing*
Hiking *Motocross / ATV*
Hunting *Mountain Biking*
Camping *Skiing*
Vehicles *Snowboarding*



First Aid for minor to severe injuries for multiple patients to provide emergency treatment until definitive care arrives



No Adventure Should Be Without the X-PAK...

From the mountains to the beach—no adventure should be without the X-PAK. NAR's Expeditionary Personnel Aid Kit (X-PAK) equips you to address the first aid needs of a small group with injuries ranging from major bleeding trauma and fractures to minor cuts and abrasions. The compact size enables optimized stowage for easy transport, ensuring that you can always have the necessary equipment to save a life within reach.

Don't go on an adventure without it...

Features:

- Treat minor to severe injuries – 2 to 3 patients
- Provide emergency treatment solutions until definitive care is available
- Compact and capable – ideal for adventure activities or extended travel

Dimensions:

- L 9 in. x W 5 in. x H 6.5 in.
- Weight: 2.72 lb

X-PAK Expeditionary Personnel Aid Kit

Item #	NSN#
82-0001	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

X-PAK Expeditionary Personnel Aid Kit

REV012318



Kit Contents :

MAJOR TRAUMA

- 2 x C-A-T® (Combat Application Tourniquet®) ORG
- 1 x 6 in. Responder ETD (Emergency Trauma Dressing)
- 1 x Compressed Gauze (4.5 in. X 4.1 yd)
- 1 x QuikClot Bleeding Control Dressing
- 2 x Emergency Survival Blanket- 52 in. x 84 in.

BANDAGES

- 2 x Burntec® Burn Dressing, 4 in. x 4 in. (10 cm x 10 cm)
- 20 x Flexible Fabric Bandage - 0.75 in. x 3 in.
- 2 x Wound Closure Strip (Pack of 10)
- 2 x Gauze Pad - 4 in. x 4 in.
- 2 x Gauze Pad - 2 in. x 2 in.
- 2 x Clear Adhesive Dressing - 4 in. x 4.75 in.
- 2 x Clear Adhesive Dressing - 2.375 in x X 2.75 in.
- 1 x Silver Bandage - 2 in. x 2 in.
- 1 x Silver Bandage - 4 in. x 4.5 in.
- 2 x Petrolatum Gauze 3 in. x 18 in.
- 1 x Flat Fold Surgical Tape – (4 strips, 2 in. x 12 in.)

PERSONAL PROTECTION AND ACCESSORIES

- 3 x Bear Claw Nitrile Trauma Gloves—LG (1 pair)
- 1 x Trauma Shears (lg) - 7.25 in.
- 1 x Permanent Marker- Blk, Small
- 1 x Penlight
- 1 x CPR Microshield Mask
- 1 x Splinter Forceps - 4.5 in.
- 1 x Luminous MED ID Patch

EYE CARE

- 2 x Polycarbonate Eye Shield (PES)

WOUND CARE

- 1x IV Syringe, 10 cc
- 1x Wound Cleaser, 4 oz

SPLINTING

- 1x Triangular Bandage - 40 in. x 40 in. x 56 in.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277