

M-FAK Mini First Aid Kit

NAR's most compact, versatile Mini First Aid Kit

The M-FAK, North American Rescue's Mini First Aid Kit, was designed for delivering immediate critical care for penetrating, blast or other traumatic injuries. At half the size of our popular T.O.R.K. platform, the M-FAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by First Responders operating in the line of duty. Treatment pack out is available in BLS (Basic Life Support), with Combat Gauze LE, and ALS (Advanced Life Support).

The M-FAK's super compact, 500 Denier nylon bag opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T® tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips to attach additional gear such as a TQ/Holder (sold separately) if needed.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

M-FAK Mini First Aid Kit

| Color | Configuration | Item # | NSN# |
|-------|----------------------|---------|------|
| BLK | Basic | 80-0494 | |
| COY | Basic | 80-0495 | |
| MTC | Basic | 80-0496 | |
| ODG | Basic | 80-0497 | |
| Gray | Basic | 80-0872 | |
| RGR | Basic | 80-0996 | |
| BLK | with Combat Gauze LE | 80-0498 | |
| COY | with Combat Gauze LE | 80-0499 | |
| MTC | with Combat Gauze LE | 80-0500 | |
| ODG | with Combat Gauze LE | 80-0501 | |
| Gray | with Combat Gauze LE | 80-0873 | |
| RGR | with Combat Gauze LE | 80-0998 | |
| BLK | Advanced | 80-0502 | |
| COY | Advanced | 80-0503 | |
| MTC | Advanced | 80-0504 | |
| ODG | Advanced | 80-0505 | |
| Gray | Advanced | 80-0874 | |
| RGR | Advanced | 80-0997 | |

Nylon Colors:

BLK

COY

MTC

ODG

GRY

RGR

Features:

- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access
- Multiple elastic loops for secure gear organization
- Vertical mount can be set to open left-to-right or right-to-left based on shooter preference
- Horizontal mount on a belt allows opening directly to your C-A-T®

REV010819

M-FAK Mini First Aid Kit

M-FAK CONTENT options:

BLS (Basic Life Support):



with COMBAT GAUZE:



ALS (Advanced Life Support):



Specifications:

- L 6 in. x W 3 in. x D 3 in.
- Weight: 13 oz



REV010819



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®