


MEDIC ECG[®] from Resolute Medical



Specifications:

Medic ECG[®] Handheld Unit:

- Leads: 12
- Electrode Compatibility: Standard
- Electrode Interface: Snap and Tab
- Computer Interface: Wired
- Frequency Range: 0.05 Hz - 150 Hz
- Sampling Rate: 8000 Hz [125 μs]
- Analysis Window: 16 seconds
- Report Generation: <3 seconds post analysis
- W 4.06 in. x H 6.26 in. x D 0.64 in.
- Weight: 1 lb



Lightest and smallest ECG designed for use in the field - ruggedized, portable, and easy to use

The Medic ECG[®] is the ideal solution for rapid field diagnosis of the condition of the heart. This 12-lead ECG is the first device of its kind designed for the field with ruggedization and portability in mind. It is extremely easy to operate and delivers fast results. The device's accuracy and quality exceeds industry standards, and industry-leading algorithms allow for automated interpretation of the results.

Device comes with Sure-Lock[®] Electrode Clips compatible with both snap and tab electrodes.

Features:

- Automated interpretation of ECG rhythms to alert the user of patient risks
- Results can be transferred through field communication channels
- Useful for early diagnosis of abnormal heart conditions before sending personnel to the field
- Versatile, inexpensive, and highly portable solution for base hospital ECG testing
- Exceeds industry standards for ECG accuracy and quality
- Supplied in muted mil spec colors with crushproof, waterproof, and dustproof case for easy and reliable storage

Medic ECG[®] Kits available in 3 options:

Lite: Medic ECG[®] unit ONLY (no tablet)

Compact: Medic ECG[®] unit with 7 in. Panasonic Toughpad[®] FZ-M1

Elite: Medic ECG[®] unit with 10 in. Panasonic Toughpad[®] FZ-G1

Medic ECG (Lite, Compact, or Elite)

Tablet Option	Series	Item #	NSN#
N/A	Lite	20-0083	
7 in. FZ-M1	Compact	20-0084	
10 in. FZ-G1	Elite	20-0085	



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION[®]

MEDIC ECG[®] from Resolute Medical

EVERY MedicECG[®] Kit contains:

- 1 x MedicECG[®] 12-lead ECG data acquisition unit
- 12 x Electrode Clips compatible with both snap and tab electrodes
- 12 x Electrodes
- 12 x Electrode Leads
- 1 x Charger
- 1 x USB Cable
- 1 x crushproof, waterproof case
- 1 x Nylon Bag



MedicECG[®] Lite
(no tablet)

Lite Kit Dimensions:

- W 13.5 in. x H 11.25 in. x D 5.5 in.
- Weight with hardcase: 8 lb
- Weight in Rucksack Transport: 1.5 lb



MedicECG[®] Compact
(with 7 in. FZ-M1 Toughpad[®] tablet)

Panasonic Toughpad[®] FZ-M1:

- W 7.9 in. x H 5.2 in. x D 0.7 in.
- Weight: 1.2 lb

Compact Kit Dimensions:

- W 13.5 in. x H 11.25 in. x D 5.5 in.
- Weight with hardcase: 9.2 lb
- Weight in Rucksack Transport: 2.5 lb



MedicECG[®] Elite
(with 10 in. FZ-G1 Toughpad[®] tablet)

Panasonic Toughpad[®] FZ-G1:

- W 10.6 in. x H 7.4 in. x D 0.8 in.
- Weight: 2.5 lb

Elite Kit Dimensions:

- W 13.5 in. x H 11.25 in. x D 5.5 in.
- Weight with hardcase: 10.5 lb
- Weight in Rucksack Transport: 3.8 lb



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION[®]