

EEL *Emergency Evacuation Litter · Non-Rigid Litters*

REV020218

Lightweight and Compact Solution to Perform Rapid Casualty Evacuation...

The EEL is a cost-effective casualty evacuation platform that features a full-size, durable, lightweight nylon construction with a load carriage capacity of 700 lb. It has eight easy to grasp carry handles spaced appropriately on the litter to accommodate various size rescue teams.

A streamlined storage pouch that attaches to a tactical-type vest via MOLLE/PALS-style connectors is included with the litter. When mounted between other mission critical tactical equipment such as ammo/chest pouches or on the small of the back, the quick deploy handle easily allows the EEL to be removed from its storage pouch to be put into use. The EEL requires minimal cube space on the limited area available on the assault vest, making it the ideal platform for tactical casualty evacuation.

Features:

- Cost-effective, non-rigid casualty evacuation platform designed for the tactical environment quick and effective casualty evacuation
- Quick deploy handle allows rapid removal from storage pouch
- Rugged, lightweight, Berry Amendment nylon construction
- Eight carry handles accommodate mission teams ranging from 2 to 8 rescuers
- Compact, small cube storage pouch with easy MOLLE/PALS style attachment to the assault vest
- Weight limit: 700 lb

Specifications:

- Pouch - L 9.75 in. x W 7.5 in. x D 2.5 in.
- Litter (Open) L 78 in. x W 30 in.
- Weight: 1 lb 14 oz
- Carrying Capacity: 700 lb

Available Nylon Pouch Colors:



Includes streamlined storage pouch with quick deploy handle.

Emergency Evacuation Litter

Color	Item#	NSN#
BLK	60-0047	
COY	60-0048	
MTC	60-0050	
ODG	60-0049	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®