

NAR QUIKLITTER™ LITE *Non-Rigid Litter*

REV2018-09-21



Lowest cost solution for safe & effective patient transport that fits in a cargo pocket!

Expanding on the success and feedback from the QuikLitter™ and QuikLitter™ XL, North American Rescue's QuikLitter™ Lite has added capability and reduced weight and size to be able to fit into a cargo pocket.

The QuikLitter™ Lite is an emergency carry litter that is designed for rapid deployment for the movement of casualties for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. Constructed of a lightweight nylon material with heavy-duty nylon webbing, the NAR QuikLitter™ Lite is rated to 750 lb and is latex free. It features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points..

Features:

- Perfect balance of size, weight, cost, and capability
- Features ten large carry handles with a load capacity tested to 750 lb
- Made of lightweight, heavy-duty nylon
- Vacuum-packed to conserve space and fit in a cargo pocket
- Latex free, water-resistant

Standard Size Specifications:

- Litter Bed: L 74 in. x W 27 in.
- Packaged: L 7 in. x W 6 in. x D 1.75 in.
- Weight: 12 oz
- Load rating: 750 lb

NAR QuikLitter™ Lite



Size	Item#	NSN#
Lite	60-0071	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®