

# TOENAIL REMOVAL TRAYS

Available individually or in a Pack of 10



## Conveniently packaged set of instruments ideal for Toenail Removal

The Toenail Removal Tray from North American Rescue is a conveniently packaged set of instruments in a half deep tray with a peel to open pouch for ease of access to the kit components. The Toenail Removal Tray meets the requirements for NSN 6515-01-630-8050 and is cost effective, single use kit with a 24 month shelf life. Available as a single unit or in a Pack of 10.

### Features:

- Conveniently packaged set of instruments for Toenail Removal
- Sterile
- Easy to open Pouch
- Cost Effective
- Disposable
- Available as a single unit or in a Pack of 10

### Kit Contents:

- 1 x Tray, Half Deep
- 1 x CSR Wrap, Kimguard 24 x 24 in.
- 1 x Kelly Forcep, Curved 5.5 in.
- 1 x Adson Forceps, 1 x 2 Teeth, 4.75 in.
- 1 x Kelly Forcep, Straight 5.5 in.
- 1 x Iris Scissors, Straight 4.5 in.
- 1 x Cuticle Nipper 4 in.
- 1 x Probe w/Eye 5.5 in.
- 10 x 12-Ply Gauze, 4 in. x 4 in.
- 1 x Kwiktear Header Bag, 10 in. x 15 in.

### Specifications:

- Packaged (Single Unit): L 9.72 in. x W 5.4 in. x D 2 in.
- Weight (Single Unit): 11 oz
- Weight (Pack of 10): 9 lb

### Toenail Removal Trays



Quantity	Item#	NSN#
Single	ZZ-0956	
Pack of 10	ZZ-0957	6515-01-630-8050



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**