

# SURGICAL SUTURE KITS

Available individually or in a Pack of 20



*Sutures not included*



## Conveniently packaged set of instruments for Surgical Sutures

The Surgical Suture Kit from North American Rescue is a conveniently packaged kit of instruments in a shallow tray with a peel to open pouch for ease of access to the kit components. The Sterile Surgical Set meets the requirements for NSN 6515-01-311-0360 and is cost effective, single use kit with a 27 month shelf life. Available as a single unit or in a Pack of 20. Sutures are not included.

### Features:

- Conveniently packaged set of surgical instruments for Surgical Sutures
- Sterile
- Easy to open Pouch
- Cost Effective
- Disposable
- Available as a single unit or in a Pack of 20

### Kit Contents:

- 1 x Tray, Shallow
- 1 x Adson Forceps w/Teeth
- 1 x Webster Needle Holder
- 1 x Iris Scissors, Straight
- 1 x Mosquito Hemostat, Curved
- 1 x Poly Bag w/Twist Tie
- 1 x Fenestrated Drape
- 1 x Alcohol Prep Swab
- 1 x PVP Prep Swab
- 1 x Cotton Filled Sponge
- 1 x Poly-backed Towel

### Specifications:

- Packaged (Single Unit): L 9 in. x W 5 in. x D 1 in.
- Weight (Single Unit): 6 oz
- Weight (Pack of 20): 7.5 lb

### Surgical Suture Kits



Quantity	Item#	NSN#
Single	ZZ-0958	
Pack of 20	ZZ-0961	6515-01-311-0360



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**