

EXPEDITIONARY CASUALTY RESPONSE KIT (ECRK)

REV020218

Lightweight platform designed to rapidly access and organize all your essential trauma management equipment

The Expeditionary Casualty Response Kit (ECRK) features silent, dirt and sand-resistant zippers that close over a large main compartment with six (6) segregated removable equipment pouches; 4 Bleed Pouches, 1 Airway Pouch, and 1 Miscellaneous Pouch. Additionally, a separate lid compartment secures the R.E.E.L. Non-Rigid Litter while two internal pouches on the internal lid compartment and two external pouches on each side of the bag accommodate additional trauma management essentials.

Features:

- Constructed of 500D Cordura® solution dyed nylon
- Sand- and water-resistant, self repairing zippers with silent 550 pull cords
- Separate lid compartment for the R.E.E.L. Non-Rigid Litter
- Six segregated removable pouches secured by hook and loop:
 - 4 Internal Bleed Pouches with Red pull handle
 - 1 Airway Pouch with Blue pull handle
 - 1 Accessories Pouch with Coyote pull handle
- 2 internal lid pouches to secure additional medical products
- 2 external side pouches on either side of the kit for additional medical products
- 2 external straps with dual adjust MILSPEC IR reduced buckles to secure the kit in a vehicle or boat

Specifications:

- Case closed: H 17 in. x W 16.5 in. x D 10 in.
- Weight: 18 lb

Option	Item#	NSN#
Kit	80-0889	
Bag Only	80-0888	

Available Color:

COY



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

EXPEDITIONARY CASUALTY RESPONSE KIT (ECRK)

REV020218

ECRK Contents:

- 6 x Combat Application Tourniquet, (C-A-T®)
- 1 x Blizzard Blanket
- 4 x Nasopharyngeal Airway, 28F with lubricant
- 4 x ARS® Needle Decompression Device, 14g, 3.25 in.
- 4 x Combat Gauze, Z-Fold Hemostatic Dressing
- 6 x Emergency Trauma Dressing, ETD™, 6 in.
- 10 x NAR Compressed Gauze, 4.5 in. x 4.1 yd
- 4 x USMC HyFin® Chest Seal Combo Pk (1 HyFin® Vent & 1 HyFin® Chest Seal)
- 6 x Pair Bear Claw™ Gloves, Large Tan
- 6 x Emergency Trauma Dressing, ETD™ ABD 8 x 10 in.
- 1 x Minor BurnTec Kit
- 1 x Trauma Shears, Large, 7.25 in.
- 6 x Combat Casualty Response Cards
- 2 x Surgical Tape 1 in.
- 1 x R.E.E.L. Non-Rigid Litter
- 10 x Antiseptic Wipes
- 1 x Combat Casualty Reference Card
- 2 x SAM Splint II
- 2 x Mini Sharpie, Black
- 3 x Elastic Bandage 6 in.
- 2 x Triangular Bandage, 45 in. x 45 in. x 63 in.
- 2 x Mini Duct Tape, 2 in. x 100 in.
- 1 x Model 7 Rescue Hook
- 6 x Polycarbonate Eye Shield
- 1 x Expeditionary Casualty Response Bag



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®