

TORK™ - Red (Tactical Operator Response Kit)

REV012918



Point-of-wounding IFAK designed for the Assault Vest or Body Armor

The Tactical Operator Response Kit™ (TORK™) by North American Rescue® was designed to meet the needs of our customers for a compact individual first aid kit which carries all the necessary equipment to address the top leading causes of potentially preventable death in a tactical environment.

The TORK™ contains essential life-saving equipment that can be used at the point of wounding for self-aid/buddy-aid to address injuries caused by penetrating or blast trauma. It's sleek, compact size allows for easy attachment using MOLLE/PALS style connectors to the assault vest/body armor or NAR's Quick Detach Panel (available for separate purchase). The TORK™ also has a rapid opener with NAR's signature Red-Tip Technology® for fast and easy, one-pull tab opening of the kit.

Features:

- Essential life-saving tools to address injuries caused by penetrating or blast trauma
- Compact individual first aid kit providing critical support for immediate Self-Aid/Buddy-Aid at the point of wounding
- Rugged Red Nylon Bag with easy MOLLE/PALS style attachment to the assault vest or body armor
- Also attaches to NAR's Quick Detach Panel (sold separately)

Specifications:

- Cased Closed:
H 7 in. x W 4.5 in. x D 2.75 in.
- Weight: 1 lb 5 oz

RED

TORK™ - Red (Tactical Operator Response Kit)

| Level | Hemostatic | Item# | NSN# |
|----------|-----------------|---------|------|
| Basic | None | 85-1805 | |
| Advanced | Combat Gauze LE | 85-1806 | |
| Advanced | ChitoGauze PRO | 85-1807 | |
| Advanced | CELOX | 85-1808 | |



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

TORK™ - Red *(Tactical Operator Response Kit)*

REV012918



BLS (Basic Life Support) Kit Contents:

- 1 x Red Nylon Bag (TORK™)
- 2 x pair Bear Claw™ Nitrile Glove, Lg.
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x C-A-T® (Combat Application Tourniquet®) Blk
- 2 x S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x 6 in. Responder ETD, Roll
- 1 x pair Trauma Shears (6.25 in.)
- 1 x Polycarbonate Eye Shield



ALS (Advanced Life Support) Kit Contents:

- 1 x Red Nylon Bag (TORK™)
- 2 x pair Bear Claw™ Nitrile Glove, Lg.
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x C-A-T® (Comba Application Tourniquet®) Blk
- 1 x Your Choice of Hemostatic: Combat Gauze LE, Celox Gauze, -or- ChitoGauze PRO
- 1 x S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x 6 in. Responder ETD, Roll
- 1 x pair Trauma Shears (6.25 in.)
- 1 x Polycarbonate Eye Shield

*Streamlined and compact...
designed to MOLLE to an assault vest or body armor*



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®