

NAR K-9 LITTER & Carrying Strap Accessories

REV051018

Designed for the safe & efficient movement of an injured canine...

The NAR K9 Litter is a lightweight, compact litter designed for the movement of injured canines. The nylon litter can be packed to store in a cargo pocket or pouch and can be quickly deployed. The litter has leg cut outs to place the canine in the litter for safe movement. The K9 Litter also has six large and four small carry handles for balanced carrying and allows for one to four person carries. Accessory options include an additional K-9 Carry Strap & 2 Carabiners to aid in securing an over the shoulder carry of an injured canine using the NAR K-9 Litter. Products may be purchased individually or in an all inclusive kit.

Features:

- Lightweight, compact litter for the movement of an injured canine
- Litter can be carried in cargo pocket or pouch
- Leg cut outs to safely carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Optional K-9 Carry Strap & 2 Carabiners for over the shoulder carry of injured K-9 also available

Dimensions:

- Litter: W 35 in. x H 48 in.
- Adjustable Carry Strap: L 9 to 54 in. x W 2 in.
- Carabiner open gate clearance: 27 mm
- Packaged: W 5.5 in. x H 10 in. x D 2.5 in.
- Weight: 1.38 lb

NAR K-9 Litter & Accessories

Options	Item #	NSN#
K-9 Litter Only	60-0074W	
Carry Strap & Carabiners	ZZ-1004	
K-9 Litter with Carry Strap & Carabiners	60-0076	

Carry Strap Accessory features:

- 2 in. wide adjustable seatbelt strap with durable plastic buckle to aid in securing an over the shoulder carry of an injured canine using the NAR K-9 Litter
- Includes 2 full-size wire gate carabiners with I-Beam construction, easy to handle, grippy shape and clean nose profile to reduce snagging



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A
MISSION®