

# LEATHER C-A-T® HOLDER *with Strap Attachment*

REV110218

## *Designed to meet the specific duty gear needs of fire and EMS personnel...*

A tourniquet needs to be accessible and readily available for one-handed deployment in under a minute, 30 seconds if possible. To meet the specific duty gear needs of fire and EMS personnel, North American Rescue now offers leather C-A-T Holders for use on a radio strap and designed to fit within organizational uniform standards. Constructed of full grain leather with high quality hook-and-loop and snaps, this premium tourniquet holder securely fastens to your radio strap with easy one-handed access to allow you to take control of extreme medical situations in a matter of seconds. It's convenient attachment allows the holder to be added to your radio strap without removing the radio pouch from the strap and is ideal for carrying your C-A-T® Gen 7 but can also accommodate other popular windlass tourniquets as well.

### Features:

- Designed to meet the specific duty gear needs of fire and EMS personnel - for use on a radio strap and designed to fit within organizational uniform standards
- Securely fastens to your radio strap with easy one-handed access to allow you to take control of extreme medical situations in a matter of seconds
- Constructed of flat black leather with high quality hook-and-loop and snaps to ensure secure placement on your radio strap to ease accessibility to your tourniquet
- Convenient attachment allows the holder to be added to your radio strap without removing the radio pouch from the strap
- Ideal for carrying the C-A-T® Gen 7 but can also accommodate other popular windlass tourniquets

### *Leather C-A-T® Holder with Strap Attachment*

Style	Item #	NSN#
Full Grain Black	30-0196	

**C-A-T® Gen 7 (sold separately)**

### Dimensions:

- L 5.75 in. x W 2 in. x D 2 in.
- Weight: 3.2 oz



***Allows the holder to be added to your radio strap without removing the radio pouch from the strap...***



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**