

NAR RESPONDER QUIKLITTER™ Non-Rigid Litter

REV2018-09-21



Lowest cost solution for safe & effective patient transport

North American Rescue's Responder QuikLitter™ is a cost effective emergency carry litter designed for rapid deployment for the movement of casualties and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements.

Constructed of a latex free water-resistant, woven polyethylene with signature red nylon webbing & stitching, the Responder QuikLitter™ is rated for up to 500 lb. and features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points.

Features:

- Cost effective emergency carry litter
- Features ten large carry handles for ease of patient handling/transport
- Load capacity tested to 500 lb
- Latex free, Water-resistant
- Signature red nylon webbing & stitching

Standard Size Specifications:

- Litter Bed: L 74 in. x W 28 in.
- Packaged: L 11 in. x W 7.75 in. x D 2 in.
- Weight: 13.3 oz
- Payload: 500 lb

NAR Responder QuikLitter™



Size	Item#	NSN#
Standard	60-0080	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®