

MINI MEDIC KIT

REV010919



Compact, multi-purpose trauma aid bag for ultimate versatility at the point of wounding...

The Mini Medic Kit is a compact trauma aid bag designed to provide various capabilities for different mission profiles. External pouches and side/bottom MOLLE allow for external tourniquet or other attachment options. Large pouch on back panel includes a drop-down pouch with abrasive resistant padded bottom allows storage of compact breaching tools. External expansion zipper allows for increased load capacity if needed. The inside of the bag has two small and one large, clear tear-away pouches for medical or other equipment supply storage. Along with two mesh inner pouches and an internal pouch on the inner lid, many storage options exist. The Mini Medic Kit provides a compact, versatile bag equipped with the essential medical equipment to treat multiple casualties.

Features:

- Compact, multi-purpose trauma aid bag with vital medical equipment for point of wounding treatment of multiple injuries
- External storage options include 2 pockets, 2 tourniquet holders & extra drop-down pouch for breaching tools or hydration pack (not included)
- Expansion zipper increases load capacity
- MOLLE capability on the sides, bottom and inside back of the bag
- Internal pouch storage includes 3 clear tear-away, 2 main compartment mesh pouches and 1 larger inner lid pouch

Specifications:

- H 14 in. x W 11 in. x D 5 in.
(Depth expandable to 7.25 in.)
- Weight: 12.5 lb

**Nylon
Colors:**

BLK

COY

ODG

RGR

Mini Medic Kits

Color	Level	Item#	NSN#
Black	Basic	80-0979	
Coyote	Basic	80-0980	
OD Green	Basic	80-0981	
Ranger Green	Basic	80-1003	
Black	Combat Gauze	80-0985	
Coyote	Combat Gauze	80-0986	
OD Green	Combat Gauze	80-0987	
Ranger Green	Combat Gauze	80-1005	
Black	Advanced	80-0982	
Coyote	Advanced	80-0983	
OD Green	Advanced	80-0984	
Ranger Green	Advanced	80-1004	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

MINI MEDIC KIT

Every **MINI MEDIC KIT** includes these core items:



PLUS these Individual Components based upon configuration of choice:



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

REV010919