

K-9 RUF KIT

REV062519

Slim profile trauma aid pack designed to carry a minimum of two complete sets of IFAK contents along with K9 medical supplies

The NAR K9 RUF Kit is a slim profile trauma aid pack designed for the K9 Handler or Tracking Team member. The K9 RUF carries a minimum of two complete sets of IFAK contents along with K9 medical supplies. Designed with two tear away internal panels that hold standard IFAK contents and two tear away clear accessory pouches for additional supplies. The internal kit also has a shock cord panel designed to hold a compact, K9 litter. The kit can be worn with a single bandolier style sling strap as well as comes with a waist strap that can be used to secure kit while tracking (waist strap removeable based on user's needs). The slim profile of the case along with the amount of contents of life-saving trauma equipment makes it a necessity for your handler, K9 and tracking team members.

Features:

- Slim profile
- Can be worn with single bandolier style sling strap as well as with a detachable waist strap
- Two tear away internal panels designed to hold standard IFAK contents
- Two additional tear away clear accessory pouches for additional equipment
- Kit contains life-saving trauma equipment and K9 specific equipment
- Internal shock cord area to hold K9 Litter
- External MOLLE on sides of bag for external tourniquet pouches
- Large panel on exterior of kit for ID patches

Specifications:

- Dimensions: H 18 in. x W 8 in. x D 2 in.
- Kit Weight: 6.3 lb
- Bag Weight: 2 lb

K-9 RUF Kit

BLK

Level	Item#
with S-Rolled Gauze	80-1027
with Combat Gauze LE	80-1026



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

K-9 RUF KIT

Can be worn with single bandolier style sling strap as well as with a detachable waist strap



PLUS+

**S-ROLLED
GAUZE**
also includes:



QTY
2 Gauze
S-Rolled

OR

**COMBAT
GAUZE LE**
also includes:



QTY
2 Combat
Gauze® LE

Every K-9 RUF KIT includes these core items:

 QTY 2 C-A-T® Tourniquet	 QTY 2 Gauze S-Rolled	 QTY 2 ETD™ 6 in. Dressing	 QTY 1 ETD™ Abdominal	 QTY 2 HyFin® Vent Compact, Twin
 QTY 2 Gauze, Petrolatum	 QTY 2 ARS® Needle, 14g	 QTY 2 Survival Blanket, Org	 QTY 1 NAR K-9 Litter	 QTY 8 Alcohol Pad
 QTY 1 Splinter Forceps, 4.5"	 QTY 1 Veterinary Tape, 2 in.	 QTY 1 Veterinary Thermometer	 QTY 1 Tissue Adhesive, Vet	 QTY 2 PES Eye Shield
 PAIR 4 Nitrile Gloves, Lg	 QTY 1 Responder Trauma Shears, Lg	 QTY 2 Permanent Marker, Large	 QTY 1 Medical Cross Patch	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®