

NAR PPE KIT

REV071219



NAR PPE Kit

Item#	NSN#
70-0007	6515-01-537-1534

Lightweight Kit, State of the Art Protection...

The NAR PPE Kit is the essential personal protection for anyone who is responsible for indigenous healthcare, detainee processing or disaster response where various hazards and contaminants may be encountered.

Our kit features the state of the art NAR Respirator and Bear Claw Nitrile Gloves in a compact lightweight kit.

The NAR PPE Kit is supplied in a vacuum sealed, rugged, easy to open package that allows for quick and effective application in the most extreme conditions.

Features:

- Certified to the NIOSH N95 standard for disposable particulate respirators
- Secure, comfortable fit with adjustable nose bridge and low breathing resistance
- Reduced packaging volume for easy storage
- Supplied with two individually wrapped pairs of Bear Claw™ Nitrile Gloves (Size: Large)

Dimensions:

- Packaged: L 3.5 in. x W 4 in. x H 1.5 in.
- Weight: 2 oz

Kit Contents:

- 1x NAR Respirator
- 4 x Bear Claw™ Nitrile Trauma Gloves (2 pr.)
- 1 x Plastic Vacuum Bag



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**