

CAT® HOLDER

REV072319

Fits Securely and Easily on the Belt...

Whether you are a member of public safety or a conscientious, situationally aware civilian who's in charge of their own Personal Security Detail, one of the most neglected Every Day Carry pieces of equipment is a combat proven windless-style tourniquet.

North American Rescue has partnered with FieldCraft Survival to fill that gap by offering a tourniquet holster for the GEN-7 CAT that is not only functional, allowing fast, consistent deployment of the tourniquet, but it's comfortable too.

This unique Inside Waist Band CAT Holder conveniently clips inside your waistband or outside on your pants belt loop, board shorts, or even on your MOLLE gear.

Features:

- Tailored specifically for your C-A-T®, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- 1000D IR Signature Reduced Nylon
- Easy-open pull elastic tab with tourniquet ID patch – "TQ"
- NATO stock numbers available for COY and ODG. Contact Customer Service for details.

Dimensions:


- L 6.25 in. x W 1.75 in. x D 1.5 in.
- Max belt width for horizontal mount: 1.5 in.
- Weight: 1.26 oz



C-A-T® Holder

Color	Item #	NSN#
BLK	30-0057	
COY	30-0044	NATO NSN AVAIL
RED	30-0128	
ODG	30-0053	NATO NSN AVAIL
MTC	30-0059	
RTC	30-0207W	
RGR	30-0208	

Nylon Colors:



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277