

READY HEAT BLANKET

REV102518



10 hours of continuous dry heat...

The ReadyHeat blanket utilized proprietary technology and all natural ingredients to deliver continuous and consistent warming. This product is part of US military protocol and carried by every US Army Medic. It can be used as an evacuation litter in an emergency situation, as it can hold a patient up to 250 lb. The Ready Heat Blanket's self-heating, 4-cell liner provides 10 hours of continuous dry heat up to 125°F (52°C). **WARNING: DO NOT PLACE DIRECTLY ON EXPOSED SKIN.**

Features:

- Blanket begins heating immediately upon opening, no activation required
- Can be used as an evacuation litter for a patient up to 250 lb
- Carried by every US Army Medic
- Approved for use on fixed wing and rotary aircraft by the US military
- Self-heating, 4-cell liner
- Provides 10 hours of continuous dry heat up to 125°F (52°C)

Specifications:

- Packaged: W 11 in. x H 13 in. x D 0.5 in.
- Blanket: 48 in. x 36 in.
- Weight: 1 lb 11 oz

Ready Heat Blanket

Item #	NSN#
ZZ-0052	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®