

# RESPONDER *Training ETDs*

REV112019



4 in. Dressing and Package

## Responder Training ETDs

Width	Item #	NSN#
4 in.	30-0145	



## ***Ideal low cost solution to help train students in the application of immediate direct pressure...***

The Responder Training ETD's from North American Rescue are non-sterile, elastic pressure bandages that incorporate the same form factor as our standard ETDs. These features include a built-in, non-adherent pad, easy-to-use hook & loop securing device, quick grip roll control system to prevent accidental unrolling if dropped and vacuum-sealed, low-cube packaging with easy-open tear notches. By becoming familiar with all of the same features, students can not only learn the procedures but also develop the muscle memory needed to complete the task more quickly and efficiently. Responder Training ETDs are available in 4 in. wide rolled configuration.

### **Features:**

- Economic, low cost solution for training in the application of immediate direct pressure
- Non-sterile, elastic pressure bandage ideal for use in training scenarios
- Same form factor as our standard NAR Responder Emergency Trauma Dressing to aid in the development of muscle memory during skills training
- Available in 4 in. width

### **Dimensions:**

#### **4 Inch ETD:**

- Bandage: L 60 in. x W 4 in.
- Pad: H 6 in. x W 4 in.
- Packaged: L 4.5 in. x W 2.5 in. x D 1.75 in.
- Weight: 2.0 oz



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**