

# BEAR CLAW™ GLOVES

REV053018

Sold on bulk box or  
in pre-packaged  
pairs of 25



## High performance textured gloves for improved grip in wet and dry situations

Bear Claw® Ultimate Nitrile gloves are made of 100% nitrile material which is free of the allergens found in standard latex gloves. They are engineered to give our users high performance by offering both maximum protection and unsurpassed tactile sensitivity.

This careful attention to detail includes textured fingertips to help with fine motor skills, beaded cuffs for quick and easy donning and a thickened palm for durability and reliability. These ambidextrous gloves are packaged 50 pairs per box (Extra large, 45 pairs per box); 25 pairs per bag. Designed for single-use only.

### Features:

- 100% nitrile gloves, textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Extended, beaded cuff for quick and easy donning
- Available in Sand color
- Form fitting, Non-sterile, Available in sizes S, M, L and XL



### Glove Specs:

- Length: 11.8 in.
- Cuff Thickness (mils): 4.7
- Palm Thickness (mils): 6.0
- Finger Thickness (mils): 8.0

### Bulk Dimensions:

- (Box) L 11 in. x W 5 in. x D 3.75 in.
- Weight: 2 lb 4 oz

### Kit Dimensions:

- (Bag) L 4.75 in. x W 6 in. x D 4 in.
- Weight: 16.1 oz

Size	Unit	Qty.	Item#	NSN#
S	Box	50 pr.	70-0009	6515-01-600-6059*
S	Kit	25 pr.	ZZ-0246	6515-01-619-5350
M	Box	50 pr.	70-0010	6515-01-600-6057
M	Kit	25 pr.	ZZ-0247	6515-01-619-5689
L	Box	50 pr.	70-0011	6515-01-600-6039
L	Kit	25 pr.	ZZ-0248	6515-01-600-0262
XL	Box	45 pr.	70-0012	6515-01-600-6052
XL	Kit	25 pr.	ZZ-0249	6515-01-619-5367

\*NATO NSN



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**