

RAPTOR IV SECURING DEVICE

REV010920



**Quickly & Efficiently
Secure IV Lines, even in
Wet or Damp Conditions...**



Engineered specifically to quickly secure an IV line in the most difficult conditions. Once vascular access is achieved, the Raptor™ secures the IV tubing to the patient without tape.

The Raptor™ reduces equipment Tactical Healthcare professionals must carry and by securing IV lines without the use of adhesive material - the Raptor™ can be used successfully on wet skin and in wet environments.

Features:

- Minimizes equipment for Tactical Healthcare Professionals
- Works well in conjunction with Saline Locks
- Secures IV tubing on wet skin
- Sold 6 per box

Specifications:

- Packaged (Individual): L 3.75 in. x W 1 in.
- Device Length (Fully Extended): L 19 in.
- Weight: 1.9 oz each
- Packaged (Box of 6): L 6 in. x W 4.8 in. x D 0.5 in.

| Unit | Item# | NSN# |
|-----------|---------|------------------|
| Pack of 6 | 30-0006 | 6515-01-523-3317 |



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®