

# EFAK EXPEDITIONARY FIRST AID KIT

REV020620



## EFAK Expeditionary First Aid Kit

Color	Item #	NSN#
BLK	80-1044	
COY	85-2533	

## Available Colors:

BLK COY

The EFAK, North American Rescue's Expeditionary First Aid Kit, was designed to be the most compact, multi-use IFAK for delivering immediate critical care for penetrating, blast or other traumatic injuries. At half the size of our popular T.O.R.K. platform, the EFAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by Expeditionary Military First Responders operating in the line of duty. Treatment pack out is the Advanced version meeting the Expeditionary Military First Responder Requirements.

The EFAK's super compact, 500 Denier nylon bag opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference. The TQ/Holder is attached to the (3) three MOLLE strips on the opposite side of the zipper opening. The EFAK also comes with an additional TQ Holder and Tourniquet for primary placement separate from the EFAK, meeting the (2) Tourniquet Expeditionary Military requirement. Also included is a loop patch on front of the bag for custom user labeling.

## Features:

- NAR's most compact, versatile Expeditionary First Aid Kit
- Designed to meet the Expeditionary Military Forces First Responder requirements for critical point-of-wounding, treating penetrating, blast or other traumatic injuries in the line of duty
- Super compact, rugged nylon platform that allows attachment both vertically and horizontally
- Treatment pack out option: Advanced
- Smaller, more compact size for ease of placement with other Military gear on the First Responder
- Half the size of our popular T.O.R.K platform




**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# EFAK EXPEDITIONARY FIRST AID KIT

REV020620



*Our most compact, half size First Aid Kit empowering Expeditionary Military First Responders with the most requested point-of-wounding equipment in the smallest cube space possible.*

## Specifications:

- L 6 in. x W 3 in. x D 3 in.
- Weight: 1 lb 3 oz

## Kit Contents:

- 1 x C-A-T® (Combat Application Tourniquet®) Black
- 1 x Tourniquet Holder (Color matches bag) externally MOLLE
- 1 x 4 in. Flat ETD™
- 2 x S-rolled Gauze (4.5 in. x 4.1 yd)
- 1 x HyFin® Vent Chest Seal, Twin Pack
- 1 x pair Bear Claw® Nitrile Trauma Gloves, Lg.
- 1 x Nasal Pharyngeal w' Lube – 28 F
- 1 x Needle Decompression (ARS) (14 GA) x 3.25"
- 1 x Eye Shield, Polycarbonate – ORG
- 1 x Hemostatic Bandage, Combat Gauze w' Z-Fold
- 1 x Tape, Duct-Mini – 2" x 100"
- 1 x Marker, Permanent – BLK, Lg
- 1 x Card, Combat Casualty Document Tool

## Separate Components:

- 1 x C-A-T® (Combat Application Tourniquet®) Black (separate - not in kit)
- 1 x Tourniquet Holder (Color matches bag) externally MOLLE



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277