

BEAR CLAW™ GLOVES

REV03062020

Sold on bulk box or
in pre-packaged
pairs of 25



High performance textured gloves for improved grip in wet and dry situations

Bear Claw® Ultimate Nitrile gloves are made of 100% nitrile material which is free of the allergens found in standard latex gloves. They are engineered to give our users high performance by offering both maximum protection and unsurpassed tactile sensitivity.

This careful attention to detail includes textured fingertips to help with fine motor skills, beaded cuffs for quick and easy donning and a thickened palm for durability and reliability. These ambidextrous gloves are packaged 50 pairs per box (Extra large, 45 pairs per box); 25 pairs per bag. Designed for single-use only.

Features:

- 100% nitrile gloves, textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Extended, beaded cuff for quick and easy donning
- Available in Sand color
- Form fitting, Non-sterile, Available in sizes S, M, L and XL

Glove Specs:

- Length: 11.8 in.
- Cuff Thickness (mils): 4.7
- Palm Thickness (mils): 6.0
- Finger Thickness (mils): 8.0

Bulk Dimensions:

- (Box) L 11 in. x W 5 in. x D 3.75 in.
- Weight: 2 lb 4 oz

Kit Dimensions:

- (Bag) L 4.75 in. x W 6 in. x D 4 in.
- Weight: 16.1 oz



| Size | Unit | Qty. | Item# | NSN# |
|------|------|--------|---------|-------------------|
| S | Box | 50 pr. | 70-0009 | 6515-01-600-6059* |
| S | Kit | 25 pr. | ZZ-0246 | 6515-01-619-5350 |
| M | Box | 50 pr. | 70-0010 | 6515-01-600-6057 |
| M | Kit | 25 pr. | ZZ-0247 | 6515-01-619-5689 |
| L | Box | 50 pr. | 70-0011 | 6515-01-600-6039 |
| L | Kit | 25 pr. | ZZ-0248 | 6515-01-600-0262 |
| XL | Box | 45 pr. | 70-0012 | 6515-01-600-6052 |
| XL | Kit | 25 pr. | ZZ-0249 | 6515-01-619-5367 |

*NATO NSN



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®