

NAR RESPIRATOR



Sealed edge to avoid open layers and adjustable nose bridge for a secure, comfortable fit.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

Front-line Respiratory Protection...

The NAR Respirator provides front-line protection for First Responders with an advanced particulate barrier.

Low breathing resistance for increased wearer comfort. Foldable design and sealed edge avoids open layers and eliminates the possibility of irritation.

Latex-free synthetic rubber head straps, adjustable nose clip and soft foam nosepiece for comfortable custom fit.

Exhalation valve reduces hot air buildup and makes it easier to breathe in hot/humid environments.

Features:

- Certified to the NIOSH N95 standard for disposable particulate respirators
- Secure, comfortable fit with adjustable nose bridge
- Minimal breathing resistance
- Reduced packaging volume for easy storage

Dimensions:

- Packaged: L 5 in. x W 4.25 in. x H 0.75 in.
- Mask size: Adult
- Weight: 0.9 oz

NAR Respirator

Item #	NSN#
70-0006	6515-01-537-1531

The NAR Respirator has an advanced particulate barrier and low-breathing resistance.

PRODUCTS WITH A
MISSION®