

# BLUE RESPONDER GLOVE KITS (Pack of 25)

REV042420



*25 conveniently  
pre-packaged pairs*



## Blue Responder Glove Kits (Pack of 25)

Size	Item#
S	ZZ-1059
M	ZZ-1060
L	ZZ-1043
XL	ZZ-1061



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

## High performance textured gloves for improved grip in wet and dry situations

Responder Nitrile Examination Gloves are made of 100% Nitrile material which is free of the allergens found in standard latex gloves. They are engineered to give our users high performance by offering both maximum protection and unsurpassed tactile sensitivity. This careful attention to detail includes textured fingertips to help with fine motor skills, beaded cuffs for quick and easy donning and a thickened palm for durability and reliability. These powder free, ambidextrous gloves come in 25 pre-packaged pairs per bag. Designed for single-use only.

### Features:

- Individually wrapped pairs for convenient rapid deployment
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Extended, beaded cuff for quick and easy donning
- Powder free, Form fitting
- Blue color
- Available in sizes S, M, L and XL
- Non-sterile, Single Use Only

### Glove Specs:

- Length: 11.8 in.
- Cuff Thickness (mils): 4.7
- Palm Thickness (mils): 6.0
- Finger Thickness (mils): 8.0

### Dimensions:

- Packaged (25 pairs):  
L 4.75 in. x W 6 in. x D 4 in.
- Individual pair:  
L 2.5 in. x W 1 in. x D .8 in.
- Weight: 16.1 oz

**PRODUCTS**  
WITH A  
**MISSION®**