

Public Access **INDIVIDUAL BLEEDING CONTROL KITS**

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control ID Patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.5 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

Essential equipment needed for persons to take action as immediate responders in helping to stop life threatening bleeding in Public Areas...



Every **INDIVIDUAL BLEEDING CONTROL KIT** contains:



PLUS these Individual Components based upon **YOUR** configuration choice:



PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KITS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	80-0452
Intermediate	Red Nylon	16.5 oz	80-0453
Advanced BCD	Red Nylon	17.4 oz	80-0523
Advanced	Red Nylon	17.4 oz	80-0454
Basic	Vacuum Sealed	11.7 oz	80-0465
Intermediate	Vacuum Sealed	13.25 oz	80-0466
Advanced BCD	Vacuum Sealed	14 oz	80-0522
Advanced	Vacuum Sealed	14 oz	80-0467

Public Access

BLEEDING CONTROL MULTI-PACKS

5-PACKS Vacuum

- Five (5) Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon tote
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced BCD & Advanced

SPECIFICATIONS

5-Pack:
H 7.5 in. x W 13 in. x D 8 in.

Each Individual Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in



PUBLIC ACCESS BLEEDING CONTROL 5-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Vacuum Sealed	5.7 lb	80-0573
Intermediate	Vacuum Sealed	5.75 lb	80-0574
Advanced BCD	Vacuum Sealed	6.5 lb	80-0576
Advanced	Vacuum Sealed	6.5 lb	80-0575

TWIN PACKS

- Contents of 2 Individual Kits
- Compact Nylon Bag designed to fit next to an AED in the cabinet

SPECIFICATIONS

Nylon Bag:
H 10 in. x W 9 in. x D 3 in.



TWIN PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	1.66 lb	80-0951
Intermediate	Red Nylon	1.8 lb	80-0952
Advanced BCD	Red Nylon	1.8 lb	80-0953

8-PACKS Nylon | Vacuum Sealed

- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon carrying case
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced & Advanced BCD

SPECIFICATIONS

8-Pack:
H 14 in. x W 11 in. x D 8 in. Each Individual Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in.

PUBLIC ACCESS BLEEDING CONTROL 8-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	8.5 lb	80-0460
Intermediate	Red Nylon	9.3 lb	80-0461
Advanced BCD	Red Nylon	9.7 lb	80-0525
Advanced	Red Nylon	9.7 lb	80-0462
Basic	Vacuum Sealed	6.8 lb	80-0468
Intermediate	Vacuum Sealed	7.6 lb	80-0469
Advanced BCD	Vacuum Sealed	8 lb	80-0524
Advanced	Vacuum Sealed	8 lb	80-0470

Every **BLEEDING CONTROL TWIN PACK** contains:



Clear Polycarbonate WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 5 & 8 Pack Stations also include an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured casualties.

SPECIFICATIONS

5 & 8 Pack Stations: H 18 in. x W 15.5 in. X D 10 in.
Empty Case Wt: 15.0 lbs (includes rack)
Large Case Item# ZZ-0822

Individual Stations: H 10 in. x W 10 in. X D 5.25 in.
Empty Case Wt: 2.39 lb
Small Case Item# ZZ-1000

8-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	80-0471
Intermediate	8	Vacuum	22.6 lb	80-0472
Advanced BCD	8	Vacuum	23 lb	80-0526
Advanced	8	Vacuum	23 lb	80-0473
Basic	8	Nylon	23.5 lb	80-0474
Intermediate	8	Nylon	24.3 lb	80-0475
Advanced BCD	8	Nylon	27.4 lb	80-0527
Advanced	8	Nylon	27.4 lb	80-0476

5-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	20.7 lb	80-0581
Intermediate	5	Vacuum	20.7 lb	80-0582
Advanced BCD	5	Vacuum	21.5 lb	80-0584
Advanced	5	Vacuum	21.5 lb	80-0583

IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	3.3 lb	80-0938
Intermediate	1	Vacuum	3.42 lb	80-0939
Advanced BCD	1	Vacuum	3.48 lb	80-0940
Basic	1	Nylon	3.3 lb	80-0935
Intermediate	1	Nylon	3.6 lb	80-0936
Advanced BCD	1	Nylon	3.7 lb	80-0937



8 Pack Station



Individual Kit Station



Small Individual Case:
ZZ-1000

CLEAR CASE
features

Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs) to provide publically accessible medical supplies to control major bleeding



5 Pack Station



Large Multi Pack Case:
ZZ-0822

- Durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Large case features wire shelf to support the contents
- External eyelet allows a secondary safety seal to easily identify if the station has been opened or tampered with
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle



5 Pack Station



8 Pack Station

METAL CASE
features



Metal Case:
ZZ-0921

Stations with Alarms are also available for special order. Call Customer Service for more details.

Semi-Recessed Metal WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or Vacuum-Sealed pouches stored within an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured patients
- Outside Dimensions: H 17.5 in. x W 17.5 in. x D 1 in.
Inside Dimensions: H 14.13 in. x W 14.13 in. x D 6.88 in.
Empty Case Wt: 10 lb (Empty Case Item# ZZ-0921)

8-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	16.8 lb	80-0861
Intermediate	8	Vacuum	17.6 lb	80-0862
Advanced BCD	8	Vacuum	18 lb	80-0863
Advanced	8	Vacuum	18 lb	80-0864
Basic	8	Nylon	18.5 lb	80-0865
Intermediate	8	Nylon	19.3 lb	80-0866
Advanced BCD	8	Nylon	15.7 lb	80-0867
Advanced	8	Nylon	15.7 lb	80-0868

5-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	15.7 lb	80-0577
Intermediate	5	Vacuum	15.7 lb	80-0578
Advanced BCD	5	Vacuum	16.5 lb	80-0580
Advanced	5	Vacuum	16.5 lb	80-0579

All Multi-Pack Wall Stations feature: Highly visible "Bleeding Control Kit" letters for fast ID Instructions for Use posted on the interior back panel



**Low Profile
Metal Cabinet:
ZZ-1199**

Designed to be mounted on a wall with a slim 4 in. clearance and featuring easy pull-to-open design for rapid public access to major bleeding control supplies

8 Pack Nylon Station

Low Profile Metal WALL STATIONS

- Low profile metal cabinet mounts on the wall with a slim 4 in. clearance
- Designed to safely store 8 Individual Bleeding Control Kits in public areas
- Kits are available in red nylon bags or Vacuum-Sealed pouches and stored in Two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon handle
- Features highly visible "Bleeding Control Kit" lettering for fast identification and convenient Instructions for Use on the interior door panel

SPECIFICATIONS

H 15.5 in. x W 22.5 in. x D 4 in.

Station Weight: 16.5 lb

Empty Cabinet Weight: 9.365 lb

Low Pro Cabinet with ID Stickers:
Item# ZZ-1199



8 Pack Vacuum Sealed Station

8-PACK LOW PROFILE METAL WALL STATIONS

LEVEL	#KITS	BAG	ITEM #
Basic	8	Vacuum	80-1083
Intermediate	8	Vacuum	80-1084
Advanced BCD	8	Vacuum	80-1085
Basic	8	Nylon	80-1080
Intermediate	8	Nylon	80-1081
Advanced BCD	8	Nylon	80-1082

Low-Pro CABINET features

- Low profile metal cabinet with 4 in. wall clearance
- Provides convenient slimline storage for medical equipment
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Magnetic closure provides quick easy access with one pull to open - No twist and pull



STB KITS Nylon or Vacuum-sealed

- Small, compact and easy to use
- Kits contain the same Tourniquet and Hemostatic Dressing technology recommended by the U.S. Military & the American College of Surgeons Committee on Trauma to control bleeding
- Available in Vacuum-sealed Pouch or Red Nylon Bag



VAC STB KIT

ITEM

85-1574

SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 8 oz

VACUUM-SEALED STB KIT contains:



SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 15 oz

NYLON STB KIT

ITEM

85-2995

NYLON STB KIT contains:



Accessories

PVC BLEEDING CONTROL PATCH

- Soft flexible PVC rubber patch for labeling your lifesaving medical gear
- Ideal for the Public Access Bleeding Control Kit or any surface with an area for hook-&-loop attachment

SPECIFICATIONS

L 4.375 in. x H 1.75 in. x D 0.125 in.
Weight: 0.275 oz

PVC PATCH

ITEM

ZZ-0796

3-WAY BLEEDING CONTROL SIGN

- Highly visible for locating lifesaving bleeding control kits immediately
- 3 way "V" shape projects from wall for easy identification from multiple directions

3-WAY SIGN

ITEM

ZZ-1117

SPECIFICATIONS

H 7 in. x W 16.5 in. x D 6.75 in.
Weight: 0.35 oz

