

TORK™ (Tactical Operator Response Kit)

REV052720



Point-of-wounding IFAK for the Assault Vest or Body Armor

The Tactical Operator Response Kit™ (TORK™) by North American Rescue® was designed to meet the needs of our customers for a compact individual first aid kit which carries all the necessary equipment to address the top leading causes of potentially preventable death in a tactical environment.

The TORK™ contains essential life-saving equipment that can be used at the point of wounding for self-aid/buddy-aid to address injuries caused by penetrating or blast trauma. It's sleek, compact size allows for easy attachment using MOLLE/PALS style connectors to the assault vest/body armor or NAR's Quick Detach Panel (available for separate purchase). The TORK™ also has a rapid opener with NAR's signature Red-Tip Technology® for fast and easy, one-pull tab opening of the kit.

Features:


- Essential life-saving tools to address injuries caused by penetrating or blast trauma
- Compact individual first aid kit providing critical support for immediate Self-Aid/Buddy-Aid at the point of wounding
- Rugged Nylon Bag with easy MOLLE/PALS style attachment to the assault vest or body armor
- Also attaches to NAR's Quick Detach Panel (sold separately)

Specifications:

- Cased Closed: H 7 in. x W 4.5 in. x D 2.75 in.
- Weight: 1 lb 5 oz



Nylon Colors:



Option	Color	Item#
Advanced w/ChitoGauze	BLK	85-0013**
Advanced w/ChitoGauze	COY	85-0014**
Advanced w/ChitoGauze	DUC	85-1115
Advanced w/ChitoGauze	MTC	85-0017**
Advanced w/ChitoGauze	ODG	85-0055**
Advanced w/ChitoGauze	RGR	85-1116
Advanced w/ChitoGauze	RED	85-1807
Advanced w/CombatGauze	BLK	85-0048
Advanced w/CombatGauze	COY	85-0049
Advanced w/CombatGauze	DUC	85-0050
Advanced w/CombatGauze	MTC	85-0052
Advanced w/CombatGauze	ODG	85-0053
Advanced w/CombatGauze	RGR	85-2302
Advanced w/CombatGauze	RED	85-1806
Basic w/CombatGauze	BLK	80-1097
Basic w/CombatGauze	COY	80-1098
Basic w/CombatGauze	DUC	80-1099
Basic w/CombatGauze	MTC	80-1100
Basic w/CombatGauze	ODG	80-1101
Basic w/CombatGauze	RGR	80-1102
Basic w/CombatGauze	RED	80-1108
Basic	BLK	85-0008
Basic	COY	85-0009
Basic	DUC	85-2300
Basic	MTC	85-0012
Basic	ODG	85-0054
Basic	RGR	85-2301
Basic	RED	85-1805

PRODUCTS
WITH A
MISSION®



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

TORK™ Tactical Operator Response Kits



QTY
1 C-A-T®
Tourniquet




QTY
1 ETD™ 6 in.
Dressing



BASIC
QTY
2 Gauze
S-Rolled



BASIC with Combat Gauze
QTY
1 Gauze
S-Rolled



BASIC with Combat Gauze
QTY
1 Combat
Gauze Z-fold



QTY
1 HyFin®
Vent. Twin Pk



QTY
1 PES
Eye Shield



PAIR
2 Nitrile
Gloves, Lg



QTY
1 Responder
Trauma Shears, Lg

Basic/Basic with Combat Gauze Contents:

- 1 x C-A-T® (Combat Application Tourniquet®) Blk
- 1 x ETD 6 in. Emergency Trauma Dressing
- 2 x S-Rolled Gauze (4.5 in. x 4.1 yd) **BASIC**

-OR-

- 1 x S-Rolled Gauze (4.5 in. x 4.1 yd) **BASIC with Combat Gauze**
- 1 x Combat Gauze Hemostatic
- 1 x HyFin® Vent Chest Seal Twin Pack
- 2 x pair Bear Claw™ Nitrile Glove, Lg.
- 1 x Pair Responder Trauma Shears (7.25 in.)

- 1 x Polycarbonate Eye Shield
- 1 x Nylon Bag (TORK™)



QTY
1 C-A-T®
Tourniquet



QTY
1 ETD™ 6 in.
Dressing



QTY
1 Gauze
S-Rolled



QTY
1 Choice of
Hemostatic



QTY
1 28F Naso
Airway & Lube



QTY
1 HyFin®
Vent. Twin Pk



QTY
1 ARS®
Needle, 14g



QTY
1 PES
Eye Shield



PAIR
2 Nitrile
Gloves, Lg



QTY
1 Responder
Trauma Shears, Lg

Advanced Contents:

- 1 x C-A-T® (Combat Application Tourniquet®) Blk
- 1 x ETD 6 in. Emergency Trauma Dressing
- 1 x S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x Your Choice of Hemostatic:
ChitoGauze PRO -or- Combat Gauze
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 2 x pair Bear Claw™ Nitrile Glove, Lg.
- 1 x Pair Responder Trauma Shears (7.25 in.)
- 1 x Polycarbonate Eye Shield
- 1 x Nylon Bag (TORK™)

**Streamlined and compact...
designed to MOLLE to an assault vest or body armor**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



TORK™ Resupply Kits

Option	Item#
Advanced w/ChitoGauze	80-1038
Advanced w/CombatGauze	80-1039
Basic w/CombatGauze	80-1103
Basic	80-1037

REV052720