

TC SPLINT *Telescoping Collapsible Splint*

REV052720

Meets the Splinting Capability Gap...

The innovative TC Splint comes complete with nylon sleeve, 6 rigid, telescoping inner poles making it easily deployable in the field, but space saving in your pack.

Lightweight and portable, the TC Splint allows for rapid deployment in tactical, field and pre-hospital settings. Size and location adjusting hook and loop straps, padding and sturdy inner poles help with patient pain and discomfort during transport.

Created to meet a capability gap in the splinting market for securing patient injuries.

Features:

- Multi-capable device for splinting fractures and/or breaks
- Capable of securely and rigidly splinting a wide variety of major and minor fractures
- Lightweight and compact for easy field carry
- Water resistant; easy to clean, fold, pack and use
- Simple to use in field or other austere environments
- 6 Easily adjustable interconnected telescoping poles
- 4 easily adjustable wide elastic bands with hook and loop securing
- Rigid and durable high grade 304 stainless steel poles
- Exterior non-absorbent 400 denier Nylon shell
- Internal non-absorbent foam padding for patient comfort
- Patent Pending

Dimensions:

- Packaged: L 8 in. x W 5 in. x D 3 in.
- Open: W 9.5 in. x L 34 in.
- Weight: 1 lb 4 oz



***Telescoping Collapsible Splint
(TC Splint)***

Item #	NSN#
50-0039	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®