

LIFEWARMER™ Quantum™ Blood & Fluid Warming System

- **Lightweight** The Quantum™ has the lightest operational weight of any portable blood & fluid warming system (Total Operational Weight: 22 oz)
- **Powerful** Warms two units of 39°F refrigerated whole blood at 100mL/min+ on a single charge
- **Compact** Fewest components to carry (three)
- **Fast** Warm-up time of only 24 seconds
- **Intuitive** The thermal transfusion/infusion set looks, handles, and sets up/primes like any conventional IV set
- **Safe** Tested safe for blood, blood products & IV solutions

- **Effective** The system is plug-and-play for ease of setup under stressful conditions such as tactical field care (only two connection points, with no plug orientation or set-up order requirements – system status constantly visible at a glance)
- **Complete** Does not require additional standard transfusion or infusion set – all other warmers do
- **Tactical** Only system FDA cleared allowing easy disable of audible alerts to protect operators in combat



Quantum Charger

Quantum Battery

Quantum Controller

TTS-B Thermal Transfusion Set - Blood (Blood, blood products & IV solutions)

TIS Thermal Infusion Set (IV solutions)

Additional Quantum accessories available at www.NARescue.com

PRODUCT WEIGHTS

Total Operational Weight: 22 oz
(Quantum™ Controller, Battery & TTS-B)

Battery: 16 oz

Controller: 3 oz

Battery Charger: 22 oz

Controller Cover: 0.23 oz

Hard Case: 3.22 lb

Thermal Transfusion Set (TTS-B): 2.7 oz

Thermal Infusion Set (TIS): 2.3 oz

*NOTE: Complete System includes 1x Controller, 1x Battery, 1x Charger & 2x TTS-B

OPERATING PARAMETERS

IV FLUIDS

- Input Temp: 68°F
- Output Temp: 100.4°F (38*+/- 2°C)
- Flow-Rate: 2 - 200mL/min (+)
- Warm-up Time: 24 seconds
- Warming Capacity with a Single Battery: 1700 mL

REV052920

BLOOD or BLOOD PRODUCTS

- Input Temp: 39°F
- Output Temp: 100.4°F (38*+/- 2°C)
- Flow-Rate: 2 - 100mL/min (+)
- Warm-up Time: 24 seconds
- Warming Capacity with Single Battery: 2 units (+)



ORDERING INFORMATION

COMPONENTS	ITEM #
Complete System*	35-0006
TTS-B (Pack of 2)	35-0005
TIS (Pack of 2)	35-0004
Controller	35-0003
Battery Charger	35-0002
Battery	35-0001
Controller Cover	35-0010
Hard Case	35-0008



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277