

# NAVAL GUN RESPONSE AID KIT

REV081420



**Nylon Color:** ODG



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

## ***Slim-line lightweight kit stores easily***

Designed to meet the rigorous requirements of the U.S. Navy for first responder emergency care equipment on board U.S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels, the Naval Gun Response Aid Kit from North American Rescue is an ergo-dynamic, low profile casualty response system featuring 2 removable equipment pouches that attach via hook-&-loop to a lightweight, durable outer bag, thus meeting the requirements for the U.S. Navy's Naval Gun Response Aid Bag.

Two internal pouches are identical & provide rapid access & organization of critical life saving equipment needed to render on-scene care at sea.

Constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/ water resistant, self-repairing zippers with silent 550 cord pulls. Both outer bag & 2 inner pouches are securable using antipilferage seals. Slim-line, lightweight kit stores easily on a wall and can be carried over the shoulder or bandolier-style.

### **Features:**

- Two internal pouches are identical & provide rapid access & organization of critical life saving equipment needed to render on-scene care at sea
- Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels
- Constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self-repairing zippers with silent 550 cord pulls.
- Lightweight, durable slimline pouch is securable using anti-pilferage seals & features a clear window to readily view contents
- Available in NAR's 11 point Berry Amendment Nylon Standard

### **Specifications:**

- Kit Closed:  
H 16 in. x W 11 in. x D 5.5 in.
- Weight: 4 lb 11 oz

### ***Naval Gun Response Aid Kit***

Color	Item#	NSN#
ODG	80-0418	6545-01-648-2690

**PRODUCTS**  
WITH A  
**MISSION®**

# NAVAL GUN RESPONSE AID KIT

REV081420

**1 x Naval Gun Response Bag**



## Each Internal Pouch contains:



**2 x Internal Pouches (Naval First Aid Box Response Kit)**

## Kit Contents:

- 1 x Naval Gun Response Bag, CCRK®
  - 2 x Naval First Aid Box Response Kit
- Each Response Kit contains:**
- 1 x Pouch (CCRK® - Naval First Aid Box Response)
  - 1 x Nasopharyngeal Airway 28F
  - 1 x Nasopharyngeal Airway 34F
  - 1 x Surgical Lubrication
  - 1 x CPR Mask, Res-Cue Key
  - 1 x C-A-T® (Combat Application Tourniquet®)
  - 1 x ETD™ 6 in. Emergency Trauma Dressing
  - 1 x ETD™ Abdominal Emergency Trauma Dressing
  - 1 x HyFin® Vent Chest Seal Twin Pack
  - 1 x Fox Eye Shield with Garter
  - 1 x Roll 2 in. Surgical Tape
  - 1 x Trauma Shears (7.25 in.)
  - 2 x SAM® Splint II
  - 1 x Elastic Wrap Bandage, 3 in.
  - 1 x Pair Bear Claw™ Nitrile Gloves, Large
  - 1 x Permanent Marker, Large



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**