

PROTECTIVE ISOLATION GOWN

REV081920



49 in. L x 36 in. W (at the shoulders)
Arms: 20 in. L x 10.5 in. W

A Must Have...

North American Rescue is proud to offer the Protective Isolation Gown. It meets FDA Level 1, 2, and 3 protective gown guidance to help mitigate the spread of communicable diseases. They are individually wrapped for cleanliness and feature front or rear break-away capability to remove.

- Evaluated and passed:

ASTM F2407-06

Liquid Penetration Test ATTCC 42

(American Association of Textile Chemists and Colorists)

- Meets FDA Level 1, 2, and 3 protective gown guidance
- Front or rear break-away capability to remove
- Individually wrapped for cleanliness
- 1.25 mil LDPE

Protective Isolation

Item #

70-0263



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®