

KN95 RESPIRATOR MASKS

REV083120



KN95 Respirator Masks - Box of 50

ITEM#	UNIT
ZZ-1308	Box of 50



More protection than standard disposable 3-Ply masks or your conventional face cover

Offering 5-filtration-layers which meet and exceed KN95 standards for protection, our KN95 Respirator face mask will keep out 95% (0.3 microns) of airborne particles for maximum protection and perform better than a standard disposable 3-Ply mask or conventional face cover. This protective face mask is lightweight and comfortable, with ear loops and a shapeable nose bridge for a flexible and secure seal. Features include 360° protection with three dimensional ventilation space for more comfortable breathing.

Features:

- 5-Layers of protection blocking at least 95% (0.3 microns) of airborne particulates
- More protection than standard disposable 3-Ply masks or your conventional face cover
- Lightweight and foldable, easy to store and carry
- 360° protection with three dimensional ventilation space for more comfortable breathing
- Equipped with metal fastening strip ensuring optimal sealing between the bridge of the nose and the mask
- Sold in a box of 50

Specifications:

- Weight: 0.32 oz (each mask)
- Box Dimensions: L 10 in. x W 4.75 in. x D 4.75 in



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®