

BLACK RESPONDER NITRILE GLOVES

Pre-packaged Pairs

REV090120

Third-party tested compliant with the EN ISO 374-5:2016 Standard for **Virus Protection against bacteria, fungi and viruses** including, but not limited to:

- ▶ **COVID-19** (Coronavirus)
- ▶ **HIV** virus
- ▶ **Ebola** virus
- ▶ **Avian influenza** virus (Bird Flu)
- ▶ **Hantavirus** [causes Hantavirus Pulmonary Syndrome (HPS)]
- ▶ **Marburg** virus
- ▶ **Rabies** virus
- ▶ **Smallpox** virus
- ▶ **Influenza** viruses (A, B and C)
- ▶ **Rotavirus**
- ▶ **HPV** virus

Third party tested for effective protection from Fentanyl as well as bacteria, fungi and viruses including COVID-19 Coronavirus!

The Medical Protective Face Shield is designed as a cost effective general splash guard to add an extra layer of personal protection from debris and chemical hazards. This professional face shield is simply designed to work in conjunction with a face mask to provide another level of PPE. The OSHA compliant clear shield is constructed of 14 mil BPA free plastic with a fully adjustable, one size fits all headband that can accommodate extra room between the face and the shield for goggles. Made in the USA.

Features:

- 25 individually wrapped pairs for convenient rapid deployment
- High performance personal protection compliant with the EN ISO 374-5:2016 Standard for Virus Protection
- Third-party tested for effective personal protection from Fentanyl and Viruses including COVID-19 Coronavirus
- Latex- and Powder- Free with decreased risk of allergies
- Biodegradable and disposable with second-skin fit & feel
- Available in 6 sizes: XS, S, M, L, XL, & XXL

Specifications:

- Gloves: Length: 9.5 in. Thickness (mils): 4.0 Cuff: Rolled
- Packaged (25 pairs): L 4.75 in. x W 6 in. x D 4 in.
- Individual pair: L 2.5 in. x W 1 in. x D .8 in.
- Weight: 16.1 oz

Black Responder PAIRS

SIZE	PKG QTY	ITEM #
XS	Kit 25 pair	ZZ-1285
S	Kit 25 pair	ZZ-1286
M	Kit 25 pair	ZZ-1287
L	Kit 25 pair	ZZ-1288
XL	Kit 25 pair	ZZ-1289
XXL	Kit 25 pair	ZZ-1290

COLOR:
Tactical
Black

**Personal Protection from
Fentanyl and
COVID-19**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®