

HEAVY DUTY TRAUMA X-SHEAR®

REV101320

A Must Have...

The XSHEAR® is a heavy duty trauma shear designed to outperform traditional trauma shears in every way. The hardened 3mm thick stainless steel blades have a comfortable grip with slip resistant inner handles. The patented curved design for gentle approach when cutting near skin is hinged by a heavy duty center bolt that will not loosen over time. The razor sharp blades have a serrated lower blade to help cut through tough material with ease.

Features:

- Hardened stainless steel blades that are twice as thick as most trauma shears, sharpened to a razor edge
- BLACK TITANIUM coating on blades provides extra durability and sleek matte black appearance
- Patented blunt tip and curved design for gentler edge near the skin
- Soft touch inner handles for comfort and slip resistance
- Heavy duty center bolt that won't loosen over time



Heavy Duty Trauma XSHEAR®

| Item # | NSN # | Handle Color | Blade Color |
|---------|-------|--------------|-------------|
| 70-0290 | | Black | Black |
| 70-0291 | | Blue | Silver |
| 70-0292 | | Orange | Black |
| 70-0293 | | Pink | Black |
| 70-0294 | | Red | Black |



Specifications:

- Dimensions: L 7.5 in. x W 3.5 in.
- Weight: 4.0 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

TACTICAL X-SHEARS HOLSTER

REV101320



A Must Have...

The XSHEAR® Heavy Duty Trauma Shears Tactical Holster securely holds the shears in place with a flex lock design which allows the shears to face either direction. The rugged injection molded ABS plastic is very durable and it's MOLLE compatibility as well as an adjustable belt clip offer flexible mounting options and orientation.

Features:

- MOLLE Compatible adjustable belt clip secured holster in place and offers flexible orientation of holster with simple Philips head screw
- Flex lock design securely holds shears in place
- Positive click when shears are secured in place
- Universal design allows shears to be inserted either direction
- Rugged injection molded ABS Plastic

Specifications:

- Dimensions: L 5.0 in. x W 3.75 in. x H 1.75 in.
- Weight: 3.2 oz

X-Shears Holster

| Item # | NSN # | Color |
|---------|-------|---------|
| 70-0295 | | ■ Black |



*Pictured With XSHEARS®
Sold Separately*



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
— WITH —
MISSION®