

# FIELD CORPSMAN KITS *Basic or Deluxe*



## Lightweight and Durable ENT Diagnostic Sets for Field Care...

The Basic Field Corpsman Kit from NAR is a lightweight, durable and portable diagnostic kit provides excellent performance and flexibility in the field.

This water resistant and highly durable kit includes a variety of high intensity lights and magnifiers for day or night examinations.

Choose from the Basic or Deluxe portable diagnostic kit.  
(Deluxe Kit contains all the items from the Basic kit plus a Mini Ophthalmoscope, additional Spare Bulb and extra Batteries).



### Field Corpsman Kits

Option	Item#	NSN#
Basic	20-0001	6515-01-458-6057
Deluxe	20-0002	6515-01-458-6178

*Provides a focused beam of bright light for superior illumination.*

PRODUCTS WITH A  
**MISSION®**



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

NOTE: Deluxe Kit is pictured above.

### Basic & Deluxe Kit Contents:

- General examination light: LED lamp to provide a focused beam of bright white light for superior illumination. Lightweight, water resistant and guaranteed unbreakable (batteries included)
- Red lens: Fits on the exam light for night use
- Fiber optic light bender/nasal head: To effectively direct light into the nasal or oral cavity. Can be used to transilluminate the sinuses
- Spare LED bulb: Extends the use of the exam light in the field
- 2 spare AAA batteries
- Soft case with mesh pouches: Secures the exam lights and accessories. Developed with padded sides and durable zipper. Two additional mesh pouches are positioned on the inside compartment for more parts or supplies
- Color instruction booklet: Contains operating instructions and part numbers

### Deluxe Kit Contents

*includes everything from the Basic Kit above PLUS:*

- Heine Mini 3000 Ophthalmoscope (batteries and spare bulb included)

### Dimensions:

- **BASIC KIT:**  
L 8 in. x W 4.25 in. x D 2.25 in.  
Weight: 6.2 oz
- **DELUXE KIT:**  
L 8 in. x W 4.25 in. x D 2.25 in.  
Weight: 11.6 oz



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION®**





Offers healthcare professionals in the field excellence and flexibility as they care for the eyes and ears of patients.

#### Ear Care:

- **Otoscopic head:** Provides bright illumination and 3X magnification for excellent visualization. Adjustable head allows examiners to adjust for personal comfort
- **Reusable adult specula:** For use with the otoscopic head
- **Disposable adult specula:** For use with the otoscopic head. Otoscope accepts standard disposable specula
- **Lighted stainless curette:** Used to remove wax from the ear or to assist in debris removal from small wounds

#### Eye Care:

- **Cobalt blue lens:** Provides UV illumination to examine the eyes for corneal abrasions or small foreign bodies. Use with or without fluorescent eye dyes
- **High resolution eye loop:** Soft nylon loop fits on the cobalt blue lens. Designed to gently sweep foreign bodies from the surface of the eye. The lighted tip illuminates work area
- **Illuminated magnifier:** Use with both the clear and cobalt blue lens. The magnifier is preferred for examining the soft tissue of the eye. Also useful for locating splinters and examining wounds, rashes and other topical irritations
- **Lighted eye magnet:** For the removal of metal foreign bodies in the eye
- **Fluor-I-Strips:** For use with cobalt blue lens for identification of corneal abrasions or to locate small particles in the eye

*NOTE: Deluxe Kit and contents are pictured.*



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION®**