

OUT-PAK KIT *(Outdoor Personal Aid Kit)*

REV121020


A compact Mini First Aid Kit...

The Outdoor Personal Aid Kit (Out-PAK) from North American Rescue is a compact Mini First Aid Kit designed as a multi-use IFAK for delivering immediate care for injuries encountered in an outdoor adventure. At half the size of our popular T.O.R.K. platform, the Out PAK is currently one of the smallest platforms available in our product offering. Despite its reduced size, this kit comes fully loaded with the essential medical equipment for Basic Life Support with the option of a Bleeding Control Dressing.

The Out-PAK's 500 Denier nylon bag opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T.® tourniquet. Also included are a loop patch on front of the bag with a PVC Medical ID Patch and (3) three MOLLE strips on the sides to attach additional gear such as a TQ/Holder (sold separately) if needed.

Features:

- NAR's most compact, versatile Mini First Aid Kit
- Contains essential medical equipment for immediate care of injuries encountered in an outdoor adventure
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access
- Multiple elastic loops for secure gear organization
- Vertical mount can be set to open left-to-right or right-to-left based on shooter preference
- Horizontal mount on a belt allows opening directly to your C-A-T.® tourniquet
- Treatment pack out option: BLS (Basic Life Support)
- Half the size of our popular T.O.R.K. platform



Dimensions:

- L 8 in. x W 3 in. x H 5.5 in.
- Weight: 1.58 lb



NORTH AMERICAN RESCUE®



www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

OUT-PAK KIT *(Outdoor Personal Aid Kit)*

REV121020

BASIC -or- **BASIC with Bleeding Control Dressing:**



Item #	NSN #	Level	Color
82-0044		Basic	Black
82-0030		Basic	Coyote
82-0045		Basic	Gray
82-0046		Basic	OD Green
82-0047		Basic	RTC
82-0048		Basic	RGR
80-1134		Basic with Bleeding Control Dressing	Black
80-1136		Basic with Bleeding Control Dressing	Coyote
80-1137		Basic with Bleeding Control Dressing	Gray
80-1138		Basic with Bleeding Control Dressing	OD Green
80-1139		Basic with Bleeding Control Dressing	RTC
80-1140		Basic with Bleeding Control Dressing	RGR



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®