

First Light<sup>®</sup> Tomahawk<sup>®</sup> MC Tactical Light (Red/Green) with MOLLE Tactical Retention System (TRS)

Item# 70-0016

Dimensions (with TRS):  
H 3.5 in. x L 3.9 in. x W 2.2 in.

Weight (with TRS): 7.26 oz

In response to the a growing, and to this point, unmet need from tactical medical professionals, NAR is proud to offer the latest innovation from First Light USA<sup>®</sup> - the industry leader in combat-specific lighting tools. The Tomahawk<sup>®</sup> Series of hands-free tactical lights are specially designed to provide maximum versatility and enhanced illumination at the point of injury. The hands-free design adds a unique set of features specifically engineered to assist the Combat Medic - from the superior chest mounting capability and full 180° adjustability of the light source to the

proprietary mix of red and green LEDs that "highlight" blood for immediate wound identification in multi-fluid situations. This not only provides Medics with a quick tool to accurately assess a casualty in low light/night time conditions, it allows them to maintain light discipline as well as a defensible position.

- **Increase Survivability** by moving a potentially targeted light source from the head to the most protected part of the Soldier's body -the chest
- **Reduces Output Signature** by utilizing a non-glaring wavelength (red/green mix); orienting the LEDs down to "shield" them from direct view, and decreasing light movement (when compared to a head-mounted light) -key factors in signature reduction
- **Maximizes Versatility** by enabling the Operator to use the light as either a hands-free or handheld device that can be used in conjunction with small arms
- **Preserves Natural Night Vision** when compared to blue, green, or white light
- **Highlights Blood** with a proprietary mix of red & green LEDs for immediate wound identification in multi-fluid situations

## Tomahawk<sup>®</sup> MC

is the official handheld lighting tool for the **U.S. Army Combat Medic Aid Bag** (MES 246C) and is Team Soldier Certified

NAR offers 2 Tomahawk<sup>®</sup> models which both include the following features:

- Tan MOLLE Tactical Retention System (TRS) for convenient chest mount option
- Rugged construction of aerospace-grade aluminum finished with mil-spec, anodized type III hard coat
- Waterproof
- Powered by 2 CR123A 3-Volt Lithium Cells
- Max White Light: 150 Lumens
- Runtime of 2+ hrs. (High White), 60+ hrs. (Low White), 100+ hrs. (Colored LEDs)
- All output modes have adjustable brightness settings & lock-out mode to prevent accidental illumination
- Double click for color provides on demand tactical lighting from the lock mode
- Covered by lifetime limited warranty

First Light<sup>®</sup> Tomahawk<sup>®</sup> NV/TC3 Tactical Light (Red/Green) plus IR with MOLLE Tactical Retention System (TRS) & UV Task Light

Item# 70-0015

Dimensions (with TRS):  
H 4.1 in. x L 3.9 in. x W 2.2 in.

Weight (with TRS): 8 oz

Blue Light

Red/Green Light

## Tomahawk<sup>®</sup> NV/TC3

also includes integrated IR illumination and a self-powered (three CR1616 coin batteries), twist-on task light that provides diagnostic ultra violet (UV) illumination of ocular injuries in low light situations



**NORTH AMERICAN RESCUE<sup>®</sup>**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277