

# ARMADILLO Medication Storage Case

REV011819

NOTE: All Black (#ZZ-0081) includes inside top, bottom and sides lined with loop (MILSPEC AA55126A Type II)

Clear Lid (#ZZ-0259) has loop lining on the bottom only



NOTE: Medications are not included.

## NAR's most compact solution for keeping your critical medications safe...

Keeping medications secure, safe, and readily available during tactical operations is critical. The case is lined with loop on the top, bottom and sides allowing for custom fitting of medication and maximum use of space.

Adhesive hook discs are provided to hook any medication how and where you want.

### Features:

- Water resistant (3 meters) and crush proof
- Stainless steel (non-rusting) hardware
- Automatic Pressure Equalization Valve
- Inside lined with loop (MILSPEC AA55126A Type II)
- Hook discs (0.5 in.) 100 ea per sheet (MILSPEC AA55126A Type II)
- Additional Hook Discs are available for separate purchase

### Dimensions:

- Exterior: L 7.95 in. x W 5.1 in. x D 2.14 in.
- Interior: L 6.45 in. x W 3.87 in. x D 1.75 in.
- Weight: 10.5 oz

### Armadillo Medication Storage Case

Color	Item#	NSN#
Black	ZZ-0081	8145-01-573-2533
Clear Lid	ZZ-0259	6545-01-594-9399*

\*NATO NSN#

Helpful Hint: During tactical operations significant changes in altitude or temperature can create a vacuum & cause the ARMADILLO to seal closed. To remedy this, insert any flat tipped device into the slot, next to the latch, on the front of the box and twist 45° to equalize the pressure within the case.



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A   
**MISSION®**