

# SAM® SOFT SHELL SPLINT

REV032218



*Custom Moldable for Immobilizing the Hand, Wrist and Forearm*

## *Exceptionally Soft and Comfortable...*

The Soft Shell™ SAM® Splint from North American Rescue is designed to be quickly applied and custom moldable for immobilizing the hand, wrist and forearm.

The exceptionally soft construction material is fabricated to wick moisture away from the skin for the ultimate in patient comfort. Since the splint is radiolucent, it can remain in place during X-rays.

Easily cleaned with cold water and mild soap, this splint air dries for extended performance in the field.

### Features:

- Quick application
- Custom moldable for immobilizing the hand, wrist and forearm
- Exceptionally soft and comfortable
- Fabricated to wick moisture away from the skin
- Cleans effortlessly with cold water and mild soap, air dries
- Fits easily into storage
- Available individually or in a bulk box of 20

### Dimensions:

- Individual Splint: L 15 in. x W 4.25 in. x 0.5 in.
- Weight: 3.2 oz

**SAM Soft Shell Splint**



Size   Unit	Item#	NSN#
15 in.   ea.	50-0008	6515-01-523-5440
15 in.   Box of 20	50-0009	6515-01-507-8074



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**