

WORKPLACE SAFETY | OSHA COMPLIANT

TRAUMA FIRST AID KITS



PUBLIC ACCESS BLEEDING CONTROL KITS | ANSI WORKPLACE FIRST AID TRAUMA KITS



COMMUNITY
PREPAREDNESS

*The first person on scene
IS the first responder...*

Are YOU Prepared?



3^{to}5 MINUTES

The Difference between Life and Death

- ▶ A victim who is bleeding can quickly die BEFORE professional emergency responders can even arrive on the scene
- ▶ Serious bleeding accounts for 35% of prehospital deaths from traumatic injuries
- ▶ Life-threatening bleeding warrants immediate intervention
- ▶ Basic bleeding control provided by bystanders can often be the difference between life and death...

*Kit shown:
Class B Trauma &
First Aid Kit-Hardcase
see page 4*

***FACT:** from 2000-2020,
45% of **ALL Active Shooter Incidents**
occurred in **Places of Commerce**

*Death from
Uncontrolled
BLEEDING
can occur in as
little as **3 to 5 Minutes***

Trauma is a national health problem and after a traumatic injury, serious bleeding is responsible for over 35% of prehospital deaths.** Causes of trauma can range from natural disasters, falls, traffic, industrial, and workplace accidents to intentional acts of violence such as shootings, stabbings and bombings. In fact, according to data collected from 2000 to 2020, the U.S. Department of Justice found that the workplace is the most likely location of an active shooter situation.*

Because these traumatic injuries often result in uncontrolled bleeding,

immediate intervention – even before emergency medical services or other first responders can arrive – is critical to avoid potential death.

North American Rescue® is leading the way to equip bystanders and initial first responders with the essential medical gear needed to take immediate action in life threatening emergencies. These lifesaving modular aid components are designed for individuals or groups in the workplace, public buildings, and populated public areas to administer immediate bleeding control at the point of wounding.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

*Active Shooter Incidents 20-Year Review, 2000-2020, Federal Bureau of Investigation, U.S. Department of Justice, Washington, D.C.

**Chambers JA, Seastedt K, Krell R, Caterson E, Levy M, Turner N. "Stop the Bleed": A U.S. Military Installation's Model for Implementation of a Rapid Hemorrhage Control Program. Mil Med. 2019 Mar 01;184(3-4):67-71. [PubMed] [Reference list] <https://academic.oup.com/milmed/article/184/3-4/67/5062295>

Workplace First Aid Kits

Class A **TFAK** (TRAUMA & FIRST AID KIT)

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2021 Standard for Class A kits**
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T® Tourniquet - the safest and most effective Tourniquet on the market!
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

TRAUMA & FIRST AID KIT - CLASS A

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0947
Soft Case	Bleeding Control Dressing	80-1045
Soft Case	CELOX™ First Aid Gauze	80-1477
Hard Case	Compressed Gauze	80-1033
Hard Case	Bleeding Control Dressing	80-1046
Hard Case	CELOX™ First Aid Gauze	80-1481

SPECIFICATIONS

Soft Case::

L 14.75 in. x H 8 in. x D 4 in.
Weight: 2.94 lb

Hard Case:

L 13 in. x H 9.5 in. x D 5 in.
Weight: 6 lb

Type IV ANSI Container only option (#80-1113) also available at NARescue.com

TRAUMA & FIRST AID KIT CLASS A:

 QTY 1 C-A-T® Tourniquet	 QTY 1 HyFin® Vent Compact Twin	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD® 4 in. Responder, Flat	 QTY 1 4" Elastic Bandage	 QTY 1 2" Elastic Bandage	 QTY 1 Triangular Bandage
 QTY 1 SAM® Splint II	 QTY 2 4" x 4" Gauze Pad (2pk)	 QTY 2 Gauze Pad, 5" x 9"	 QTY 16 1" x 3" Bandage	 QTY 1 BurnTec 4"x4" Dressing	 QTY 1 CPR Microshield Mask	 QTY 10 Antiseptic Towelette
 QTY 10 Antibiotic Ointment	 QTY 1 Surgical Tape, 1 in.	 QTY 1 6"x9" Instant Cold Pack	 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 Sterile Oval Eye Pad	 PAIR 2 Responder Nitrile Gloves, Lg	 QTY 1 Responder Shears, Lg
 QTY 10 Hand Sanitizer	 QTY 10 Safetec Burn Gel	 QTY 1 Survival Blanket	 QTY 1 Permanent Marker, Large	 QTY 1 First Aid Pocket Guide	 QTY 1 Instructions for Use Bleeding Control	 QTY 1 Kit Contents Inventory Card

PLUS+
one of these:



TFAK RESUPPLY - CLASS A

DRESSING	ITEM #
Compressed Gauze	80-1153
Bleeding Control Dressing	80-1155
CELOX™ First Aid Gauze	80-1479

Container not included with ReSupply Kits



Workplace First Aid Kits

Class B **TFAK** (TRAUMA & FIRST AID KIT)

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2021 Standard for Class B kits**
- Same supplies as the Class A Kit but with larger quantities to accommodate larger workplaces
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

SPECIFICATIONS

Soft Case:

L 14.75 in. x H 8 in. x D 7.5 in.
Weight: 4.87 lb

Hard Case:

L 13.5 in. x H 11 in. x D 6 in.
Weight: 9 lb

TRAUMA & FIRST AID KIT - CLASS B

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0948
Soft Case	Bleeding Control Dressing	80-1047
Soft Case	CELOX™ First Aid Gauze	80-1478
Hard Case	Compressed Gauze	80-1034
Hard Case	Bleeding Control Dressing	80-1048
Hard Case	CELOX™ First Aid Gauze	80-1482

Type IV ANSI Container only option (#80-1114) also available at NARescue.com

TRAUMA & FIRST AID KIT CLASS B:



 QTY 1 C-A-T™ Tourniquet	 QTY 1 HyFin® Vent Compact Twin	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD™ 4 in. Responder, Flat	 QTY 1 4" Elastic Bandage	 QTY 2 2" Elastic Bandage	 QTY 2 Triangular Bandage
 QTY 1 SAM® Splint II	 QTY 4 4" x 4" Gauze Pad (2pk)	 QTY 4 Gauze Pad, 5" x 9"	 QTY 50 1" x 3" Bandage	 QTY 2 BurnTec 4" x 4" Dressing	 QTY 1 CPR Microshield Mask	 QTY 50 Antiseptic Towelette
 QTY 25 Antibiotic Ointment	 QTY 2 Surgical Tape, 1 in.	 QTY 2 6" x 9" Instant Cold Pack	 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 Sterile Oval Eye Pad	 PAIR 4 Responder Nitrile Gloves, Lg	 QTY 1 Responder Shears, Lg
 QTY 20 Hand Sanitizer	 QTY 25 Safetec Burn Gel	 QTY 1 Survival Blanket	 QTY 1 Permanent Marker, Large	 QTY 1 First Aid Pocket Guide	 QTY 1 Instructions, Bleeding Control	 QTY 1 Kit Contents Inventory Card
	 QTY 1 Kit Contents Inventory Card					

TFAK RESUPPLY- CLASS B

DRESSING	ITEM #
Compressed Gauze	80-1154
Bleeding Control Dressing	80-1156
CELOX™ First Aid Gauze	80-1480

Container not included with ReSupply Kits

PLUS+ one of these:

 QTY 1 Responder Compressed Gauze	- OR -	 QTY 1 Bleeding Control Dressing	- OR -	 QTY 1 5" CELOX Hemostatic Dressing
--	--------	--	--------	---

LOGGING FIRST AID & TRAUMA KIT *Hard Case*

- Rugged case designed for use in the harshest of conditions
- Modular system featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels (**Meets Type IV ANSI Container**)
- Fast open dual-latch closure system allows for rapid access while keeping contents secure
- Able to treat the most common job site injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market
- Designed with mobility in mind, the kit can be easily carried to the patient's side

SPECIFICATIONS

L 13.5 in. x H 11 in. x D 6 in.
Weight: 5.66 lb

LOGGING FIRST AID KIT & TRAUMA KIT - HARD CASE

OPTION	ITEM #
Basic	80-1190
with Bleeding Control Dressing	80-1191



LOGGING FIRST AID & TRAUMA KIT HARD CASE:

<p>QTY 1 C-A-T® Tourniquet</p>	<p>80-1190 QTY 2 Responder Compressed Gauze</p>	-OR-	<p>80-1191 QTY 1 Responder Compressed Gauze</p>	<p>QTY 1 Bleeding Control Dressing</p>
<p>QTY 1 ETD® 4 in. Responder</p>	<p>QTY 1 MicroShield CPR Face Mask</p>	<p>QTY 2 Elastic Wrap 2 in. Bandage</p>	<p>QTY 2 Gauze Pads 4"x4" (2pk)</p>	<p>QTY 2 Gauze Pads 8"x10" (2pk)</p>
<p>QTY 16 Bandage 1 in. x 3 in.</p>	<p>QTY 1 BurnTec 4"x4" Dressing</p>	<p>QTY 1 Cold Compress</p>	<p>QTY 1 Surgical Tape, 2" Mini</p>	<p>QTY 2 Triangular Bandage</p>
<p>QTY 1 SAM® Splint II</p>	<p>QTY 1 Eyewash Solution 1 oz</p>	<p>QTY 2 PES Eye Shield</p>	<p>QTY 10 Antiseptic Towelette</p>	<p>QTY 1 Forceps, Splinter, 4.5 in.</p>
<p>QTY 10 Sting & Bite Towelette</p>	<p>QTY 1 SAM® Splint II</p>	<p>QTY 1 Eyewash Solution 1 oz</p>	<p>QTY 2 PES Eye Shield</p>	<p>QTY 10 Antiseptic Towelette</p>
<p>QTY 1 Forceps, Splinter, 4.5 in.</p>	<p>QTY 1 SAM® Splint II</p>	<p>QTY 1 Eyewash Solution 1 oz</p>	<p>QTY 2 PES Eye Shield</p>	<p>QTY 10 Antiseptic Towelette</p>
<p>QTY 10 Sting & Bite Towelette</p>	<p>QTY 1 SAM® Splint II</p>	<p>QTY 1 Eyewash Solution 1 oz</p>	<p>QTY 2 PES Eye Shield</p>	<p>QTY 10 Antiseptic Towelette</p>
<p>QTY 1 Inventory Contents Card</p>	<p>QTY 1 Instructions for Use</p>	<p>QTY 1 Just in Time Instructions</p>	<p>QTY 1 Card, Inspection</p>	

LOGGING FIRST AID & TRAUMA KIT SOFT CASE:



LOGGING FIRST AID & TRAUMA KIT Soft Case

- Designed as a portable First Aid Kit for the workplace that **exceeds OSHA Standard Number: 1910.266 App A - First-aid Kits** (Mandatory) for the logging industry
- Compact size allows for storage in vehicles and machinery
- Rugged nylon constructed bag keeps contents secure
- Durable 550 cord pulls
- Able to treat the most common job site injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market
- Designed with mobility in mind, the kit can be easily carried to the patient's side

LOGGING FIRST AID KIT & TRAUMA KIT - SOFT CASE

OPTION	ITEM #
Basic	80-1188
with Bleeding Control Dressing	80-1189

SPECIFICATIONS

L 14.75 in. x H 8 in. x D 4 in.
Weight: 3.03 lb



Individual Kits | Public Access Bleeding Control

Public Access INDIVIDUAL BLEEDING CONTROL KITS

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control ID Patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.5 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

*Essential equipment needed
for persons to take action
as immediate responders in
helping to stop life threatening
bleeding in Public Areas...*



PVC BLEEDING CONTROL PATCH

- Soft flexible PVC rubber patch for labeling your lifesaving medical gear
- Ideal for the Public Access Bleeding Control Kit or any surface with an area for hook-&-loop attachment

SPECS

L 4.375 in. x H 1.75 in.
x D 0.125 in.
Weight: 0.275 oz

PATCH

ITEM

ZZ-0796



Every **INDIVIDUAL BLEEDING CONTROL KIT** contains:



PLUS these Individual Components based upon YOUR configuration choice:




PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KITS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	80-0452
Intermediate	Red Nylon	16.5 oz	80-0453
Advanced BCD	Red Nylon	17.4 oz	80-0523
Advanced	Red Nylon	17.4 oz	80-0454
Basic	Vacuum Sealed	11.7 oz	80-0465
Intermediate	Vacuum Sealed	13.25 oz	80-0466
Advanced BCD	Vacuum Sealed	14 oz	80-0522
Advanced	Vacuum Sealed	14 oz	80-0467

Public Access Bleeding Control | Multi-Packs

**All Intermediate/
Advanced Public
Access Bleeding
Control Kits and
Stations are Texas
House Bill 496
Compliant**



PUBLIC ACCESS BLEEDING CONTROL 8-PACK NYLON



PUBLIC ACCESS BLEEDING CONTROL 8-PACK VACUUM



Public Access

BLEEDING CONTROL 8-PACKS

Nylon | Vacuum-Sealed

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon carrying case
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced & Advanced BCD

SPECIFICATIONS

8-Pack: H 14 in. x W 11 in. x D 8 in.

Each Individual Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in.

PUBLIC ACCESS BLEEDING CONTROL 8-PACKS

LEVEL	BAG	WGT	ITEM #
Basic	Red Nylon	8.5 lb	80-0460
Intermediate	Red Nylon	9.3 lb	80-0461
Advanced BCD	Red Nylon	9.7 lb	80-0525
Advanced	Red Nylon	9.7 lb	80-0462
Basic	Vacuum Sealed	6.8 lb	80-0468
Intermediate	Vacuum Sealed	7.6 lb	80-0469
Advanced BCD	Vacuum Sealed	8 lb	80-0524
Advanced	Vacuum Sealed	8 lb	80-0470

Public Access BLEEDING CONTROL TWIN PACKS *Nylon*

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Contents of Two (2) Public Access Individual Bleeding Control Kits in one convenient nylon bag
- Packaged in a compact, ruggedized nylon bag designed to fit next to an AED in the cabinet
- Available in Basic, Intermediate, and Advanced BCD configurations

PUBLIC ACCESS BLEEDING CONTROL TWIN PACK

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	1.66 lb	80-0951
Intermediate	Red Nylon	1.8 lb	80-0952
Advanced BCD	Red Nylon	1.8 lb	80-0953

SPECIFICATIONS

Nylon Bag:
H 10 in. x W 9 in.
x D 3 in.



Every **BLEEDING CONTROL TWIN PACK** contains:



PLUS + these Individual Components based upon YOUR configuration choice:



Slimline profile ideal for fitting next to an AED in the cabinet!

Trauma & First Aid ESSENTIAL BLEEDING CONTROL KIT

Every **ESSENTIAL BLEEDING CONTROL KIT** contains:



PLUS+ one of these:



- Compact and Convenient Kit
- Highly mobile and easily deployed
- Hypothermia prevention included
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market



SPECIFICATIONS

L 12 in. x H 7.5 in.
Weight: 0.68 lbs

TRAUMA & FIRST AID ESSENTIAL BLEEDING CONTROL KIT

OPTION	ITEM #
Compressed Gauze	80-1474
Bleeding Control Dressing	80-1475
CELOX™ First Aid Gauze	80-1476



Public Access BLEEDING CONTROL 6-PACKS *Vacuum-Sealed*

- Designed to provide bystanders and initial first responders with quick and easy access to essential medical equipment for stopping life-threatening bleeding
- Includes an easy-to-open zippered tote containing six (6) Vacuum Sealed Individual Bleeding Control Kits packaged in ruggedized vacuum sealed pouches and one (1) NAR® QuikLitter™ LITE for moving injured casualties
- Available in Basic, Intermediate, Advanced care configurations


SPECIFICATIONS

6-Pack:

H 7.5 in. x W 13 in. x D 8 in.

Each Individual
Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in

PUBLIC ACCESS BLEEDING CONTROL 6-PACK VACUUM



PUBLIC ACCESS BLEEDING CONTROL 6-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Vacuum Sealed	5.7 lb	80-0573
Intermediate	Vacuum Sealed	5.75 lb	80-0574
Adv.	Vac Vacuum Sealed	6.5 lb	80-0575
Adv. BCD	Vacuum Sealed	6.5 lb	80-0576
Adv. w/ CELOX	Vacuum Sealed	6.5 lb	80-1571

All Int./Adv. Public Access Bleeding Control Kits and Stations are
Texas House Bill 496 Compliant

Mini First Aid Kit M-FAK®



- Personal First Aid Kit with the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

M-FAK


SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.

Weight: 13 oz

LEVEL	ITEM #
Basic	80-0890
with BCD	80-0891

M-FAK Contents:



Clear Polycarbonate WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 6 & 8 Pack Stations also include an easy to-open nylon carrier and one (1) Responder QuikLitter™ (6 Packs contain QuikLitter™ Lite) for moving injured casualties*

SPECIFICATIONS

6 & 8 Pack Stations: H 18 in. x W 15.5 in. X D 10 in.
Empty Case Wt: 15 lb (includes rack)
Large Case Item# ZZ-0822

Individual Stations: H 10 in. x W 10 in. X D 5.25 in.
Empty Case Wt: 2.39 lb
Small Case Item# ZZ-1000

8-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	80-0471
Intermediate	8	Vacuum	22.6 lb	80-0472
Adv.	8	Vacuum	23 lb	80-0473
Adv. BCD	8	Vacuum	23 lb	80-0526
Adv. w/ CELOX	8	Vacuum	23 lb	80-1575
Basic	8	Nylon	23.5 lb	80-0474
Intermediate	8	Nylon	24.3 lb	80-0475
Adv.	8	Nylon	27.4 lb	80-0476
Adv. BCD	8	Nylon	27.4 lb	80-0527
Adv. w/ CELOX	8	Nylon	27.4 lb	80-1576

6-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	6	Vacuum	20.7 lb	80-0581
Intermediate	6	Vacuum	20.7 lb	80-0582
Adv.	6	Vacuum	21.5 lb	80-0583
Adv. BCD	6	Vacuum	21.5 lb	80-0584
Adv. w/ CELOX	6	Vacuum	21.5 lb	80-1574

IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	3.3 lb	80-0938
Intermediate	1	Vacuum	3.42 lb	80-0939
Adv. BCD	1	Vacuum	3.48 lb	80-0940
Adv. w/ CELOX	1	Vacuum	3.48 lb	80-1580
Basic	1	Nylon	3.3 lb	80-0935
Intermediate	1	Nylon	3.6 lb	80-0936
Adv. BCD	1	Nylon	3.7 lb	80-0937
Adv. w/ CELOX	1	Nylon	3.7 lb	80-1581

Alarm (#ZZ-1075) is also available for the Clear Polycarbonate Wall Cabinet

All Intermediate/Advanced Public Access Bleeding Control Kits and Stations are Texas House Bill 496 Compliant



8 Pack Station



Individual Kit Station



Small Individual Clear Case: ZZ-1000



6 Pack Station



Large Multi-Pack Clear Case: ZZ-0822

Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs) to provide publically accessible medical supplies to control major bleeding

CLEAR CASE
features

- Durable, clear polycarbonate cabinet provides both physical and environmental protection for medical equipment
- Large case features wire shelf to support the contents
- External eyelet allows a secondary, tamper-evident safety seal
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle



6 Pack Station



8 Pack Station

METAL CASE
features

- Mounts permanently inside the wall with a low profile 1 in. trim
- Aluminum cabinet with a white finish and clear plexiglass window
- Stainless steel continuous hinge & chrome plated handle designed to open with one quick pull

Semi-Recessed Metal Case: ZZ-0921



Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.

Semi-Recessed Metal WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or Vacuum-Sealed pouches stored within an easy to-open nylon carrier and one (1) Responder QuikLitter™ (6 Packs contain QuikLitter™ Lite) for moving injured patients*
- Outside Dimensions: H 17.5 in. x W 17.5 in. x D 1 in.
Inside Dimensions: H 14.13 in. x W 14.13 in. x D 6.88 in.
Empty Case Wt: 10 lb (Empty Case Item# ZZ-0921)

8-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	16.8 lb	80-0861
Intermediate	8	Vacuum	17.6 lb	80-0862
Adv.	8	Vacuum	18 lb	80-0864
Adv. BCD	8	Vacuum	18 lb	80-0863
Adv. w. CELOX	8	Vacuum	18 lb	80-1578
Basic	8	Nylon	18.5 lb	80-0865
Intermediate	8	Nylon	19.3 lb	80-0866
Adv.	8	Nylon	15.7 lb	80-0868
Adv. BCD	8	Nylon	15.7 lb	80-0867
Adv. w/ CELOX	8	Nylon	15.7 lb	80-1579

6-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	6	Vacuum	15.7 lb	80-0577
Intermediate	6	Vacuum	15.7 lb	80-0578
Adv.	6	Vacuum	16.5 lb	80-0579
Adv. BCD	6	Vacuum	16.5 lb	80-0580
Adv. w/ CELOX	6	Vacuum	16.5 lb	80-0577

*All Multi-Pack Wall Stations feature:
Highly visible "Bleeding Control Kit" letters for fast ID
Instructions for Use posted on the interior back panel*



**Low Profile
Metal Cabinet:
ZZ-1199**

Designed to be mounted on a wall with a slim 4 in. clearance and featuring easy pull-to-open design for rapid public access to major bleeding control supplies

Low Profile Metal WALL STATIONS

8 Pack Nylon Station

- Low profile metal cabinet mounts on the wall with a slim 4 in. clearance
- Designed to safely store 8 Individual Bleeding Control Kits in public areas
- Kits are available in red nylon bags or Vacuum-Sealed pouches and stored in Two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon handle
- Features highly visible "Bleeding Control Kit" lettering for fast identification and convenient Instructions for Use on the interior door panel

SPECIFICATIONS

H 15.5 in. x W 22.5 in. X D 4 in.

Empty Cabinet Weight: 9.365 lb

Low Pro Cabinet with ID Stickers:

Item# ZZ-1199

All Intermediate/Advanced Public Access Bleeding Control Kits and Stations are Texas House Bill 496 Compliant



8 Pack Vacuum Sealed Station

8-PACK LOW PROFILE METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	17.49 lb	80-1083
Intermediate	8	Vacuum	18.19 lb	80-1084
Advanced BCD	8	Vacuum	18.69 lb	80-1085
Basic	8	Nylon	17.49 lb	80-1080
Intermediate	8	Nylon	18.19 lb	80-1081
Advanced BCD	8	Nylon	18.69 lb	80-1082

Low-Pro CABINET features

- Low profile metal cabinet with 4 in. wall clearance
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Provides convenient slimline storage for medical equipment
- Magnetic closure provides quick easy access with one pull to open - No twist and pull

Alarm (#ZZ-1075) is also available for the Low Profile Metal Cabinet

Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.

- Small, compact and easy to use
- Kits contain the same Tourniquet and Hemostatic Dressing technology recommended by the U.S. Military & the American College of Surgeons Committee on Trauma to control bleeding
- Available in Vacuum-sealed Pouch or Red Nylon Bag

STB KITS

Nylon or Vacuum-Sealed



VAC STB KIT

ITEM

85-1574

SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 8 oz.

VACUUM-SEALED STB KIT contains:



SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.


Weight: 15 oz

NYLON STB KIT

ITEM

85-2995

NYLON STB KIT contains:



Cabinet Security ALARMS

- Alarm compatible with Clear Polycarbonate and Low Profile Metal Cabinets
- Loud 120 decibel alarm sounds when magnet separates from base more than 0.5 in.
- Easy and convenient to install
- Requires no wiring
- Mounting hardware and instructions included
- Designed for indoor use only



SPECIFICATIONS

W 163 in. x H 2.41 in. x D 0.72 in.

Weight: 6.4 oz each

CABINET SECURITY ALARMS

ITEM#	UNIT
ZZ-1075	Pack of 2

Accessories

3-WAY SIGN

Bleeding Control



- Highly visible for locating lifesaving bleeding control kits immediately
- 3 way "V" shape projects from wall for easy identification from multiple directions

SPECIFICATIONS

H 7 in. x W 16.5 in. x D 6.75 in.

Weight: 0.35 oz

3-WAY SIGN

ITEM

ZZ-1117

Personal First Aid Kits

I-ROK® Individual Rugged Outdoor Kit

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding
- Pelican R40 ruck case is crushproof, dustproof and waterproof up to 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- Available in 4 hardcase colors and 2 sizes with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions



Every I-ROK contains the following items:



SPECIFICATIONS

Int: L 7.6 in. x W 4.7 in. x D 1.9 in.
Ext: L 9.8 in. x W 6.1 in. x D 2.9 in.
Weight: 2.48 lb.




I-ROK®

COLOR	LEVEL	ITEM #
Black	Basic	82-0031
OD Green	Basic	82-0033
Orange	Basic	82-0034
Black	w/BCD	82-0035
OD Green	w/BCD	82-0037
Orange	w/BCD	82-0038
Black	w/CELOX	82-0066
OD Green	w/CELOX	82-0068
Orange	w/CELOX	82-0069

I-ROK® XL

- Extra supplies to accommodate multiple injuries
- Large, Pelican R60 Ruck Case w/ removable tray

I-ROK XL includes these addl. items:



SPECIFICATIONS

Int: L 8.9 in. x W 5.4 in. x D 3 in.
Ext: L 11.1 in. x W 6.9 in. x D 3.9 in.
Weight: 3 lb.



I-ROK® XL

COLOR	LEVEL	ITEM #
Black	Basic	82-0049
OD Green	Basic	82-0051
Orange	Basic	82-0052
Black	w/BCD	82-0053
OD Green	w/BCD	82-0055
Orange	w/BCD	82-0056
Black	w/CELOX	82-0070
OD Green	w/CELOX	82-0072
Orange	w/CELOX	82-0073

Public Access **INDIVIDUAL BLEEDING CONTROL TRAINER KITS**

- The appropriate equipment to teach individuals the necessary life-saving skills of bleeding control
- Corresponds with NAR®'s single Public Access Bleeding Control Kits in Basic, Intermediate, or Advanced
- Packaged in a compact, blue nylon bag with PVC Bleeding Control Patch

Every **INDIVIDUAL TRAINER KIT** contains:



SPECIFICATIONS
H 7 in. x W 5 in. x D 2.5 in.

PLUS these Individual Components based upon **YOUR** configuration choice:



BLEEDING CONTROL INDIVIDUAL TRAINER KITS

LEVEL	WGT	ITEM #
Basic	15 oz	80-0901
Intermediate	16.5 oz	80-0902
Advanced	16.5 oz	80-0903



STB TRAINER contains:



STB TRAINER

- Designed to provide the essential training equipment to teach how to control bleeding
- Contains training items corresponding with the necessary tools recommended by the American College of Surgeons Committee on Trauma to equip bystanders to help save lives in the instance of life threatening bleeding

SPECIFICATIONS
H 7 in. x W 5 in. x D 2.5 in.
Weight: 8 oz

STB TRAINER

ITEM #
85-1880

Multi-User Skills Training Kits | Public Access Bleeding Control

BASIC BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



INTERMEDIATE BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



ADVANCED BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



BLEEDING CONTROL TRAINING SKILLS KITS

- Provides the appropriate equipment for up to 20 students to learn the necessary lifesaving skills of bleeding control
- Packaged in a self-contained, wheeled, rigid case with internal dividers for easy organization and inventory
- Skill set training addresses extremity bleeding, wound packing, pressure bandages, chest injuries and patient transport
- Corresponds with NAR®'s Public Access Bleeding Control Kit Levels: Basic, Intermediate & Advanced
- Available with or without Wound Mgt. Simulator

SPECIFICATIONS

Hard Storage Case:
W 32 in.x H 21 in. x D 12.5 in.

BLEEDING CONTROL TRAINING SKILLS KITS

LEVEL	WGT	ITEM #
Basic	40 lb	80-0881
Intermediate	42 lb	80-0879
Advanced	42 lb	80-0877
Advanced w/CELOX	42 lb	80-1583
Basic Wound Mgt. Simulator	44 lb	80-0880
Intermediate Wound Mgt. Simulator	48 lb	80-0878
Adv. Wound Mgt. Simulator	48 lb	80-0876
Adv. w/CELOX Wound Mgt. Simulator	48 lb	80-1584



Bleeding Control Training Aids

COMBAT GAUZE® Mouflage Trainer



- Thoroughly instruct students on the proper application of Z-Fold Combat Gauze®

- Completely inert (no hemostatic properties) - training purposes only

SPECIFICATIONS

Packaged: H 6 in. x W 4.5 in. x D 1.25 in.
Gauze: L 4 yd x W 3 in.
Weight: 0.75 oz.


CG TRAINER

ITEM

30-0063

C-A-T® TRAINER

Combat Application Tourniquet



- Fully functioning C-A-T® tourniquet colored blue to identify it as a training device
- Enables instructors and trainees hands-on training & tourniquet application practice
- Providing a distinct visual means to distinguish training devices from personal/operational life-saving devices

C-A-T® TRAINER

ITEM

30-0033

SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in.

Weight: 2.7 oz

HYFIN® TRAINER

Twin Pack



- Unique adhesive gel formulation does not damage mannequins or leave excessive residue behind
- Packaged with two training chest seals in one unique perforated package allowing students to open only one seal at a time

SPECIFICATIONS

Packaged:
L 7.5 in. x W 4.5 in. x D 0.25 in.
Chest Seal: H 6 in. x W 6 in. x D 0.1 in.
Weight: 2.5 oz
U.S. Patent: 7,504,549

HYFIN® TRAINER

ITEM

10-0032

WOUND MANAGEMENT SIMULATOR

- Rugged part-task trainer for learning essential life-saving skills of tourniquet application, wound packing and applying pressure bandages.
- Multiple wounding patterns to enhance training objectives
- Simulates severe laceration, penetrating trauma (entrance and exit wounds), along with the ability to apply tourniquets
- Proprietary silicon that is rugged and durable and will not break down due to environmental exposure
- Can be used in the classroom or in field training scenarios



Portable training aid to help learn wound management techniques to control bleeding

WOUND MANAGEMENT SIMULATOR

SKIN TONE

ITEM

Light

93-0069

Dark

93-0410

SPECIFICATIONS

L 15.38 in. x Diameter 6 in.
Weight: 10.38 lb

Tourniquets | Hemostatic | Gauze

C-A-T® Combat Application Tourniquet Gen 7

The Fastest, Safest, Most Effective Prehospital Field Tourniquet

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC

online video available

scan to view>



C-A-T TOURNIQUET

COLOR	STYLE	ITEM #
Black	Tactical	30-0001
Orange	Rescue	30-0023
Blue	Trainer	30-0033

SPECIFICATIONS

Packaged:
H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 3.4 oz.

NAR® WOUND PACKING GAUZE

NAR® WOUND PACKING GAUZE

ITEM

30-0054

- Ultra compressed, soft, white 1-ply gauze
- Includes embedded X-Ray detectable line
- Z-fold packaging for controlled application as a complete unit or in a continuous linear feed
- Used as a packing gauze for wounds, as a backing gauze for hemostatic agents, with a compression bandage to control hemorrhage, or for minor wound bandaging

SPECIFICATIONS

Gauze: W 4.5 in. x L 4.1 yd
Packaged: L 3 in. x W 2 in. x D 1 in.
Weight: 12 oz

CELOX™ FIRST AID Z-FOLD

- Z-Fold Configuration
- Fast-acting technology
- Advanced hemostatic agent combined with the simplicity of gauze
- Easy to apply on superficial wounds, minor cuts, and abrasions

SPECIFICATIONS

Gauze: W 3 in. x L 5 ft.
Packaged: L 5.5 in. x W 4.5 in.
Weight: 0.1 lbs

CELOX™ First Aid

ITEM

30-0300

QUIKCLOT® BLEEDING CONTROL DRESSING

- Advanced hemostatic agent combined with the simplicity of gauze
- Utilizes the same inorganic component as QuikClot Combat Gauze® used by the US Military

SPECIFICATIONS

Z-fold Gauze L 4 yd x W 3 in.
Pkg: H 5 in. x W 5.5 in. x D 0.625 in.
Weight: 0.8 oz
Roll Gauze L 4 ft x W 3 in.
Pkg: H 5 in. x W 3.5 in. x D 0.5 in.
Weight: 0.5 oz

QUIKCLOT BCD

STYLE	ITEM #
Z-Fold	30-0161
Roll	30-0164

NAR® Responder COMPRESSED GAUZE

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption

SPECIFICATIONS

Gauze: W 4.5 in. x L 4.1 yd
Packaged:
L 3 in. x W 2.5 in. x D 1 in.
Weight: 1.2 oz

NAR® RESPONDER COMPRESSED GAUZE

ITEM

30-0182



NAR® S-ROLLED GAUZE

- Unique inner dispenser, allows controlled linear feed application
- Premium 100% cotton

SPECIFICATIONS

Gauze: W 4.5 in. x L 4.1 yd
Packaged:
L 3 in. x W 2.5 in. x D 1 in.
Weight: 1.2 oz

NAR® S-ROLLED GAUZE

ITEM

30-0003



RESPONDER ETD® Emergency Trauma Dressing

- Sterile pressure bandages that deliver extremely effective pressure and excellent absorption
- Suitable for applying immediate direct pressure
- Equipped with a non-adherent pad
- Easy-to-use hook & loop securing device with no confusing pressure bars or hooks
- Simple application with the Quick grip roll control system which prevents the dressing from unrolling if dropped.
- Rugged vacuum-sealed, low-cube packaging with easy-open tear notches



RESPONDER EMERGENCY TRAUMA DRESSINGS

STYLE	BANDAGE	PAD	WGT	PACKAGE	ITEM #
4 in. Flat	W 4 in. x L 60 in.	4x6	2 oz	L 4.25 in. x W 2.75 in. x D 0.88 in.	30-0171
4 in. Flat (Mini)	W 4 in. x L 42 in.	4x6	2.4 oz	L 5 in. x W 3.5 in. x D 0.75 in.	30-0179
6 in. Rolled	W 6 in. x L 60 in.	6x6	2.9 oz	L 6 in. x W 2.5 in. x D 2.1 in.	30-0168

RESPONDER Nitrile Gloves



- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves

SPECIFICATIONS

Bag (25 pair): L 4.75 in. x W 6 in. x D 4 in.
Weight: 1 lb 1 oz

RESPONDER NITRILE GLOVES

SIZE	PKG QTY	ITEM #
S	Bag 25 pair	ZZ-1059
M	Bag 25 pair	ZZ-1060
L	Bag 25 pair	ZZ-1043
XL	Bag 25 pair	ZZ-1061

*also available in black at NARescue.com

HYFIN® VENT COMPACT Twin Pack

- Two Chest Seals packaged approximately 25% smaller than our standard HyFin®
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary


SPECIFICATIONS

Device (excludes tab):
4.75 in. x 4.75 in.
Wgt: 1.55 oz
Packaged:
H 3.88 in x W 6.6 in x D 0.13 in.
U.S. Patent: 7,504,549

HYFIN VENT COMPACT

ITEM

10-0042



Accessories | Hypothermia Prevention | Evacuation

NAR® SURVIVAL BLANKET

Aluminized, reflective PET to help retain up to 90% of your body heat

SURVIVAL BLANKET

ITEM #	NSN #
ZZ-0038	4240-01-697-8869

SPECIFICATIONS

Unfolded: L 84 in. x W 52 in.
Weight: 1.7 oz



NAR® Responder QUIKLITTER

- Multi-use, low cost disposable evacuation platform for effective emergency casualty transport & transfer
- Constructed of a latex free water-resistant, woven polyethylene with signature red nylon webbing & stitching
- Ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points



NAR® RESPONDER QUIKLITTER

ITEM

60-0080

SPECIFICATIONS

Litter Bed: L 74 in. x W 28 in. | Weight: 13.3 oz
Payload 500 lb.

Emergency SURVIVAL WRAP

- Lightweight, durable, laminated polyethylene
- Reversible Silver/Orange or OD Green



SURVIVAL WRAP

COLOR	ITEM #
Orange/Silver	ZZ-0834
OD Green/Silver	ZZ-1390

SPECIFICATIONS

Packaged: L 5 in. x W 3.75 in.
Open Wrap: L 96 in. x W 60 in.
Weight: 3.7 oz

HeavyDuty Trauma **XSHEAR®**

- Hardened stainless steel blades that are twice as thick as most trauma shears, sharpened to a razor edge
- Patented blunt tip and curved design for gentler edge near the skin
- Soft touch inner handles for comfort and slip resistance
- Heavy duty center bolt that won't loosen over time

XSHEAR

HANDLE	BLADE	ITEM #
Black	Black	70-0290
Blue	Silver	70-0291
Orange	Black	70-0292
Pink	Black	70-0293
Red	Black	70-0294
Coyote	Black	70-0415
Green	Silver	70-0424

SPECIFICATIONS

L7.5 in x 3.5 in. x 0.118 in.
Cutting area length: 2 in.
Weight: 4 oz.



Shears & Holsters are sold separately

Rigid **XSHEAR® HOLSTER**

- Custom designed for the XSHEAR® Heavy Duty Trauma Shears
- MOLLE Compatible adjustable belt clip secured holster in place and offers flexible orientation of holster with simple Philips head screw
- Flex lock design securely holds shears in place
- Positive click when shears are secured in place
- Universal design allows shears to be inserted either direction
- Rugged injection molded ABS Plastic

SPECIFICATIONS

L: 5 in. x W 3.75 in. x D 1.75 in.
Weight: 3.2 oz

RIGID XSHEAR® HOLSTER

ITEM

70-0295



Soft **XSHEAR® HOLSTER**

- Designed for XSHEARS but compatible with most shears
- Durable Spandex over soft foam with nylon interior
- Powder coated heavy steel belt clip
- Ambidextrous carry

SOFT XSHEAR® HOLSTER

ITEM

70-0314

SPECIFICATIONS

L 7 in. x W 3.5 in. x D 0.5 in.
Weight: 3.2 oz



SALTY BRITCHES® SKIN BARRIER

- Offers soothing, long-lasting barrier against skin irritation from windburn, friction, and chafing
- Waterproof
- Won't rub off

SPECIFICATIONS

Tubes

Packaged: W 2.5 in. x L 4 in. | Wt: 0.15 lbs

Single Use (10 pk)

Packaged: W 3 in. x L 5 in. | Wt: 0.10 lbs

SALTY BRITHCES®

OPTION	ITEM #
2 oz. Skin Barrier	75-0014
2 oz. Winter Relief	75-0015
0.1 oz. Single Use*	75-0017



*Pack of 10

Responder 7.25 in. **TRAUMA SHEARS**



Large, premium-grade Trauma Shears designed for safety and durability while cutting through tough materials

RESPONDER TRAUMA SHEARS

ITEM

ZZ-1019

SPECIFICATIONS

L 7.25 in. x W 3.7 in. x D 0.4 in.
Weight: 2 oz



COMMUNITY
PREPAREDNESS



**Redefining
FIRST AID**

NORTH AMERICAN RESCUE

NARescue.com • 888.689.6277

PRODUCTS WITH A MISSION®



BE THE DIFFERENCE®

www.NARescue.com



PUBLIC ACCESS BLEEDING CONTROL KITS | ANSI WORKPLACE FIRST AID TRAUMA KITS

REV010824