

WORKPLACE SAFETY

OSHA COMPLIANCE

TRAUMA FIRST AID KITS



PUBLIC ACCESS BLEEDING CONTROL KITS | ANSI WORKPLACE FIRST AID TRAUMA KITS



COMMUNITY PREPAREDNESS



*The first person on scene
IS the first responder...*

Are YOU Prepared?



3 to 5 MINUTES

*Kit shown:
Class B Trauma &
First Aid Kit-Hard Case
see page 5*

The Difference between Life and Death

- ▶ A victim who is bleeding can quickly die BEFORE professional emergency responders can even arrive on the scene
- ▶ Serious bleeding accounts for 35% of prehospital deaths from traumatic injuries
- ▶ Life-threatening bleeding warrants immediate intervention
- ▶ Basic bleeding control provided by bystanders can often be the difference between life and death...

***FACT:** from 2000-2020,
45% of **ALL Active Shooter Incidents**
 occurred in **Places of Commerce**



*Death from
 Uncontrolled
BLEEDING
 can occur in as
 little as **3 to 5 Minutes***

Trauma is a national health problem and after a traumatic injury, serious bleeding is responsible for over 35% of prehospital deaths.** Causes of trauma can range from natural disasters, falls, traffic, industrial, and workplace accidents to intentional acts of violence such as shootings, stabbings and bombings. In fact, according to data collected from 2000 to 2020, the U.S. Department of Justice found that the workplace is the most likely location of an active shooter situation.*

Because these traumatic injuries often result in uncontrolled bleeding,

immediate intervention – even before emergency medical services or other first responders can arrive – is critical to avoid potential death. **North American Rescue®** is leading the way to equip bystanders and initial first responders with the essential medical gear needed to take immediate action in life threatening emergencies. These lifesaving modular aid components are designed for individuals or groups in the workplace, public buildings, and populated public areas to administer immediate bleeding control at the point of wounding.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

* Active Shooter Incidents 20-Year Review, 2000-2020, Federal Bureau of Investigation, U.S. Department of Justice, Washington, D.C.

**Chambers JA, Seastedt K, Krell R, Caterson E, Levy M, Turner N. "Stop the Bleed": A U.S. Military Installation's Model for Implementation of a Rapid Hemorrhage Control Program. Mil Med. 2019 Mar 01;184(3-4):67-71. [PubMed] [Reference list] <https://academic.oup.com/milmed/article/184/3-4/67/5062295>

Workplace First Aid Kits

Class A **TFAK** (TRAUMA & FIRST AID KIT)

- Portable, Bleeding Control & First Aid Kit designed for the workplace, **exceeding ANSI/ISEA Z308.1-2021 Standard for Class A kits**
- Able to **treat the most common workplace injuries** with the added benefit of being able to **control major bleeding** with the C-A-T® Tourniquet - the safest and most effective Pre-Hospital Field Tourniquet on the market!
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

TRAUMA & FIRST AID KIT - CLASS A

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0947
Soft Case	Bleeding Control Dressing	80-1045
Soft Case	CELOX™ First Aid Gauze	80-1477
Hard Case	Compressed Gauze	80-1033
Hard Case	Bleeding Control Dressing	80-1046
Hard Case	CELOX™ First Aid Gauze	80-1481

SPECIFICATIONS

Soft Case::

L 14.75 in. x H 8 in. x D 4 in.
Weight: 2.94 lb

Hard Case::


H 6 in. x W 13.5 in. x D 11.5 in.
Weight: 6 lb

Type IV ANSI Container only option (#80-1113) also available at NARescue.com

TRAUMA & FIRST AID KIT CLASS A:

 QTY 1 C-A-T® Tourniquet	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD® 4 in. Responder, Flat	 QTY 1 Face Shield Microshield Mask	 QTY 1 HyFin® Vent Compact, Twin Pk	 QTY 1 NAR® Survival Blanket	 QTY 1 Triangular Bandage
 QTY 1 SAM® Splint II	 QTY 1 4" Elastic Wrap Bandage	 QTY 1 2 in. Elastic Wrap Bandage	 QTY 2 4" x 4" Gauze Pad 2 PK	 QTY 2 5" x 9" Gauze Pad	 QTY 1 BurnTec® 4" x 4" Dressing	 QTY 10 Safetec Burn Gel
 QTY 10 MED NAP Antiseptic Towelette	 QTY 10 Triple Antibiotic Ointment	 QTY 10 Hand Sanitizer	 QTY 1 COLDSTAR 6" x 9" Cold Pack	 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 Sterile Oval Eye Pad	 QTY 1 Surgical Tape, 1 in.
 QTY 16 1" x 3" Bandage	 PAIR 2 Responder Blue Gloves, LG	 QTY 1 Responder Shears, Lg	 QTY 1 Permanent Marker, Large	 QTY 1 First Aid Pocket Guide	 QTY 1 Instructions Bleeding Control	 QTY 1 Kit Contents Inventory Card
 QTY 1 Inspection Card						

PLUS+
one of these:



RESUPPLY - CLASS A

DRESSING	ITEM #
Compressed Gauze	80-1153
Bleeding Control Dressing	80-1155
CELOX™ First Aid Gauze	80-1479

Container not included with ReSupply Kits



Workplace First Aid Kits

Class B **TFAK** (TRAUMA & FIRST AID KIT)

- Portable, Bleeding Control and First Aid Kit designed for the workplace, **exceeding ANSI/ISEA Z308.1-2021 Standard for Class B kits**
- Same supplies as the Class A Kit but with larger quantities to accommodate larger workplaces
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

SPECIFICATIONS

Soft Case:

L 14.75 in. x H 8 in. x D 7.5 in.
Weight: 4.87 lb

Hard Case:

H 6 in. x W 13.5 in. x D 11.5 in.
Weight: 9 lb

TRAUMA & FIRST AID KIT - CLASS B

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0948
Soft Case	Bleeding Control Dressing	80-1047
Soft Case	CELOX™ First Aid Gauze	80-1478
Hard Case	Compressed Gauze	80-1034
Hard Case	Bleeding Control Dressing	80-1048
Hard Case	CELOX™ First Aid Gauze	80-1482

Type IV ANSI Container only option (#80-1114) also available at NARescue.com

Class A/B resupply modules also available at NARescue.com



TRAUMA & FIRST AID KIT CLASS B:

 QTY 1 C-A-T® Tourniquet	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD® 4 in. Responder, Flat	 QTY 1 CPR Microshield Mask	 QTY 1 HyFin® Vent Compact Twin	 QTY 1 NAR® Survival Blanket	 QTY 2 Triangular Bandage
 QTY 1 SAM® Splint II	 QTY 1 4" Elastic Wrap Bandage	 QTY 2 2" Elastic Wrap Bandage	 QTY 4 4" x 4" Gauze Pad 2 PK	 QTY 4 5" x 9" Gauze Pad	 QTY 2 BurnTec® 4" x 4" Dressing	 QTY 25 Safetec Burn Gel
 QTY 50 Antiseptic Towelette	 QTY 25 Triple Antibiotic Ointment	 QTY 20 Hand Sanitizer	 QTY 2 6" x 9" Cold Pack	 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 Sterile Oval Eye Pad	 QTY 2 Surgical Tape, 1 in.
 QTY 50 1" x 3" Bandage	 PAIR 4 Responder Blue Gloves, LG	 QTY 1 Responder Shears, Lg	 QTY 1 Permanent Marker, Large	 QTY 1 First Aid Pocket Guide	 QTY 1 Instructions Bleeding Control	 QTY 1 Kit Contents Inventory Card
					 QTY 1 Inspection Card	

RESUPPLY- CLASS B

DRESSING	ITEM #
Compressed Gauze	80-1154
Bleeding Control Dressing	80-1156
CELOX™ First Aid Gauze	80-1480

Container not included with ReSupply Kits



PLUS+
one of these:

w/ Comp. Gauze QTY 1 Responder Compressed Gauze	- OR -	w/ BCD QTY 1 Bleeding Control Dressing	- OR -	w/ CELOX™ QTY 1 5' CELOX Hemostatic Dressing
--	--------	---	--------	---

Workplace First Aid Kits

LOGGING FIRST AID & TRAUMA KIT *Hard Case*

- Rugged case designed for use in the harshest of conditions
- Modular system featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels (**Meets Type IV ANSI Container**)
- Fast open dual-latch closure system allows for rapid access while keeping contents secure
- Able to **treat the most common job site injuries** with the added benefit of being able to **control major bleeding** with the C-A-T® Tourniquet - the fastest, safest, most effective pre-hospital field tourniquet on the market
- Designed with mobility in mind, the kit can be easily carried to the patient's side

LOGGING FIRST AID & TRAUMA KIT - HARD CASE

SPECIFICATIONS

H 6 in. x W 13.5 in. x D 11.5 in.
Weight: 5.66 lb

OPTION	ITEM #
Basic	80-1190
with Bleeding Control Dressing	80-1191

LOGGING FIRST AID & TRAUMA KIT HARD CASE:

 QTY 1 C-A-T® Tourniquet	80-1190  QTY 2 Responder Compressed Gauze	-OR- 80-1191  QTY 1 Responder Compressed Gauze	 QTY 1 QuikClot® BCD
 QTY 1 ETD® 4 in. Responder, Flat	 QTY 2 8" x 10" Gauze Pad 2 PK	 QTY 2 4" x 4" Gauze Pad 2 PK	 QTY 2 2 in. Elastic Wrap Bandage
 QTY 2 Triangular Bandage	 QTY 16 1" x 3" Bandage	 QTY 1 BurnTec® 4" x 4" Dressing	 QTY 1 SAM® Splint II
 PAIR 2 Responder Blue Gloves, LG	 QTY 10 Antiseptic Towelette	 QTY 10 Sting & Bite Antiseptic	 QTY 1 Microshield CPR Mask
 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 PES Eye Shield	 QTY 1 Instant Cold Compress	 QTY 1 2" Mini Surgical Tape
 QTY 1 Splinter Forceps, 4.5"	 QTY 1 Responder Shears, Lg	 QTY 1 NAR® Survival Blanket	 QTY 1 Instructions for Use PUBLIC ACCESS BLEEDING CONTROL
 QTY 1 First Aid Pocket Guide	 QTY 1 Inspection/ Content Card		



Workplace First Aid Kits

LOGGING FIRST AID & TRAUMA KIT SOFT CASE:



LOGGING FIRST AID & TRAUMA KIT Soft Case


- Portable Bleeding Control and First Aid Kit designed for the workplace that **exceeds OSHA Standard Number: 1910.266 App A - First-aid Kits** (Mandatory) for the logging industry
- Compact, rugged nylon case keeps contents secure while stored in vehicles and machinery
- Able to **treat the most common job site injuries** with the added benefit of being able to **control major bleeding** with the C-A-T® - the fastest, safest, most effective pre-hospital field tourniquet on the market
- Designed with mobility in mind, the kit can be easily carried to the patient's side

SPECIFICATIONS

L 14.75 in. x H 8 in. x D 4 in.
Weight: 3.03 lb

LOGGING FIRST AID & TRAUMA KIT - SOFT CASE

OPTION	ITEM #
Basic	80-1188
with BCD	80-1189



BLOODBORNE PATHOGEN

PERSONAL PROTECTION KIT

- PPE in an ultra-portable, re-sealable pouch
- Complies with Title 29 of the Code of Federal Regulations at 29 CFR 1910.1030.

BLOODBORNE PATHOGEN PP KIT

ITEM


70-0426

SPECIFICATIONS

Dimensions: L 8 in. x W 5.5 in. x D 1.5 in. | Weight: 0.35 oz



BLOODBORNE PATHOGEN PP KIT contains:



MINE SAFETY AND HEALTH FIRST AID KIT



- First Aid Kit that exceeds **MSHA's (Mine Safety & Health Administration) § 75.1713-7** requirements for first aid equipment in mining environments.
- Able to **treat the most common workplace injuries** with the added benefit of being able to control major bleeding with the
- C-A-T® Tourniquet, the **fastest, safest, most effective** prehospital field tourniquet in the world
- Kit contains **bleeding control, splinting and evacuation** equipment
- Designed so the kit can be **easily transported** to the patient's side

SPECIFICATIONS

External: H 32.5 in. x W 21 in. x D 13 in.

Internal: H 29 in. x W 18 in. x D 10 in.













Weight: 53.6 lbs.

MINE SAFETY & HEALTH FIRST AID KIT

ITEM

80-1791

MINE SAFETY & HEALTH FIRST AID KIT:

 QTY 2 C-A-T® Tourniquet	 QTY 2 ETD® 4 in. Responder, Flat	 QTY 2 NAR® Gauze, Compressed	 QTY 8 4" Elastic Wrap Bandage	 QTY 8 2" Elastic Wrap Bandage
 QTY 24 Triangular Bandage	 QTY 2 SAM® Splint II	 QTY 2 Wool Blanket	 QTY 4 BurnTec® 4"x4" Dressing	 QTY 12 Ammonia Inhalant
 QTY 1 Foilie Ointment, 1 oz.	 QTY 100 1" x 3" Bandage	 QTY 12 1" Coban Wrap Bandage	 QTY 2 Responder Shears, Lg	 QTY 1 MedEVAC Quad-Fold Litter
			 QTY 1 6'x6' Thermal Resistant Blanket	

ZOLL® AED AUTOMATIC DEFIBRILLATORS



AED 3 Contents

- 1x ZOLL® AED 3 Automatic Defibrillator
- 1x CPR Uni-Padz™ Electrodes with Accessory Kit
 - » 1x Scissors
 - » 1x Disposable Razor
 - » 1x Face Barrier/Shield
 - » 1x Antimicrobial wipe
 - » 1x Dry Towel
 - » 1x ZOLL® AED 3 Battery Pack

ALL ZOLL AED AUTOMATIC DEFIBRILLATORS: INCLUDE MEDICAL PRESCRIPTION / AUTHORIZATION WITH PURCHASE

- Provide detailed CPR guidance that supports rescuers with integrated text, voice, and visual prompts
- Easy to replace 123 batteries
- Longer-lasting 5-year battery and electrode pads reduce AED device maintenance requirements
- **Device Classification:** Class III and internally powered per EN60601-1

ZOLL AED 3:

- **AED 3:** RapidShock® Analysis enables short rhythm analysis for continuous care and more life-saving CPR
- **AED 3:** Pediatric rescue is made easier with ZOLL's unique CPR Uni-padz® electrodes and a child-mode setting
- **AED 3:** Equipped with optional Wi-Fi connectivity, allowing seamless communication with AED program management systems through your wireless network
- **AED 3 Design Standards:** Meets applicable requirements of EN 60601-1, IEC 60601-1-11, IEC 60601-2-4

ZOLL AED PLUS:

- **AED Plus:** Includes ZOLL's CPR-D-padz® electrodes, which feature a unique one-piece, pull-strip design to ensure fast and accurate placement
- **AED Plus:** Fully Automatic AED
- **AED Plus Design Standards:** Meets applicable requirements of UL 2601, AAMI DF-39, IEC 601-2-4, EN60601-1, & IEC60601-1-2

ZOLL AED AUTO DEFIBRILLATORS

OPTION	ITEM #
ZOLL AED 3	75-0024
ZOLL AED PLUS	75-0019

Find **ZOLL®** carrying
case and other
AED accessories at
NARescue.com



SPECIFICATIONS

ZOLL AED 3:

L 5 in. x W 9.3 in. x D 9.7 in.
Weight: 5 lb 8 oz. (w/ battery)

ZOLL AED PLUS:

L 5.25 in. x W 9.5 in. x D 11.5 in.
Weight: 6 lb 11.2 oz. (w/ battery)



AED Plus Contents

- 1x ZOLL® AED Plus Automatic Defibrillator
- 1x Carrying Case - Black
- 1x CPR-D-Padz® Electrodes with Accessory Kit
 - » 1x Scissors
 - » 1x Disposable Razor
 - » 1x Face Barrier/Shield
 - » 1x Antimicrobial wipe
 - » 1x Dry Towel
 - » 1x Type 123A Photo Flash lithium manganese dioxide batteries (10 pk)



- The latest cooling technology for effective Human Thermo-regulation
- Polar Skin Mini-Mesh,™ a unique smart textile with a hydrofunctional polymer that allows for liquid absorptive capacity 350% greater than base weight
- ArcticAqua,™ a proprietary solution that prevents rigidity at temperatures well below 0°C while providing the requisite moisture to support evaporative, convective and conductive cooling
- Adheres to current U.S. Army and National Athletic Trainers' Association protocols for heat illness treatment
- Protracted direct cooling of high blood flow regions

COMPREHENSIVE COOLING SYSTEM

- Highly-effective portable cooling system with C.O.R.E. technology for the treatment of heat-related illnesses
- Scientifically proven to reduce core temperature by .09°F/minute in 100°F heat when used according to protocol

COMPREHENSIVE COOLING

ITEM

ZZ-0674

SPECIFICATIONS

L 16 in. x W 12.5 in. x H 3.25 in.
Weight: 9 lbs



COMPLETE CARE PACK

- Portable, highly insulated option
- Contains one Polar Skin Comprehensive Cooling System and one heavy-duty, leak & sweat proof cooler that keeps ice from melting for 24 hours in +100°F heat

COMPLETE CARE PACK

ITEM

80-0451

SPECIFICATIONS

W 14 in. x L 12 in. x D 7 in.
Weight: 10 lb 3.2 oz



COMPLETE CARE BACKPACK

Available in either Single Patient (one comprehensive cooling system) or Double Patient (two cooling systems) configurations

COMPLETE CARE BACKPACK

SIZE	ITEM #
Single	80-0543
Double	80-0544

SPECIFICATIONS

W 19.5 in. x L 15 in. x D 7.5 in.
Single Weight: 12 lb 13 oz
Double Weight: 22 lb



COOLERS *Available in Standard, Large and Backpack*

Standard size accommodates one Polar Skin™ Comprehensive Cooling System and Large or Backpack Size accommodates two cooling systems

POLAR SKIN COOLERS

COLOR	SIZE	ITEM #
Red	Standard	ZZ-0682
Black	Standard	ZZ-0759
Black	Backpack	ZZ-0794



SPECIFICATIONS

Standard: L 14 in. x W 12 in. x H 7 in.
Weight: 2.2 lb.
Backpack: L 19.5 in. x W 15 in. x D 7.5 in.
Weight: 3.8 lb.



ICE SHEETS *Standard*

- Principle cooling component of Polar Skin™ cooling systems, featuring C.O.R.E. technology
- Adheres to current National Athletic Trainers' Assoc and U.S. Army protocols for heat illness treatment
- Sheet Size: W 36 in. x L 30 in.



STANDARD ICE SHEETS

UNIT	ITEM #	NSN #
Single Pack	ZZ-0668	6532-01-653-6054
Double Pack	ZZ-0669	----
Triple Pack	ZZ-0755	----

SPECIFICATIONS

Single: L 4.5 in. x W 7.75 in. x H 0.5 in. Wgt: 9.6 oz
 Double: L 9.5 in. x W 8 in. x H 0.5 in. Wgt: 1 lb 6.4 oz
 Triple: L 12.5 in. x W 7.5 in. x H 0.5 in. Wgt: 2 lb 3.2 oz

ICE SHEETS *Large*

Same characteristics as the standard size Ice Sheets, but in a larger 36 in. x 60 in. size

LARGE ICE SHEETS

SPECIFICATIONS

ITEM #

L 8.5 in. x W 5.5 in. x H 1.5 in.
 Weight: 1 lb 8 oz

ZZ-0758

C2E ICE PACK BUNDLE

- Four cooling pads featuring C.O.R.E. Technology
- Facilitates Conductive, Convective, and Evaporative (C2E) cooling of high blood flow (large artery) regions of the body



C2E ICE PACK

SPECIFICATIONS

ITEM #

Indiv. Gel Pk: L 15 in. x W 11 in. x D 1.5 in.
 Weight: 4.6 lb

ZZ-0675

COOLING COLLAR

- Two gel packs encased in integrated pouches
- Designed to align with the carotid artery

COOLING COLLAR

SPECIFICATIONS

ITEM #

Indiv. Gel Pk: W 6 in. x L 10 in.
 Pkg: L 8 in. x W 9.5 in. x H 1.3 in.
 Weight: 2.7 lb

ZZ-0670

Burn Treatment



BurnTec® Hydrogel Dressings

- Sterile non-stick hydrogel dressings for treating 1st, 2nd & 3rd-degree burns (thermal, chemical, electric)
- Also treats other skin injuries, such as sunburns, abrasions, skin irritations, bruises and swelling (insect bites)
- Applied dressings provide a cool, soothing effect for up to 24 hours
- Easy to handle without spilling or damaging the hydrogel
- Permeable to water vapor and oxygen but impermeable to bacteria
- Biocompatible and hypo-allergenic; Reduces risk of hypertrophic scars

BURNTec DRESSINGS

SIZE	ITEM #	NSN #
2 in. x 2 in.	30-0114	6510-01-631-1913
2.5 in. x 5 in.	30-0115	6510-01-631-2080
4 in. x 4 in.	30-0110	6510-01-631-1883
5 in. x 5 in.	30-0116	6510-01-631-2138
5 in. x 10 in.	30-0117	6510-01-631-1907
8 in. x 8 in.	30-0118	6510-01-631-1897
8.5 in. x 11 in.	30-0119	6510-01-631-1890
4 in. x 16 in.	30-0111	6510-01-631-2124
8 in. x 16 in.	30-0112	6510-01-626-2990
16 in. x 24 in.	30-0120	6510-01-631-1886
10 in. x 10 in. Face	30-0121	6510-01-627-1834
12 in. x 16 in. Face	30-0113	6510-01-627-1855
16 in. x 24 in. Face	30-0122	----



BurnTec® EMERGENCY BURN KIT

- Easy to handle without spilling or damaging the hydrogel
- Does not leave a mess like other burn care products
- Cool, soothing effect for up to 24 hours
- Impermeable to bacteria
- Biocompatible and hypoallergenic
- Reduces risk of hypertrophic scars
- Award winning technology

Recommended Uses

- First, second and third-degree burns (thermal, chemical, electric)
- Sunburns
- Abrasions and skin irritations
- Bruises and swelling (insect bites)
- Other skin injuries



SPECIFICATIONS

Small Kit

Dimensions: L 8 in. x W 14 in. x D 3 in.
Weight: 3 lbs 3.2 oz

Large Kit

Dimensions: L 10 in. x W 20 in. x D 3 in.
Weight: 4 lbs. 16 oz.

BURNTec EMERGENCY BURN KIT - SMALL



BURNTec EMERGENCY BURN KIT

OPTION	ITEM #	NSN #
Small Kit	80-1661	----
Large Kit	80-1662	----
Small Case	30-0150	----

BURNTec EMERGENCY BURN KIT - LARGE



CRITICAL BURN KIT

- Easy to handle without spilling or damaging the hydrogel
- Does not leave a mess like other burn care products
- Cool, soothing effect for up to 24 hours
- Impermeable to bacteria
- Biocompatible and hypoallergenic
- Reduces risk of hypertrophic scars
- Award winning technology

BURNTec CRITICAL BURN KIT



SPECIFICATIONS

Dimensions: L 8 in. x W 14 in. x D 3 in.
Weight: 3 lbs. 3.2 oz



online video
available

◀scan to view

BURNTec CRITICAL BURN KIT

OPTION	ITEM #	NSN #
Kit	80-1663	----
Small Case	30-0150	----

BurnTec®

BURN DRESSING KIT

- Convenient array of BurnTec® dressings in a handy carrying case
- Kit contains (1) each of the following dressings:

2 in. x 2 in.	5 in. x 10 in.	8 in. x 16 in.
2.5 in. x 5 in.	8 in. x 8 in.	10 in. x 10 in. (face)
4 in. x 4 in.	4 in. x 16 in.	12 in. x 16 in. (face)
5 in. x 5 in.	8.5 in. x 11 in.	

BURNTec BURN KIT

OPTION	ITEM #	NSN #
Kit	80-0361	6515-01-630-9177
Bag Only	30-0123	6510-01-626-6079

SPECIFICATIONS

L 12 in. x W 15.5 in. x D 3.5 in.
Weight: 6 lb 4 oz



BurnTec® MINOR BURN DRESSING KIT

- Variety of BurnTec® dressings in a handy carrying case
- Kit contains the following:
Two 4 in. x 4 in. dressings
Two 4 in. x 16 in. dressings
One 8.5 in. x 11 in. dressing

BURNTec MINOR BURN KIT

OPTION	ITEM #	NSN #
Kit	80-0415	6510-01-644-5427
Bag Only	30-0150	----

SPECIFICATIONS

L 14 in. x W 8.5 in. x D 2.5 in.
Weight: 3 lbs. 16 oz.



Personal First Aid Kits

I-ROK® Individual Rugged Outdoor Kit

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding
- Pelican R40 ruck case is crushproof, dustproof and waterproof up to 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- Available in 4 hard case colors and 2 sizes with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions



I-ROK®

COLOR	LEVEL	ITEM #
Black	Basic	82-0031
OD Green	Basic	82-0033
Black	w/BCD	82-0035
OD Green	w/BCD	82-0037
Black	w/CELOX	82-0066
OD Green	w/CELOX	82-0068

SPECIFICATIONS

Int: L 7.6 in. x W 4.7 in. x D 1.9 in.
Ext: L 9.8 in. x W 6.1 in. x D 2.9 in.
Weight: 2.48 lb.

Every I-ROK contains the following items:

 QTY 1 C-A-T® Tourniquet	 QTY 1 Responder Compressed Gauze	 QTY 1 Choice of Hemostatic	 QTY 1 NAR® Mini Responder 4" ETD®	 QTY 1 HyFin® Vent Compact, Twin Pk				
 QTY 5 Butterfly Bandage	 QTY 5 1" x 3" Bandage	 QTY 3 Safetec Burn Gel	 QTY 1 Microshield CPR Mask	 QTY 3 Antibiotic Ointment	 QTY 2 Sting & Bite Antiseptic	 QTY 1 NAR® Survival Blanket	 QTY 3 Hand Sanitizer	 QTY 3 Salty Britches® Chafing Ointment
 PAIR 1 Responder Blue Gloves, LG	 QTY 1 NAR® Mini Duct Tape, 2 in.	 QTY 1 Permanent Marker, Large	 QTY 1 Instructions for Use: BLEEDING CONTROL					

I-ROK® XL

- Extra supplies to accommodate multiple injuries
- Large, Pelican R60 Ruck Case w/ removable tray

I-ROK XL includes these addl. items:

 QTY 5 Butterfly Bandage	 QTY 5 1" x 3" Bandage	 QTY 2 Safetec Burn Gel						
 QTY 1 BurnTec® 4"x4" Dressing	 QTY 2 Antibiotic Ointment	 QTY 1 Sting & Bite Antiseptic	 QTY 2 Hand Sanitizer	 QTY 2 4" x 4" Gauze Pad 2 PK	 QTY 1 Triangular Bandage	 PAIR 1 Responder Blue Gloves, LG	 QTY 1 Splinter Forceps, 4.5"	 QTY 1 NAR® Trauma Shears, Lg

I-ROK® XL

COLOR	LEVEL	ITEM #
Black	Basic	82-0049
OD Green	Basic	82-0051
Black	w/BCD	82-0053
OD Green	w/BCD	82-0055
Black	w/CELOX	82-0070
OD Green	w/CELOX	82-0072

SPECIFICATIONS

Int: L 8.9 in. x W 5.4 in. x D 3 in.
Ext: L 11.1 in. x W 6.9 in. x D 3.9 in.
Weight: 3 lb.



Personal First Aid Kits

Mini First Aid Kit

M-FAK®

RED



Other colors available at
NARescue.com

SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz

M-FAK is compatible with
NAR® Quick Detach Panel

M-FAK

LEVEL	ITEM #
Basic	80-0890
with BCD	80-0891

- NAR's most compact, versatile Mini First Aid Kit
- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries
- Rugged nylon platform allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

M-FAK Contents:



PLUS+
one of these:



DOOR PANEL KIT

- Slim design to fit in vehicle door panel storage bin
- Equipped to treat major bleeding as well as minor injuries, burns, lacerations, etc.

SPECIFICATIONS

H 7 in. x W 7 in. x D 15 in.
Weight: 1 lb 6 oz

DOOR PANEL KIT

LEVEL	ITEM #
Basic	80-1612
with BCD	80-1614
with CELOX™	80-1616

Also available in black at **NARescue.com**



HEADREST IFAK

COLOR	LEVEL	ITEM #
Black	Basic	80-1091
Black	with BCD	80-1151
Coyote	Basic	80-1150
Coyote	with BCD	80-1152

SPECIFICATIONS

H 5.25 in. x W 7.25 in.
x D 4.25 in.
Weight: 1 lb. 9.9 oz.

- Conveniently mounts to most vehicle headrests



HEADREST IFAK



Public Access **INDIVIDUAL BLEEDING CONTROL KITS**

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control ID Patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.5 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.



**Essential equipment
needed for persons to
take action as immediate
responders in helping
to stop life threatening
bleeding in Public Areas...**

Every **INDIVIDUAL BLEEDING CONTROL KIT** contains:




**PLUS these
Individual
Components
based
upon YOUR
configuration
choice:**



PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KITS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	80-0452
Intermediate	Red Nylon	16.5 oz	80-0453
Advanced BCD	Red Nylon	17.4 oz	80-0523
Advanced	Red Nylon	17.4 oz	80-0454
Basic	Vacuum Sealed	11.7 oz	80-0465
Intermediate	Vacuum Sealed	13.25 oz	80-0466
Advanced BCD	Vacuum Sealed	14 oz	80-0522
Advanced	Vacuum Sealed	14 oz	80-0467

**All Intermediate/
Advanced Public
Access Bleeding
Control Kits and
Stations are Texas
House Bill 496
Compliant**



**PUBLIC ACCESS
BLEEDING CONTROL 8-PACK NYLON**



**PUBLIC ACCESS
BLEEDING CONTROL 8-PACK VACUUM**



Public Access

BLEEDING CONTROL 8-PACKS

Nylon | Vacuum-Sealed

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon carrying case
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced, & Advanced BCD

SPECIFICATIONS

8-Pack: H 14 in. x W 11 in. x D 8 in.

Each Individual Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in.

PUBLIC ACCESS BLEEDING CONTROL 8-PACKS

LEVEL	BAG	WGT	ITEM #
Basic	Red Nylon	8.5 lb	80-0460
Intermediate	Red Nylon	9.3 lb	80-0461
Advanced BCD	Red Nylon	9.7 lb	80-0525
Advanced	Red Nylon	9.7 lb	80-0462
Basic	Vacuum Sealed	6.8 lb	80-0468
Intermediate	Vacuum Sealed	7.6 lb	80-0469
Advanced BCD	Vacuum Sealed	8 lb	80-0524
Advanced	Vacuum Sealed	8 lb	80-0470



Public Access **BLEEDING CONTROL 6-PACKS** Vacuum-Sealed

- Designed to provide bystanders and initial first responders with quick and easy access to essential medical equipment for stopping life-threatening bleeding
- Includes an easy-to-open zippered tote containing six (6) Vacuum Sealed Individual Bleeding Control Kits packaged in ruggedized vacuum sealed pouches and one (1) NAR® QuikLitter™ LITE for moving injured casualties
- Available in Basic, Intermediate, and Advanced care configurations

PUBLIC ACCESS BLEEDING CONTROL 6-PACK VACUUM



PUBLIC ACCESS BLEEDING CONTROL 6-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Vacuum Sealed	5.7 lb	80-0573
Intermediate	Vacuum Sealed	5.75 lb	80-0574
Adv.	Vac Vacuum Sealed	6.5 lb	80-0575
Adv. BCD	Vacuum Sealed	6.5 lb	80-0576
Adv. w/ CELOX	Vacuum Sealed	6.5 lb	80-1571

All Int./Adv. Public Access Bleeding Control Kits and Stations are
Texas House Bill 496 Compliant

SPECIFICATIONS

6-Pack:

H 7.5 in. x W 13 in. x D 8 in..

Each Individual
Bleeding Control Kit:

H 7 in. x W 5 in. x D 2.5 in

Trauma & First Aid

ESSENTIAL BLEEDING CONTROL KIT

Every **ESSENTIAL BLEEDING CONTROL KIT** contains:



PLUS+
one of these:



- Compact and Convenient Kit
- Highly mobile and easily deployed
- Hypothermia prevention included
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market



SPECIFICATIONS

L 12 in. x H 7.5 in.

Weight: 0.68 lbs

TRAUMA & FIRST AID ESSENTIAL BLEEDING CONTROL KIT

OPTION	ITEM #
Compressed Gauze	80-1474
Bleeding Control Dressing	80-1475
CELOX™ First Aid Gauze	80-1476

Clear Polycarbonate WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 6 & 8 Pack Stations also include an easy to-open nylon carrier and one (1) Responder QuikLitter™ (6 Packs contain QuikLitter™ Lite) for moving injured casualties*

SPECIFICATIONS

6 & 8 Pack Stations: H 18 in. x W 15.5 in. X D 10 in.
Empty Case Wt: 15 lb (includes rack)
Large Case Item# ZZ-0822

Individual Stations: H 10 in. x W 10 in. X D 5.25 in.
Empty Case Wt: 2.39 lb
Small Case Item# ZZ-1000

8-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	80-0471
Intermediate	8	Vacuum	22.6 lb	80-0472
Adv.	8	Vacuum	23 lb	80-0473
Adv. BCD	8	Vacuum	23 lb	80-0526
Adv. w/ CELOX	8	Vacuum	23 lb	80-1575
Basic	8	Nylon	23.5 lb	80-0474
Intermediate	8	Nylon	24.3 lb	80-0475
Adv.	8	Nylon	27.4 lb	80-0476
Adv. BCD	8	Nylon	27.4 lb	80-0527
Adv. w/ CELOX	8	Nylon	27.4 lb	80-1576

6-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	6	Vacuum	20.7 lb	80-0581
Intermediate	6	Vacuum	20.7 lb	80-0582
Adv.	6	Vacuum	21.5 lb	80-0583
Adv. BCD	6	Vacuum	21.5 lb	80-0584
Adv. w/ CELOX	6	Vacuum	21.5 lb	80-1574

IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	3.3 lb	80-0938
Intermediate	1	Vacuum	3.42 lb	80-0939
Adv. BCD	1	Vacuum	3.48 lb	80-0940
Adv. w/ CELOX	1	Vacuum	3.48 lb	80-1580
Basic	1	Nylon	3.3 lb	80-0935
Intermediate	1	Nylon	3.6 lb	80-0936
Adv. BCD	1	Nylon	3.7 lb	80-0937
Adv. w/ CELOX	1	Nylon	3.7 lb	80-1581

Alarm (#ZZ-1075) is also available for the Clear Polycarbonate Wall Cabinet

All Intermediate/Advanced Public Access Bleeding Control Kits and Stations are **Texas House Bill 496 Compliant**



8 Pack Station



Individual Kit Station



Small Individual Clear Case: ZZ-1000

Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs) to provide publically accessible medical supplies to control major bleeding



6 Pack Station



Large Multi-Pack Clear Case: ZZ-0822

CLEAR CASE
features

- Durable, clear polycarbonate cabinet provides both physical and environmental protection for medical equipment
- Large case features wire shelf to support the contents
- External eyelet allows a secondary, tamper-evident safety seal
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle



**Low Profile
Metal Cabinet:
ZZ-1199**

Designed to be mounted on a wall with a slim 4 in. clearance and featuring easy pull-to-open design for rapid public access to major bleeding control supplies

Low Profile Metal WALL STATIONS

8 Pack Nylon Station

- Low profile metal cabinet mounts on the wall with a slim 4 in. clearance
- Designed to safely store 8 Individual Bleeding Control Kits in public areas
- Kits are available in red nylon bags or Vacuum-Sealed pouches and stored in Two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon handle
- Features highly visible "Bleeding Control Kit" lettering for fast identification and convenient Instructions for Use on the interior door panel

SPECIFICATIONS

H 15.5 in. x W 22.5 in. X D 4 in.

Empty Cabinet Weight: 9.365 lb

Low Pro Cabinet with ID Stickers:

Item# ZZ-1199

All Intermediate/Advanced Public Access Bleeding Control Kits and Stations are Texas House Bill 496 Compliant



8 Pack Vacuum Sealed Station

8-PACK LOW PROFILE METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	17.49 lb	80-1083
Intermediate	8	Vacuum	18.19 lb	80-1084
Advanced BCD	8	Vacuum	18.69 lb	80-1085
Basic	8	Nylon	17.49 lb	80-1080
Intermediate	8	Nylon	18.19 lb	80-1081
Advanced BCD	8	Nylon	18.69 lb	80-1082

Low-Pro CABINET features

- Low profile metal cabinet with 4 in. wall clearance
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Provides convenient slimline storage for medical equipment
- Magnetic closure provides quick easy access with one pull to open - No twist and pull

Alarm (#ZZ-1075) is also available for the Low Profile Metal Cabinet

Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.



*6 Packs
contain one
(1) QuikLitter™
Lite.
8 Packs
contain one
(1) Responder
QuikLitter™

6 Pack Station

Semi-Recessed Metal WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or Vacuum-Sealed pouches stored within an easy to-open nylon carrier and one (1) Responder QuikLitter™ (6 Packs contain QuikLitter™ Lite) for moving injured patients*
- Outside Dimensions: H 16 in. x L 17 in. x D 3.25 in.
Inside Dimensions: H 14 in. x W 15 in. x D 7.25 in.
Empty Case Wt: 10 lb (Empty Case Item# ZZ-0921)

8-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	16.8 lb	80-0861
Intermediate	8	Vacuum	20.6 lb	80-0862
Adv.	8	Vacuum	21 lb	80-0864
Adv. BCD	8	Vacuum	21 lb	80-0863
Adv. w. CELOX	8	Vacuum	21 lb	80-1578
Basic	8	Nylon	21.5 lb	80-0865
Intermediate	8	Nylon	22.3 lb	80-0866
Adv.	8	Nylon	21.5 lb	80-0868
Adv. BCD	8	Nylon	21 lb	80-0867
Adv. w/ CELOX	8	Nylon	21 lb	80-1579

6-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	6	Vacuum	18.7 lb	80-0577
Intermediate	6	Vacuum	19.5 lb	80-0578
Adv.	6	Vacuum	19.5 lb	80-0579
Adv. BCD	6	Vacuum	19.5 lb	80-0580
Adv. w/ CELOX	6	Vacuum	19.5 lb	80-0577


*All Multi-Pack Wall Stations feature:
Highly visible "Bleeding Control Kit" letters for fast ID
Instructions for Use posted on the interior back panel*



8 Pack Station

METAL CASE features

- Mounts permanently inside the wall with a low profile 1 in. trim
- Aluminum cabinet with a white finish and clear plexiglass window
- Stainless steel continuous hinge & chrome plated handle designed to open with one quick pull



Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.

Accessories | Public Access Bleeding Control

3-WAY SIGN Bleeding Control



SPECIFICATIONS

H 7 in. x W 16.5 in. x D 6.75 in.
Weight: 0.35 oz

- Highly visible for locating lifesaving bleeding control kits immediately
- 3 way "V" shape projects from wall for easy identification from multiple directions

3-WAY SIGN

ITEM

ZZ-1117

Cabinet Security ALARMS

- Alarm compatible with Clear Polycarbonate and Low Profile Metal Cabinets
- Loud 120 decibel alarm sounds when magnet separates from base more than 0.5 in.
- Easy and convenient to install
- Requires no wiring
- Mounting hardware included
- Designed for indoor use only



SPECIFICATIONS

W 1.63 in. x H 2.41 in. x D 0.72 in.
Weight: 6.4 oz each


CABINET SECURITY ALARMS

ITEM#	UNIT
ZZ-1075	Pack of 2

Bleeding Control Training Aids

C-A-T® TRAINER

Combat
Application
Tourniquet



- Fully functioning C-A-T® tourniquet colored blue to identify it as a training device
- Enables instructors and trainees hands-on training & tourniquet application practice
- Providing a distinct visual means to distinguish training devices from personal/operational life-saving devices

C-A-T® TRAINER

ITEM

30-0033

SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in.
Weight: 2.7 oz

HYFIN® TRAINER

Twin Pack



- Unique adhesive gel formulation does not damage mannequins or leave excessive residue behind
- Packaged with two training chest seals in one unique perforated package allowing students to open only one seal at a time

SPECIFICATIONS

Packaged: L 7.5 in. x W 4.5 in. x D 0.25 in.
Chest Seal: H 6 in. x W 6 in. x D 0.1 in.
Weight: 2.5 oz
U.S. Patent: 7,504,549

HYFIN® TRAINER

ITEM

10-0032

WOUND MANAGEMENT SIMULATOR

- Rugged part-task trainer for learning essential life-saving skills of tourniquet application, wound packing and applying pressure bandages.
- Multiple wounding patterns to enhance training objectives
- Simulates severe laceration, penetrating trauma (entrance and exit wounds), along with the ability to apply tourniquets
- Proprietary silicon that is rugged and durable and will not break down due to environmental exposure
- Can be used in the classroom or in field training scenarios

SPECIFICATIONS

L 15.38 in. x Diameter 6 in. | Weight: 10.8 lb

WOUND MGMT SIMULATOR

SKIN TONE	ITEM #
Light	93-0069
Deep	93-0410



MAKE THE **RIGHT CHOICE!**

THE **C-A-T®** TOURNIQUET IS THE
**FASTEST, SAFEST, MOST
EFFECTIVE** PREHOSPITAL FIELD
TOURNIQUET IN THE **WORLD!**



1 **Single Routing Buckle**

- Allows for fast application which decreases blood loss
- Simplifies effective slack removal

Dimensions:

- Device: L 37.5 in. | Weight: 2.7 oz
- Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.

Free-moving Internal Band 7

- Distributes pressure across TQ width, reducing pressure spikes on the TQ band edges as seen with competitor devices
- Minimizes friction during application



2 **Windlass Strap**

- NEW welding process features strong adherence of strap to clip (Released 2025)
- Secures windlass rod and prevents dislodgment



3 **Windlass Clip**

- Rapid windlass lock from either side (ambidextrous)
- Durable construction keeps windlass secure



4 **Windlass Rod**

- Large diameter for strength
- Aggressive ribbing for slip-free grip
- Twist for additional tightening of the tourniquet to stop bleeding



5 **Stabilization Plate**

- Reinforced beveled stabilization plate minimizes discomfort and provides stable pressure under windlass

6 **Red Tip Technology®** RED TIP

- Features NAR's signature of premium quality: red elliptical tip
- Rounded tip easily feeds through buckle

C-A-T TOURNIQUET

COLOR	STYLE	ITEM #
Black	Tactical	30-0001
Orange	Rescue	30-0023
Blue	Trainer	30-0033

SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 2.7 oz.



Find **C-A-T®**
videos, info,
and more at
NARescue.com

DON'T RISK LIVES WITH COUNTERFEITS!

- The **Combat Application Tourniquet®** has **saved more lives** than **all other** prehospital tourniquets **in the world** combined
- **Meets & exceeds ANSI/OSHA requirements** for self-applied, arterial **tourniquets** with a windlass or ratchet
- **Scientifically proven** 100% effective & trusted by the **U.S. Military**
- **Made in the U.S.A.**



CHOOSE THE GENUINE C-A-T® TOURNIQUET

NAR® WOUND PACKING GAUZE

ITEM

30-0054

NAR® WOUND PACKING GAUZE

- Ultra compressed, soft, white 1-ply gauze
- Includes embedded X-Ray detectable line
- Z-fold packaging for controlled application as a complete unit or in a continuous linear feed
- Used as a packing gauze for wounds, as a backing gauze for hemostatic agents, with a compression bandage to control hemorrhage, or for minor wound bandaging

SPECIFICATIONS

Gauze: W 4.5 in. x L 4.1 yd
Packaged: L 3 in. x W 2 in. x D 1 in.
Weight: 12 oz

NAR® Responder COMPRESSED GAUZE

- The compressed, rolled gauze is ideal for wound bandaging and wrapping of injuries
- Unique crinkle weave ensures excellent fluid absorption
- Ultra-compact size that expands to cover larger wound areas
- Used as backing gauze for hemostatic agents, in conjunction with a compression bandage to control hemorrhage

NAR® RESPONDER COMPRESSED GAUZE

ITEM

30-0182

SPECIFICATIONS

Gauze: W 4.5 in. x L 4.1 yd
Packaged: L 3 in. x W 2.5 in. x D 1 in.
Weight: 12 oz



CELOX™ FIRST AID Z-FOLD HEMOSTATIC DRESSING

- Z-Fold Configuration
- Fast-acting technology
- Advanced hemostatic agent combined with the simplicity of gauze
- Easy to apply on superficial wounds, minor cuts, and abrasions

CELOX™ First Aid

ITEM

30-0300

SPECIFICATIONS

Gauze: W 3 in. x L 5 ft.
Packaged: L 5.5 in. x W 4.5 in.
Weight: 0.1 lbs

QUIKCLOT® BLEEDING CONTROL DRESSING

- Advanced hemostatic agent combined with the simplicity of gauze
- Utilizes the same inorganic component as QuikClot Combat Gauze® used by the US Military

SPECIFICATIONS

Z-fold Gauze L 4 yd x W 3 in.
Pkg: H 5 in. x W 5.5 in. x D 0.625 in.
Weight: 0.8 oz

Roll Gauze L 4 ft x W 3 in.
Pkg: H 5 in. x W 3.5 in. x D 0.5 in.
Weight: 0.5 oz

QUIKCLOT BCD

STYLE	ITEM #
Z-Fold	30-0161
Roll	30-0164



NAR® RESPONDER ETD® Emergency Trauma Dressing

- Sterile pressure bandages that deliver extremely effective pressure and excellent absorption
- Suitable for applying immediate direct pressure
- Equipped with a non-adherent pad
- Easy-to-use hook & loop securing device (30-0179 also features C-Clip securing device)
- Simple application with the Quick grip roll control system which prevents the dressing from unrolling
- Rugged vacuum-sealed, low-cube packaging with easy-open tear notches



RESPONDER EMERGENCY TRAUMA DRESSINGS

STYLE	BDG.	PAD	WGT	PACKAGE	ITEM #
4in.Flat	4 x 60	4x6	2 oz	L 4.25 in. x W 2.75 in. x D 0.88 in.	30-0171
4in.(Mini)	4 x 42	4x6	2.4 oz	L 5 in. x W 3.5 in. x D 0.75 in.	30-0179
4in.Roll	4 x 60	4x6	2 oz	L 4.5 in. x W 2.5 in. x D 1.75 in.	30-0167
6in. Flat	6 x 56	6x6	3.31 oz	L 6.25 in. x W 3.75 in. x D 0.88 in.	30-0172
6in. Roll	6 x 60	6x6	2.9 oz	L 6 in. x W 2.5 in. x D 2.1 in.	30-0168

RESPONDER Nitrile Gloves

- High performance personal protection compliant with the EN ISO 374-5:2016 Standard for Virus Protection
- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Latex-free and Powder Free

SPECIFICATIONS

Bag (25 pair): L 4.75 in. x W 6 in. x D 4 in.
Weight: 1 lb 1 oz

RESPONDER NITRILE GLOVES

SIZE	PKG QTY	ITEM #
S	Bag 25 pair	ZZ-1059
M	Bag 25 pair	ZZ-1060
L	Bag 25 pair	ZZ-1043
XL	Bag 25 pair	ZZ-1061



HYFIN® VENT COMPACT Twin Pack

- Smaller packaging, approximately 25% smaller than our standard HyFin® Vent
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Adheres to sweaty or hairy chests
- Each chest seal also includes a gauze pad to wipe the wound surface prior to application.
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary


SPECIFICATIONS

Device (excludes tab): 4.75 in. x 4.75 in.
Weight: 1.55 oz
Packaged: H 3.88 in x W 6.6 in x D 0.13 in.
U.S. Patent: 7,504,549

HYFIN VENT COMPACT

ITEM #

10-0042



Accessories | Hypothermia Prevention | Evacuation

NAR® SURVIVAL BLANKET

- Lightweight, durable space blanket made of a combination of mylar and polyethylene, perfect for regulating body temperature
- Reflective silver color helps to put radiant heat back into the body
 - Wind and waterproof, further aiding in retaining warmth
 - 56 in. x 84 in. wrap for full body coverage



SPECIFICATIONS

Packaged: L 5 in. x W 3 in.
Blanket: L 84 in. x W 56 in.
Weight: 1.84 oz.

SURVIVAL BLANKET

ITEM #

NSN #

ZZ-0038

4240-01-697-8869

NAR® Responder QUIKLITTER

- Multi-purpose, low cost, disposable evacuation platform
- Woven polyethylene
- Latex free & water-resistant
- Red nylon webbing & stitching
- Ten (10) carry handles

NAR® RESPONDER QUIKLITTER

ITEM #

60-0080

SPECIFICATIONS

Litter Bed: L 74 in. x W 28 in.
Weight: 13.3 oz
Payload 500 lb.



LEATHERMAN® RAPTOR EMERGENCY MEDICAL MULTI-TOOL



High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears
2. Strap Cutter
3. Oxygen Tank Wrench
4. Ring Cutter
5. Carbide Glass Breaker
6. 5 cm Ruler
7. Holster with Clip

SPECS

Shear Length: 1.9 in.
Closed Length: 5 in.
Weight: 5.8 oz

HANDLE COLORS:



RAPTOR MULTI-TOOL

HANDLE	ITEM #
Black	70-0014
Blue/Black	70-0264
Red/Black	70-0266
Yellow/Blue	70-0326

Heavy Duty Trauma XSHEAR®

- Hardened stainless steel blades that are twice as thick as most trauma shears
- Patented blunt tip for gentler edge near the skin
- Soft touch inner handles

XSHEAR®

HANDLE	BLADE	ITEM #
Black	Black	70-0290
Blue	Silver	70-0291
Orange	Black	70-0292
Pink	Black	70-0293
Red	Black	70-0294
Coyote	Black	70-0415
Green	Silver	70-0424



SPECIFICATIONS

L7.5 in x 3.5 in. x 0.118 in.
Cutting area length: 2 in.
Weight: 4 oz



Rigid XSHEAR® HOLSTER

- Custom designed for the XSHEAR® Heavy Duty Trauma Shears
- MOLLE Compatible adjustable clip
- Flex lock design securely holds shears in place with a positive click
- Universal design allows shears to be inserted either direction
- Rugged injection molded ABS Plastic

RIGID XSHEAR® HOLSTER

ITEM

70-0295

SPECIFICATIONS

L 5 in. x W 3.75 in. x D 1.75 in. | Weight: 3.2 oz

Trauma kit WALL MOUNT HOOK

- Designed to store Class A or B Trauma and First Aid Kits (sold separately) in areas where rapid access is needed
- 5/16 in. thick, heavy duty, zinc plated steel hook
- Easily mounts to any wall

SPECIFICATIONS

Hook Dimensions:
H 2.87 in. x W 3 in. x D .75 in
Mounting Plate Dimensions:
H 3 in. x W 5 in. x D 0.12 in.
Weight: 1 lb.

WALL MOUNT

ITEM

ZZ-1819



SALTY BRITCHES® SKIN BARRIER

- Offers soothing, long-lasting barrier against skin irritation from windburn, friction, and chafing
- Waterproof
- Won't rub off

SPECIFICATIONS

Tubes

Packaged: W 2.5 in. x L 4 in. | Wt: 0.15 lbs

Single Use (10 pk)

Packaged: W 3 in. x L 5 in. | Wt: 0.10 lbs

SALTY BRITHCES®

OPTION	ITEM #
2 oz. Skin Barrier	75-0014
2 oz. Winter Relief	75-0015
0.1 oz. Single Use*	75-0017



*Pack of 10



COMMUNITY
PREPAREDNESS



Redefining FIRST AID

NORTH AMERICAN RESCUE

NARescue.com • 888.689.6277

TRAUMA
AND FIRST AID KIT

Class A
Bleeding Control Dressing

80-1045

CLASS A



TRAUMA &
FIRST AID KIT

ANSI/OSHA COMPLIANT

ANSI Z391.1-2015

23-0863

North American Rescue, LLC • 355 Tenthred Ct. • Denver, CO 80602 • USA

www.NARescue.com • 888-689-6277

REV0424



PELICAN
1400 CASE

PRODUCTS WITH A MISSION™

BE THE DIFFERENCE™

www.NARescue.com

Sign up to get updates and special offers when you...

#JOIN THE MISSION



PUBLIC ACCESS BLEEDING CONTROL KITS | ANSI WORKPLACE FIRST AID TRAUMA KITS

Volume 4 • Number 3 • REV041625