

TOURNIQUETS | FIRST AID KITS | TRAUMA KITS | RESCUE TOOLS | LITTERS

LE

LAW ENFORCEMENT
Catalog



Since 1996, North American Rescue® has focused on **increasing survivability** on both the battlefield and the homefront. With our innovative medical/rescue equipment, NAR® continues to lead the way by delivering the **tailored, mission-critical healthcare solutions** needed by today's first responders to counter the rapidly evolving threats of the 21st century. See what a lifesaving difference our combat proven make in the prehospital environment!



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A MISSION®



BE THE DIFFERENCE™

www.NARescue.com

CONTENTS

TOURNIQUETS	3-4
FIRST AID KITS	4 -12
TRAUMA KITS	13-16
EXTRAS	17

DISCOUNT PRICING

Contact Customer Service to have your account upgraded to discount pricing. Shop your pricing on our full catalog of products, available online 24-7 at www.NARescue.com

888.689.6277





CUSTOMIZABLE SOLUTIONS

BUILT TO MEET YOUR REQUIREMENTS

If you haven't found exactly what you need, our team of **subject matter experts** at North American Rescue® are standing by to guide you to the **ideal lifesaving solution** tailored to meet all your specifications.

Whether it's the essential gear needed for individuals or customized kits for multiple casualties, a **North American Rescue® Specialty Kit** is your answer to the precise aid kit you and your organization need to **increase survivability** should injury occur.



Call or Email us **TODAY**
to discover all your
options for the tailored
lifesaving solution
YOU can live with!



Explore all the possibilities at
NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

High Quality Tactical Nylon - OR - Vacuum Bag

Massive Hemorrhage | Tourniquets

C-A-T® Combat Application Tourniquet Gen 7

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC



online video available

< scan to view



**The FASTEST, SAFEST,
Most EFFECTIVE
Prehospital Field Tourniquet!**



COMBAT APPLICATION TOURNIQUET

SPECIFICATIONS

Packaged:
H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 3.04 oz.

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972

RIGID C-A-T® TQ CASE *with Shirt Shield*

- Allows carriage of a Gen 7 C-A-T® Tourniquet on your duty/patrol/pistol belt (Tourniquets sold separately)
- Shirt Shield protects uniform shirts from undue friction wear
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally, or to a Safariland Drop-leg Shroud
- Injection Molded Nylon Case stands up to rugged military usage
- Cross medical identifier patch included

RIGID C-A-T TQ CASE

COLOR	ITEM #
Black	80-0944
Black, Basket Weave	80-0962
Black, Polished	80-0963
Coyote	80-0945
Ranger Green	80-0946
Navy	80-1161

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in.
Weight: 5.36 oz

Explore other TQ carrier
options available at
NARescue.com

RIGID C-A-T® INSIDE-THE-WAISTBAND TOURNIQUET CASE

- Keeps your C-A-T® Tourniquet close at hand for rapid deployment
- Holds the C-A-T® Tourniquet (Gen 7) utilizing friction (Tourniquet sold separately)
- Designed with built-in Shirt Shield that creates a barrier between your body, shirt, and the exposed portion of your tourniquet
- Injection Molded Case to stand up to rugged LE/MIL usage
- Drainage hole at bottom of case
- Offers a range of adjustment
- Available Mounting Options – 1.5" and 1.75" IWB Overhook (included) designed, manufactured, and assembled in the USA.

RIGID CAT TQ CASE

COLOR	ITEM #
Black	80-1488
Gray	80-1489

SPECIFICATIONS

Packaged: L 5 in. x W 3.25 in. x D 1.625 in.
Weight: 3.0 oz

LEATHER C-A-T® HOLDERS



- Designed to meet the specific duty gear needs of law enforcement officers - for use on daily duty belt to fit within organizational uniform standards
- Securely fastens to your duty belt with easy one-handed access to allow you to take control of extreme medical situations in seconds
- Constructed of black full grain leather in different finishes, attachment options, and snap colors

- High quality snaps ensure secure placement on your duty belt and easy accessibility to your tourniquet
- Ideal for the C-A-T® Gen 7 but also accommodates other popular windlass TQs
- Attachment options: Clip, Loop or Hook & Loop Strap attachment on select finishes and snap colors



◀ **with Strap Attachment**
with hook and loop allows direct attachment to a radio strap.

Rigid Case MOLLE ATTACHMENT

ITEM

ZZ-1696

SPECIFICATIONS

H 4 in. x W 1 in. x D 0.5 in.
Weight: 0.4 oz.



SPECIFICATIONS

L 5.75 in. x W 2 in. x D 2 in.
Weight: 3.2 oz.

UNIFORM DUTY BELT LEATHER TQ HOLDER

FINISH	ATTACHMENT	SNAP	ITEM #
High Gloss	Clip	Silver	30-0194
Full Grain	Clip	Silver	30-0193
Basket Weave	Loop	Silver	30-0195
Full Grain	Loop	Silver	30-0197
Basket Weave	Strap Attach	Silver	30-0228
High Gloss	Strap Attach	Silver	30-0229

Nylon C-A-T® HOLDER

- Tailored specifically for your C-A-T®, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- Easy-open pull elastic tab with tourniquet ID patch - "TQ"

NYLON TQ HOLDER

COLOR	ITEM #
Black	30-0057
Coyote	30-0044
Multi-Cam	30-0059
Red	30-0128
Ranger Green	30-0208
RealTree®	30-0207
OD Green	30-0053



SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 2.17 oz

LASER CUT C-A-T® HOLDER

- Constructed with laser cut laminated Cordura® nylon
- Easy Pull Elastic Cover
- PALS/MOLLE-style connections
- Accommodates GEN 6 and GEN 7 C-A-T® Tourniquets
- Weighs less than traditional TQ holders

SPECIFICATIONS

L 7.25 in. x W 2.75 in. x D 1.5 in.
Weight: 1.29 oz

LASER CUT TQ HOLDER

COLOR	ITEM #
Black	30-0297
Coyote	30-0298
Ranger Green	30-0299
OCP	30-0306
Red	30-0551



IFAK Individual First Aid Kits



Individual Patrol Officer Kit **IPOK**

- Individual hemorrhage control kit
- Compact, vacuum-sealed & packaged to fit into a pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance

SPECIFICATIONS

H 4 in. x W 6.5 in. x D 2.75 in.

Weight: 7.4 oz



IPOK:



IPOK - INDIVIDUAL PATROL OFFICER KIT

HEMOSTATIC/GAUZE	ITEM#
Wound Packing Gauze	80-0167
Combat Gauze, LE	80-0168

ETAK

Enhanced Trauma Aid Kit

Designed to pair with a separate C-A-T® Tourniquet & holder (not included) for a great option to add more capability to previously owned kit

- Compact IFAK platform accessible from a belt, duty belt, assault vest/body armor or any equipment with PALS/MOLLE attachment points.
- Featuring multiple interior elastic loops keep your kit organized and secure
- Exterior tear away panel is easily removable for added mobility and access to your kit
- MED patch included

Every ETAK includes these core items:




PLUS+

BASIC also includes:



with **COMBAT GAUZE** also includes:

ADVANCED also includes:



Bag only option
also available at
NARescue.com

ETAK ENHANCED TRAUMA AID KIT

COLOR	LEVEL	ITEM #
Black	Basic	80-1168
Coyote	Basic	80-1169
Ranger Green	Basic	80-1171
MultiCam	Basic	80-1170
Black	Combat Gauze LE	80-1172
Coyote	Combat Gauze LE	80-1173
Ranger Green	Combat Gauze LE	80-1175
MultiCam	Combat Gauze LE	80-1174
Black	Advanced	80-1176
Coyote	Advanced	80-1177
Ranger Green	Advanced	80-1179
MultiCam	Advanced	80-1178

SPECIFICATIONS

H 5.75 in. x W 4 in. x D 1.75 in.

Weight: 0.65 lb

IPOK HOLDER

- Protects your medical kit from the elements
- Enclosure can be secured by Fastex snap buckle, hook-and-loop or both



IPOK HOLDER

COLOR	ITEM#
Black	80-0205T
Coyote	80-0207
OD Green	80-0254T
Red	80-0358T

SPECIFICATIONS

H 6.5 in. x W 6 in. x D 3.5 in.
Weight: 4.1 oz

VISOR MED POUCH

- Securely mounts to a vehicle visor
- Orange webbing handle to grab for easy deployment
- Pouch with easy-open buckle



VISOR MED POUCH

COLOR	ITEM#
Black	80-0269
Coyote	80-0270
Red	80-0360

SPECIFICATIONS

H 6.5 in. x W 5.5 in. x D 3 in.
Weight: 5.4 oz

HEADREST IFAK

- Designed to provide personnel with an easily accessible, quick to deploy medical pouch that conveniently mounts to a vehicle's headrest
- Fully loaded with emergency medical equipment for point-of-wounding treatment until definitive care is available
- Durable nylon construction features multiple interior elastic loops allowing for the contents to be organized and secured
- Easy to grab side handles allow for easy deployment
- IFAK pouch only is also available
- The pouch design allows for the contents of an individual first aid kit or vacuum sealed kit



SPECIFICATIONS

H 5.25 in. x
W 7.25 in. x D 4.25 in.
Weight: 1.63 lb



HEADREST IFAK



HEADREST IFAK

COLOR	OPTIONS	ITEM #
Black	Basic	80-1091
Coyote	Basic	80-1150
Black	Combat Gauze LE	80-1151
Coyote	Combat Gauze LE	80-1152
Black	Bag Only	80-1090T
Coyote	Bag Only	80-1149T

IFAK Individual First Aid Kits

M-FAK Mini First Aid Kit - LE

- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

Bag only
option also
available at
NARescue.com

SPECS

H 6 in. x W 3 in.
x D 3 in.
Weight: 13 oz



LE M-FAK | ROO® M-FAK

BASIC or BASIC with Combat Gauze:



LE M-FAK | ROO® M-FAK ADVANCED:



Compact Officer Response Emergency Kit **CORE KIT**

Designed to go where you go, this compact medical kit will fit in a BDU pocket, vest pouch or Go-Bag and provides personnel with essential life-saving bleeding control equipment that can be kept on their person. The CORE Kit is packaged in a compact, rugged, clear vacuum-sealed bag with easy-to-open tear notches.

- Compact design with the essential life-saving equipment needed to address moderate and severe traumatic bleeding
- Contains the #1 CoTCCC Recommended Prehospital Tourniquet and Hemostatic Gauze
- The Best Tourniquet paired with the Best Hemostatic Gauze and an effective pressure bandage at the Best Price
- Meets the Tactical Emergency Casualty Care (TECC) Guidelines for hemorrhage control

CORE KIT:



SPECS

H 7 in. x W 3.5 in. x D 2.25 in. | Weight: 6.4 oz

CORE KIT

ITEM

85-0834

ROO® M-FAK Mini First Aid Kit

Rapid one handed access to your tourniquet without ever having to open the full kit!

Contains the same packout options as the M-FAK:
Basic, Basic with Combat Gauze® LE and Advanced

ROO M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM #
Black	Basic	80-1049
Coyote	Basic	80-1050
OD Green	Basic	80-1051
Ranger Green	Basic	80-1052
Black	Combat Gauze LE	80-1053
Coyote	Combat Gauze LE	80-1054
OD Green	Combat Gauze LE	80-1055
Ranger Green	Combat Gauze LE	80-1056
Black	Advanced	80-1057
Coyote	Advanced	80-1058
OD Green	Advanced	80-1059
Ranger Green	Advanced	80-1060



L-FAK Lumbar First Aid Kit

- Durable 500-denier exterior carrier
- Compact design to reduce space needed for life-saving equipment
- Designed for wear in the small of the back
- Multiple interior elastic loops for securing gear
- Rapid access to medical gear
- C-A-T® tourniquet sold separately




Every L-FAK includes these core items:



PLUS+ BASIC
also includes:



COMBAT GAUZE®
also includes:



ADVANCED also includes:



SPECIFICATIONS

Pkgd - H 4.75 in. x W 5.75 in.
Weight: 0.88 lbs
Adv. Weight: 0.94 lbs

LUMBAR FIRST AID KIT

COLOR	LEVEL	ITEM #
Black	Basic	85-3627
Coyote	Basic	85-3633
Ranger	Basic	85-3630
Black	w/ Combat Gauze™	85-3628
Coyote	w/ Combat Gauze™	85-3634
Ranger	w/ Combat Gauze™	85-3631
Black	Advanced	85-3629
Coyote	Advanced	85-3635
Ranger	Advanced	85-3632

Bag only option also available at NARescue.com

IFAK Individual First Aid Kits



EDC ANKLE TRAUMA KIT

OPTIONS	ITEM #
Wound Packing	80-1157
BCD	80-1158
Carrier Only	30-0294T

EDC ANKLE TRAUMA KIT

- Lightweight, low profile, versatile carrier for essential trauma and EDC equipment with quick, single-handed, ambidextrous access
- Separate concealed compartment for essential EDC items or ID
- Designed to accommodate a C-A-T® Tourniquet and gauze/bandage
- Advanced materials that provide user comfort and long-term wearability
- Individual hook and loop pouches for each carried component
- One size fits most with convenient hook and loop adjustability
- Universal magazine pocket

EVERY DAY CARRY ANKLE TRAUMA KIT



SPECIFICATIONS

L 17.75 in. x H 4 in. x D 0.25 in.
Weight: 0.657 lb

REFLEX™ SYSTEM ROLL & CARRIER

- Primary construction of heavy-duty nylon laminate allows rapid deployment of medical supplies
- Features High Abrasion-Resistant Neoprene Kevlar® (HANK™) composite grab handles for long-term durability
- Can mount vertically or horizontally on belts 1.5 in.- 2.25 in., and horizontally on MOLLE
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack

- "Mini" MOLLE allows for precise fit on belts as narrow as 1.5 in. or two-row MOLLE panels
- Features loop panels on med roll and carrier for label patches (Medical Cross Patch included)
- Two-piece system designed for ambidextrous access
- Medical contents sold separately

REFLEX SYSTEM ROLL & CARRIER

COLOR	ITEM #
Black	80-1117
LE Blue	80-1118
Coyote	80-1119
Gray	80-1120
Multi-Cam Blk	80-1121
Multi-Cam	80-1122
OD Green	80-1123



SPECIFICATIONS

L 8 in. x W 3 in. x H 5.5 in.
Weight: .58 lb




REFLEX SYSTEM IFAK

COLOR	LEVEL	WGT	ITEM #
Black	Basic	1.58 lb	80-1127
LE Blue	Basic	1.58 lb	80-1128
Coyote	Basic	1.58 lb	80-1129
Gray	Basic	1.58 lb	80-1130
Multi-Cam Black	Basic	1.58 lb	80-1131
Multi-Cam	Basic	1.58 lb	80-1132
OD Green	Basic	1.58 lb	80-1133
Black	Combat Gauze	1.6 lb	80-1134
LE Blue	Combat Gauze	1.6 lb	80-1135
Coyote	Combat Gauze	1.6 lb	80-1136
Gray	Combat Gauze	1.6 lb	80-1137
Multi-Cam Black	Combat Gauze	1.6 lb	80-1138
Multi-Cam	Combat Gauze	1.6 lb	80-1139
OD Green	Combat Gauze	1.6 lb	80-1140
Black	Advanced	1.61 lb	80-1141
LE Blue	Advanced	1.61 lb	80-1142
Coyote	Advanced	1.61 lb	80-1143
Gray	Advanced	1.61 lb	80-1144
Multi-Cam Black	Advanced	1.61 lb	80-1145
Multi-Cam	Advanced	1.61 lb	80-1146
OD Green	Advanced	1.61 lb	80-1147

REFLEX™ IFAK SYSTEM

- Developed with direct input from active-duty medical personnel, the ReFlex™ System IFAK comes fully loaded with essential medical equipment to treat life threatening injuries in the field of operation
- Organized to put all medical supplies within easy reach
- Heavy-duty nylon laminate allows for rapid deployment
- HANK™ composite grab handles for long-term durability
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack
- "Mini" MOLLE allows for precise vertical/horizontal mount on belts 1.5 in. - 2.25 in.
- Mounts horizontally on MOLLE



SPECIFICATIONS

L 8 in. x W 3 in. x H 5.5 in. | Weight: 1.6 lb

Every REFLEX SYSTEM IFAK includes these core items:



PLUS+ BASIC

also includes:





COMBAT GAUZE

also includes:



ADVANCED also includes:



ReFlex™ System VEHICLE MOUNT

- Attaches ReFlex™ Med Roll to the back of a vehicle headrest
- All-nylon construction with built-in elastic tourniquet loop



REFLEX MOUNT

COLOR	ITEM #
Coyote	ZZ-1369
Gray	ZZ-1370

SPECIFICATIONS

L 5.5 in. x W 0.125 in. x H 6 in.

Weight: 0.21 lb



RIG SERIES EAGLE IFAK

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access

RIG SERIES EAGLE IFAK

COLOR	LEVEL	ITEM #
● Black	Basic	80-0509
● Coyote	Basic	80-0510
● Multi-Cam	Basic	80-0511
● Ranger Green	Basic	80-0991
● Black	Combat Gauze	80-0517
● Coyote	Combat Gauze	80-0518
● Multi-Cam	Combat Gauze	80-0519
● Ranger Green	Combat Gauze	80-0993
● Black	Advanced	80-0513
● Coyote	Advanced	80-0514
● Multi-Cam	Advanced	80-0515
● Ranger Green	Advanced	80-0992

RX ONLY

Bag only option
also available at
NARescue.com

RIG SERIES EAGLE BASIC:




RIG SERIES EAGLE with COMBAT GAUZE LE:



SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
Kit Packed Weight: 118 lb
Empty Bag Weight: 7 oz

RIG SERIES EAGLE ADVANCED:



Individual First Aid Kits IFAK

T.O.R.K. TACTICAL OPERATOR RESPONSE KIT

COLOR	LEVEL/HEMOSTATIC	ITEM #
Black	Basic	85-0008
Coyote	Basic	85-0009
Digital Camo	Basic	85-2300
Multi-Cam	Basic	85-0012
Ranger Green	Basic	85-2301
Red	Basic	85-1805
MTC Alpine	Basic	80-1633
Black	Basic w/Combat Gauze	80-1097
Coyote	Basic w/Combat Gauze	80-1098
Digital Camo	Basic w/Combat Gauze	80-1099
Multi-Cam	Basic w/Combat Gauze	80-1100
Ranger Green	Basic w/Combat Gauze	80-1102
Red	Basic w/Combat Gauze	80-1108
MTC Alpine	Basic w/Combat Gauze	80-1635
Black	Advanced w/Combat Gauze	85-0048
Coyote	Advanced w/Combat Gauze	85-0049
Digital Camo	Advanced w/Combat Gauze	85-0050
Multi-Cam	Advanced w/Combat Gauze	85-0052
Ranger Green	Advanced w/Combat Gauze	85-2302
Red	Advanced w/Combat Gauze	85-1806
MTC Alpine	Advanced w/Combat Gauze	80-1636
7 options online	Advanced w/ChitoGauze	NARescue.com

T.O.R.K. Tactical Operator Response Kit

- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid injuries caused by penetrating or blast trauma

T.O.R.K. ReSUPPLY KITS*

CONTENT OPTIONS	ITEM #
Basic	80-1037
Basic Combat Gauze	80-1103
Advanced Combat Gauze	80-1039
Advanced ChitoGauze	80-1038

*ReSupply Kits include contents only. Nylon bag is sold separately.



T.O.R.K. BASIC/BASIC with Combat Gauze®:



SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz

T.O.R.K. ADVANCED:



Bag only options also available at NARescue.com



RAID KIT V2

Releasable Abdominally Integrated Drop kit



Individual First Aid Kits IFAK



RAID KIT V2

COLOR	ITEM#	OPTION	NSN#
BLK	80-1545	Basic	----
COY	80-1546	Basic	----
MTC	80-1547	Basic	----
RGR	80-1548	Basic	----
BLK	80-1549	Basic w/ CG	----
COY	80-1550	Basic w/ CG	----
MTC	80-1551	Basic w/ CG	----
RGR	80-1552	Basic w/ CG	----
BLK	80-1625	Basic w/ CELOX™	----
COY	80-1626	Basic w/ CELOX™	----
MTC	80-1627	Basic w/ CELOX™	----
RGR	80-1628	Basic w/ CELOX™	----
BLK	80-1553	Advanced	----
COY	80-1554	Advanced	----
MTC	80-1555	Advanced	----
RGR	80-1556	Advanced	----
BLK	80-1629	Adv. w/ CELOX™	----
COY	80-1630	Adv. w/ CELOX™	----
MTC	80-1631	Adv. w/ CELOX™	----
RGR	80-1632	Adv. w/ CELOX™	----

Provides several options for carrying essential life-saving equipment



Every RAID Kit V2 includes these core items:



- RAID V2 pouch's unique take on the industry standard hanger combined with the utility of a fanny pack allows user to decide how to wear it based on their needs during self-aid/buddy-aid
- Designed to meet the need for an individual first aid kit with multiple carry options
- First Spear Tube makes reattaching the pouch as easy as removal

PLUS+

BASIC
also includes:

OR

BASIC with COMBAT GAUZE®
also includes:

OR

BASIC with CELOX™ GAUZE
also includes:



OR

ADV. with COMBAT GAUZE®
also includes:

OR

ADV. with CELOX™ GAUZE
also includes:

AND



SPECIFICATIONS

H 8 in. x W 7 in. x D 2 in.
Weight: 3 lbs 3.2 oz.



MINI MEDIC KIT

- Compact, multi-purpose trauma aid bag with vital medical equipment for point of wounding treatment of multiple injuries
- Two external pockets with two external tourniquet holders
- One external pouch with additional drop-down pouch for use with breaching tools or storage for a hydration pack
- Expansion zipper to increase load capacity
- External MOLLE on sides and bottom of bag / Internal MOLLE rows on inside back of bag
- Three clear tear away pouches on the inside main compartment
- Two internal mesh pouches and one larger internal pouch on inner lid
- Internal tape roll strap

Every MINI MEDIC KIT includes these core items:



PLUS+ BASIC
also includes:

COMBAT GAUZE
also includes:

ADVANCED
also includes:



Bag only option
also available at
NARescue.com

Inside Back
Compartment for
stowing extra
equipment

Breaching
Equipment
Sold
Separately

MINI MEDIC KIT


COLOR	LEVEL	ITEM #
Black	Basic	80-0979
Coyote	Basic	80-0980
OD Green	Basic	80-0981
Ranger Green	Basic	80-1003
MultiCam	Basic	80-1231
Black	Combat Gauze	80-0985
Coyote	Combat Gauze	80-0986
OD Green	Combat Gauze	80-0987
Ranger Green	Combat Gauze	80-1005
MultiCam	Combat Gauze	80-1232
Black	Advanced	80-0982
Coyote	Advanced	80-0983
OD Green	Advanced	80-0984
Ranger Green	Advanced	80-1004
MultiCam	Advanced	80-1233

Patrol **VEHICLE TRAUMA KIT**

- Provide personnel with an easily accessible, quick to deploy medical kit from your Patrol Vehicle
- Modular storage system comes organized with an array of medical products ranging from basic first aid to more intermediate trauma care
- One side covers all your agency needs to administer basic first aid and community assistance
- Second side offers a first time combination of Basic First Aid and Basic Trauma medical items
- Exterior includes MOLLE-style system to accommodate extra equipment as needed



PATROL VEHICLE TRAUMA KIT:



PATROL VEHICLE TRAUMA

LEVEL	ITEM #
Basic	80-0542
Combat Gauze LE	80-0859

SPECIFICATIONS

H 7 in. x W 12 in. x D 5 in.
Weight: 6 lb



Trauma Kits

RANGE TRAUMA AID KIT *Hard Case*

- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches
- Labeled internal pouches for ease of use

RANGE TRAUMA AID

ITEM

85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.

Weight: 8 lb 12 oz



RANGE TRAUMA AID KIT HARD CASE:



School Resource Officer Kits

SRO CRISIS RESPONSE KIT

- Equip School Resource Officers with the essential life-saving tools to begin treating injured victims until medical personnel arrive
- Empowers an SRO (School Resource Officer) to administer self-aid in the event he or she is injured or provide immediate care to injured casualties during a school incident
- Multiple carry options

SRO CRISIS RESPONSE KIT:



SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.

Weight: 5 lb 16 oz

SRO CRISIS RESPONSE KIT

ITEM

85-0417




Bag only option
also available at
NARescue.com

I-ROK Individual Rugged Outdoor Kit

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding
- Storage case is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- Available in 3 hardcase colors and 2 sizes with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions




Every I-ROK contains the following items:



SPECIFICATIONS

H 7 in. x W 4.5 in.
x D 2.5 in.
Weight: 2.48 lb



I-ROK



COLOR	LEVEL	ITEM #
Black	Basic	82-0031
OD Green	Basic	82-0033
Orange	Basic	82-0034
Black	w/BCD	82-0035
OD Green	w/BCD	82-0037
Orange	w/BCD	82-0038
Black	w/CELOX	82-0066
OD Green	w/CELOX	82-0068
Orange	w/CELOX	82-0069

I-ROK XL

SPECIFICATIONS

H 8.9 in. x W 5.4 in. x D 3 in.
Weight: 3 lb

I-ROK XL includes these extra items:



I-ROK XL

COLOR	LEVEL	ITEM #
Black	Basic	82-0049
OD Green	Basic	82-0051
Orange	Basic	82-0052
Black	w/BCD	82-0053
OD Green	w/BCD	82-0055
Orange	w/BCD	82-0056
Black	w/CELOX	82-0070
OD Green	w/CELOX	82-0072
Orange	w/CELOX	82-0073



SPECS

Shear Length: 1.9 in.
Closed Length: 5 in.
Weight: 5.8 oz

Leatherman® RAPTOR® RESPONSE SHEARS

1. Stainless Steel folding shears
2. Integrated ring cutter
3. 5 cm ruler
4. Oxygen tank wrench
5. Lanyard connection in handle
6. Removable pocket clip

RAPTOR RESPONSE SHEARS

HANDLE	ITEM #
● Cement	70-0328
● Crimson	70-0329
● Navy	70-0330



SPECS

Length: Closed 4.3 in.
Open 7.25 in.
Shear Length: 1.9 in.
Width: 1.7 in | 4.3 cm
Overall Thickness: .85 in.
Weight: 5.5 oz

Benchmade® 7 HOOK PERSONAL SAFETY CUTTER

- Saves time – more than 1.5 minutes faster than traditional trauma shears when cutting 2 articles of clothing.
- Extensively tested and widely used by the US military for medical and egress purposes
- Maintains sharpness for hundreds of cuts (up to 500+ cuts through seatbelt test material)
- 440C cutting surface and opening with handle to provide a more compact carry
- Includes a MOLLE-compatible carrying sheath



HOOK SAFETY CUTTER

ITEM #	STYLE
ZZ-0652	7 Hook

SPECS

Length: 4.3 in.
Weight: 1.5 oz

NAR® QUIKLITTER LITE

- Perfect balance of size, weight, cost, and capability
- Features ten large carry handles with a load capacity tested to 750 lb
- Made of lightweight, heavy-duty nylon
- Vacuum-packed to conserve space and fit in a cargo pocket
- Latex-free & water-resistant

QUIKLITTER

SIZE	ITEM #
LITE*	60-0071

*Nylon Construction

SPECIFICATIONS

Payload 750 lb.
Litter Bed: L 74 in. x W 27 in.
Weight: 12 oz



UNO EMERGENCY SPLINT

- Modular Design
- Stepless Adjustment
- Small and Lightweight
- Compatible with TCCC Equipment
- Pack of 3 includes additional ankle strap for stabilizing the ankle and enables the device to provide traction



SPECIFICATIONS

Packaged - L 21.50 in. x W 10.55 in. x D 1.50 in.
Unpackaged - L 15.75 in. x W 2.17 in. x D 2.36 in.
Weight: 0.91 lbs (sgl)
2.27 lbs (3 pk)

UNO EMERGENCY SPLINT

ITEM #	OPTION
50-0045	Single
50-0046	Pack of 3



NORTH AMERICAN RESCUE®

Education and Training

North American Rescue Education and Training (NAR® Training)

has roots that run deep and wide in the world of pre-hospital medicine and technical rescue. Our team was built from the ground up on the vision and experience of seasoned military medics seeking to fill the training gaps commonly encountered in the operational environment.

Since then, we have continued to broaden our vision to include other invaluable first responder perspectives such as pre-hospital healthcare providers in civilian, law enforcement, emergency medical services and community public safety initiatives.

CORE Competencies

- Continuing Medical Education (CMEs)
- Prolonged Casualty Care
- Austere Medicine
- Technical Rescue
- Emergency Medical Services
- Critical Care Flight Medicine



ONLINE




INSTRUCTOR



HANDS ON



IMMERSIVE



Training Locations:

- NAR Training
- Training Partners & Affiliates

ACCREDITATIONS:



AUTHORIZED
TRAINING
SITE

AUTHORIZED
TRAINING
CENTER



American
Heart
Association®


Now SAME DAY Cards Issued At End Of Course!



Explore
training
opportunities
online

NAReducation.com ➤ **Phone:** 888.689.6277 ➤ **Email:** NARTraining@NARescue.com

CAT[®] COMBAT APPLICATION TOURNIQUET



**THE FASTEST,
SAFEST, MOST
EFFECTIVE
PREHOSPITAL FIELD
TOURNIQUET**

TOURNIQUETS | FIRST AID KITS | TRAUMA KITS | RESCUE TOOLS | LITTERS



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

35 Tedwall Court | Greer, SC 29650 | USA

PRODUCTS WITH A MISSION[®]



BE THE DIFFERENCE[®]

www.NARescue.com

REV010424 | Volume 3 | Number 1