NAR® ELITE LITTER™ - RESPONDER LITE

Evacuation Litter for Incident Tactical Extraction

Expanding on the success and feedback from the ELITE Litter™ - Responder and ELITE Litter™ - Responder XL Litters, North American Rescue's ELITE (Evacuation Litter for Incident Tactical Extraction) Litter - Responder LITE has added capability and reduced weight and size to be able to fit into a cargo pocket. The ELITE Litter™ - Responder LITE is an emergency carry litter that is designed for rapid deployment for the movement of casualties for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements. Constructed of a heavy-duty polypropylene, the NAR® ELITE Litter™ - Responder LITE is rated to 750 lb and is latex free. It features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points. Check out the entire series of NAR® ELITE Litter™ products for the right features based on your requirements.

FEATURES/SPECIFICATIONS

- Cost effective emergency carry litter
- Features ten large carry handles with a load capacity tested to 750 lb. (340 kg)
- Constructed of a heavy-duty polypropylene
- Vacuum-packed to conserve space and fit in a cargo pocket
- Latex free, water-resistant
- Grey nylon webbing & stitching

DIMENSIONS

- Litter Bed: L 74 in. x W 28 in.
- Packaged: L 7 in. x W 6 in. x D 1.75 in.
- Weight: 13 oz
- Load rating: 750 lb

Learn more* about this and other *Products with a Mission*® at *www.NARescue.com





NAR® ELITE LITTER™ - RESPONDER LITE

ITEM #
60-0071

** Explore other soft litters and patient evacuation solutions at* **www.NARescue.com**