NARCO POUCH WITH MEDICATION PANEL

The NARCO Pouch is a dedicated carry and protection system for the NARCO Medication Panel. It keeps your organized medication setup contained, shielded from the environment, and mounted where you can reach it, so the panel workflow stays intact from storage to point of care.

The NARCO Medication Panel is a rugged, quick-access panel built for organization and rapid deployment in the most critical environments. Built for compatibility with our Armadillo Medication Storage Cases (sold separately), the NARCO Panels are a lightweight addition for improved organization and efficiency.

Pouch and Medication Panel is manufactured by Blue Tide Innovations.



Learn more* about this and other *Products with a Mission*® at *www.NARescue.com



POUCH FEATURES/SPECIFICATIONS

- **Purpose-built housing**. Designed specifically to fit and protect the NARCO Medication Panel with precision
- **Durable, water-resistant shell.** 420D TPU construction delivers rugged resilience while keeping a lightweight, streamlined profile for operational efficiency
- **YKK AquaGuard zippers.** Provide a clean seal that adds an extra layer of water-resistant protection
- **Full clamshell access.** Opens completely so the entire panel and contents are visible and reachable
- **Laser-cut MOLLE rear.** Mounts modularly to packs, vests, and platforms
- **Dual-sided loop Velcro interior.** Secures panels and accessories on both faces

PANEL FEATURES/SPECIFICATIONS

- Keeps critical medications secure during tactical operations
- Rapid-access tabs
- Two-panel system (QuickDraw and ADMIN)
- Berry Amendment compliant
- Medications not included with panels
- Compatible with Armadillo Medication Storage Cases (sold separately)

DIMENSIONS

Kit

- Total Dimensions: H 7 in. x L 4.5 in. x W 2 in.
- Total Weight: 0.4 lb.

Pouch

- Pouch Weight: 0.23 lb. (3.6 oz.)

Panels

- Panel Dimensions: L 6.4 in. x W 3.8 in.
- Panel Weight: 0.18 lb. (2.8 oz.)

NARCO POUCH W/ PANELS

ITEM #	NSN#
80-2016	----

˙ Explore more NARCO products at www.narescue.com

NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277